# C-REACTIVE PROTEIN HIGH SENSITIVE

Cat. No.	Pack Name	Packaging (Content)
XSYS0084	CRP-HS	R1: 2 × 40 ml (Buffer), R2: 2 × 11 ml (Latex)

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of C-Reactive Protein in serum and plasma.

## CLINICAL SIGNIFICANCE

### Diagnostic Implications

C-Reactive Protein (CRP) is an acute marker of inflammatory processes. In case of an acute inflammation the concentration of CRP increases and decreases more quickly than the red cell sedimentation rate.

The increase of CRP occurs in a non-specific way in different kinds of tissular aggression, as for example in infectious states, rheumatoid arthritis, myocard infarct, malignant tumour, etc.

Although not diagnostic it is very useful for following-up and monitoring such illnesses, as well as for differential diagnosis in certain cases.

Routinely available immunochemical assay methods for CRP have limited sensitivity, and until recently, CRP concentrations below 10 mg/l could not be measured precisely, leading to the wide spread adoption of this value as the upper limit of the health-associated reference range. This is satisfactory for most purposes in general medicine.

However, in neonatal pediatric practice, a high sensitive CRP immunoassay shows that health-associated reference values are below 1–2 mg/l and that any rise above such values is associated with serious disease, usually bacterial infection.

More recently, application of sensitive CRP assays to studies of adult cardiovascular disease has revealed important prognostic relationships between modest increase of CRP and the occurrence, progression, and thrombo-occlusive complications of atherosclerosis. We therefore developed a high sensitive CRP assay with a detection limit around 0.13 mg/l and a high measuring range (0–140 mg/l HS CRP).

## METHODOLOGY

Measurement of antigen-antibody reaction by the end-point method.

## REAGENT COMPOSITION

### R1 (Buffer)

Sodium chloride	9 g/l
Detergent	0.1 %
Sodium azide	0.09 %

### R2 (Latex)

Glycine buffer pH 8.42	
Rabbit anti-human CRP	sensitized latex 0.20 %
Sodium azide	0.09 %

## REAGENT PREPARATION

Liquid reagents, ready to use.

## STABILITY AND STORAGE

The reagents are stable until expiry date when kept at 2–8 °C. Stability in the instrument is at least 4 weeks if contamination is avoided. Do not freeze.

## WARNING AND PRECAUTIONS

- For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents of the kit are not classified as dangerous.

## SAMPLE COLLECTION

Use fresh serum. If the test can not be carried out on the same day, the serum may be stored at 2–8 °C for 48 hours. If stored for a longer period, the sample should be frozen.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance accurately at 340 nm may be used.
- Analyser specific consumables such as sample cups.
- Controls.
- Saline (9 g/l NaCl)

## ASSAY PROCEDURE

Refer to the assay parameters for details.

## CALIBRATION

Blank: Saline

Cat. No.:	Pack name	Content
XSYS0088	CRP-HS CAL	1 × 1 ml

**Calibraton curve:** generate a 6 point calibration curve by diluting the calibrator 1:32, 1:16, 1:8, 1:4, 1:2 and undiluted in saline.

## Calibration frequency

Calibration verification: Not necessary.

## QUALITY CONTROL

For quality control use

Cat. No.:	Pack name	Content
BLT20034	MULTICON L1	1 × 1 ml
BLT20035	MULTICON L2	1 × 1 ml
BLT20014	CRP CON L	1 × 1 ml

## CALCULATION

Results are calculated automatically by the instrument.

## EXPECTED VALUES

### Reference Values

Less than 1.0 mg/l = Low Risk for CVD  
1.0–2.9 mg/l = Intermediate Risk for CVD  
Greater than 3.0 mg/l = High Risk for CVD

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

## Accuracy

A comparison between XL-Systems CRP (y) and a commercially available test (x) using 40 samples gave following results:

$y = 1.069x - 0.29$  mg/l  $r = 0.999$   
**Measuring Range:** 0.73–160 mg/l  
**Detection Limit:** 0.24 mg/l  
**Hook Effect:** No risk

## PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/l)	SD (mg/l)	CV (%)
Sample 1	75.26	0.86	1.14
Sample 2	35.50	0.79	2.21

Inter-assay precision Run to run (n=20)	Mean (mg/l)	SD (mg/l)	CV (%)
Sample 1	4.20	0.19	4.63
Sample 2	29.85	0.78	2.62
Sample 3	46.13	0.94	2.04

**Specificity:** Monospecific  
**Interferences:** No interference for: Hemoglobin (10 g/l), rheumatic factor (560 IU/ml), Heparin (50 mg/dl), Bilirubin (40 mg/dl), Triglyceride (500 mg/dl)

**Limitations:** None  
**Stability at 4 °C:** At least 3 years after production

## WASTE MANAGEMENT

Please refer to local legal requirements.

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**ASSAY PARAMETERS (conventional units)**

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
<b>Test Details</b>						
Test	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS
Test Code	70	70	70	70	70	70
Report Name	C-Reactive Protein HS	C-Reactive Protein HS	C-Reactive Protein HS	C-Reactive Protein HS	C-Reactive Protein HS	C-Reactive Protein HS
Unit	mg/dl	mg/dl	mg/dl	mg/dl	mg/dl	mg/dl
Decimal Places	2	2	2	2	2	2
Wavelength-Primary	600	600	600	600	600	600
Wavelength-Secondary	0	0	0	0	0	0
Assay type	2-Point	2-Point	2-Point	2-Point	2-Point	2-Point
Curve type	5P Calibration Logit - Log	5P Calibration Logit - Log	5P Calibration Logit - Log	5P Calibration Logit - Log	5P Calibration Logit - Log	5P Calibration Logit - Log
M1 Start	16	16	12	24	10	16
M1 End	16	16	12	24	10	16
M2 Start	34	36	51	63	31	34
M2 End	34	36	51	63	31	34
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	NA	NA	NA	NA	NA	NA
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	NA	NA	NA	NA	NA	NA
Technical Maximum	NA	NA	NA	NA	NA	NA
<b>Y=aX+b</b>						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	NA	NA	NA	NA	NA	NA
Reagent Abs Max	0	0	0	0	0	0
Auto Rerun	No	No	No	No	No	No
Total Reagents	2	2	2	2	2	2
Reagent R1	CRPHS R1	CRPHS R1	CRPHS R1	CRPHS R1	CRPHS R1	CRPHS R1
Reagent R2	CRPHS R2	CRPHS R2	CRPHS R2	CRPHS R2	CRPHS R2	CRPHS R2
Reagent R3	NA	NA	NA	NA	NA	NA
<b>Test Volumes</b>						
Test	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
<b>Sample Volumes</b>						
Normal	2	2	2	2	2	2
Dilution Ratio	1	1	1	1	1	1
Increase	4	4	4	4	4	4
Dilution Ratio	1	1	1	1	1	1
Decrease	2	2	2	2	2	2
Dilution Ratio	10	10	10	10	10	10
Standard volume	2	2	2	2	2	2
<b>Reagent Volumes and Stirrer speed</b>						
RGT-1 Volume	250	250	250	250	250	250
R1 Stirrer Speed	High	High	NA	High	High	High
RGT-2 Volume	60	60	60	60	60	60
R2 Stirrer Speed	High	High	NA	High	High	High
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
<b>Reference Ranges</b>						
Test	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
<b>Category Male</b>						
Normal-Lower Limit	0	0	0	0	0	0
Normal-Upper Limit	0.3	0.3	0.3	0.3	0.3	0.3
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
<b>Category Female</b>						
Normal-Lower Limit	0	0	0	0	0	0
Normal-Upper Limit	0.3	0.3	0.3	0.3	0.3	0.3
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
<b>Revision Number</b>						
Revision	<A-100- CRPHS-2 26.06.2015>	<A-200- CRPHS-2 26.06.2015>	<A-300/600- CRPHS-2 26.06.2015>	<A-640- CRPHS-2 26.06.2015>	<A-1000- CRPHS-2 26.06.2015>	<A-180- CRPHS-2 26.06.2015>

**ASSAY PARAMETERS (SI units)**

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
<b>Test Details</b>						
Test	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS
Test Code	70	70	70	70	70	70
Report Name	C-Reactive Protein HS	C-Reactive Protein HS	C-Reactive Protein HS	C-Reactive Protein HS	C-Reactive Protein HS	C-Reactive Protein HS
Unit	mg/l	mg/l	mg/l	mg/l	mg/l	mg/l
Decimal Places	1	1	1	1	1	1
Wavelength-Primary	600	600	600	600	600	600
Wavelength-Secondary	0	0	0	0	0	0
Assay type	2-Point	2-Point	2-Point	2-Point	2-Point	2-Point
Curve type	5P Calibration Logit - Log	5P Calibration Logit - Log	5P Calibration Logit - Log	5P Calibration Logit - Log	5P Calibration Logit - Log	5P Calibration Logit - Log
M1 Start	16	16	12	24	10	16
M1 End	16	16	12	24	10	16
M2 Start	34	36	51	63	31	34
M2 End	34	36	51	63	31	34
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	NA	NA	NA	NA	NA	NA
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	NA	NA	NA	NA	NA	NA
Technical Maximum	NA	NA	NA	NA	NA	NA
<b>Y=aX+b</b>						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	NA	NA	NA	NA	NA	NA
Reagent Abs Max	0	0	0	0	0	0
Auto Rerun	No	No	No	No	No	No
Total Reagents	2	2	2	2	2	2
Reagent R1	CRPHS R1	CRPHS R1	CRPHS R1	CRPHS R1	CRPHS R1	CRPHS R1
Reagent R2	CRPHS R2	CRPHS R2	CRPHS R2	CRPHS R2	CRPHS R2	CRPHS R2
Reagent R3	NA	NA	NA	NA	NA	NA
<b>Test Volumes</b>						
Test	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
<b>Sample Volumes</b>						
Normal	2	2	2	2	2	2
Dilution Ratio	1	1	1	1	1	1
Increase	4	4	4	4	4	4
Dilution Ratio	1	1	1	1	1	1
Decrease	2	2	2	2	2	2
Dilution Ratio	10	10	10	10	10	10
Standard volume	2	2	2	2	2	2
<b>Reagent Volumes and Stirrer speed</b>						
RGT-1 Volume	250	250	250	250	250	250
R1 Stirrer Speed	High	High	NA	High	High	High
RGT-2 Volume	60	60	60	60	60	60
R2 Stirrer Speed	High	High	NA	High	High	High
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
<b>Reference Ranges</b>						
Test	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS	CRPHS
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
<b>Category Male</b>						
Normal-Lower Limit	0	0	0	0	0	0
Normal-Upper Limit	3	3	3	3	3	3
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
<b>Category Female</b>						
Normal-Lower Limit	0	0	0	0	0	0
Normal-Upper Limit	3	3	3	3	3	3
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
<b>Revision Number</b>						
Revision	<ASI-100- CRPHS-2 26.06.2015>	<ASI-200- CRPHS-2 26.06.2015>	<ASI-300/600- CRPHS-2 26.06.2015>	<ASI-640- CRPHS-2 26.06.2015>	<ASI-1000- CRPHS-2 26.06.2015>	<ASI-180- CRPHS-2 26.06.2015>

**REFERENCES / ЛИТЕРАТУРА / ЛІТЕРАТУРА / LITERATURA / LITERATÚRA / REFERENCIAS**

1. Claus DR, Osmand AP, Gewurz H. Radioimmunoassay of human C-reactive protein and levels in normal sera. J Lab Clin Med 1976; 87: 120-128
2. Wasunna A, Whitelaw A, Gallimore R, Hawkins PN, Pepys MB. C-reactive protein and bacterial infection in preterm infants. Eur J Pediatr 1990; 149:424-427
3. Heinrich J, Schulte H, Schönfeld R, Köhler E, Assmann G. Association of variables of coagulation, fibrinolysis and acute-phase with atherosclerosis in coronary and peripheral arteries and those arteries supplying the brain. Thromb Haemostas 1995; 73: 374-379
4. Grau AJ, Buggle F, Becher H, Werle E, Hacke W. The association of leukocyte count, fibrinogen and C-reactive protein with vascular risk factors and ischemic vascular diseases. Thromb Res 1996; 82:245-255
5. Dixon JS, Bird HA, Sitton NG et al. C-Reactive protein in the serial assessment of disease activity in rheumatoid arthritis. Scand J Rheum 1984; 13: 39-44
6. Hind CRH, and Pepys MB. The role of serum C-Reactive protein (CRP) measurement in clinical practice. Int Med 1984; 5: 112-151
7. Hanson LÅ, Wadsworth Ch. Das C-reaktive Protein und sein diagnostischer Wert, insbesondere bei infektionen. Laboratoriumblätter 1979; 29: 58-68
8. Macy EM, Hayes TE, Tracy RP. Variability in the measurement of C-reactive protein in healthy subjects: implications reference intervals and epidemiological applications. Clin Chem 1997; 43: 52-58

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