HbA1c 2R

Cat. No.	Pack Name	Packaging (Content)
XSYS0096	HbA1c 2R	R1: 2 × 21 ml, R2: 2 × 8 ml, R3: 3 × 50 ml

EN

IVD

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of glycated hemoglobin HbA1c in whole blood on photometric systems.

CLINICAL SIGNIFICANCE

Hemoglobin A1c (HbA1c) is a glycated hemoglobin which is formed by the non-enzymatic reaction of glucose with native hemoglobin. This process runs continuously throughout the circulatory life of the red cell (average life time 100–120 days). The rate of glycation is directly proportional to the concentration of glucose in the blood. The blood level of HbA1c represents the average blood glucose level over the preceding 6 to 8 weeks (due to the kinetics of erythrocyte turnover this period is more affected by the blood glucose level than the preceding weeks). Therefore, HbA1c is suitable for retrospective long-term monitoring of blood glucose concentration in individuals with diabetes mellitus. Clinical studies have shown that lowering of HbA1c level can help to prevent or delay the incidence of late diabetic complications.

As the amount of HbA1c also depends on the total quantity of hemoglobin the reported HbA1c value is indicated as a percentage of the total hemoglobin concentration. Falsely low values (low HbA1c despite high blood glucose) may occur in people with conditions with shortened red blood cell survival (hemolytic diseases) or significant recent blood loss (higher fraction of young erythrocytes). Falsely high values (high HbA1c despite normal blood glucose) have been reported in iron deficiency anemia (high proportion of old erythrocytes). These circumstances have to be considered in clinical interpretation of HbA1c values.

PRINCIPLE

Particle enhanced immunoturbidimetric test.

HbA1c is determined directly without measurement of total hemoglobin.

This test is based on latex immunoagglutination where HbA1c in the test sample is absorbed onto latex particles, and then Cross-linked anti-HbA1c is added to form an Antigen-Antibody reaction. The measured absorbance is proportional to the HbA1c bound to particles, which in turn is proportional to the percentage of HbA1c in the sample.

STANDARDIZATION

The assay is standardized according to the approved IFCC reference method [3]. Calibration according to DCCT/NGSP is also possible. Corresponding calibrator values are listed in the package insert of the calibrator set HbA1c liquid. DCCT/NGSP and IFCC values show a linear relationship and can therefore be calculated from each other. Also new IFCC recommended units mmol/mol (mmol HbA1c / mol Hb) can be easily calculated from IFCC values. Recalculation equations are following:

IFCC = (NGSP – 2.15) / 0.915

NGSP = 0.915 × IFCC + 2.15

mmol/mol = 10 × IFCC

IFCC: International Federation of Clinical Chemistry [3,4]

DCCT: Diabetes Control and Complications Trial [5]

NGSP: National Glycohemoglobin Standardization Program [6]

REAGENT COMPOSITION

R1:	Latex	0.1 %
R2:	Anti-HbA1c Cross linked anti-human hemoglobin HbA1c mouse monoclonal antibody	
	Buffer	
	NaCl	2 %
R3:	Hemolyzing solution	

REAGENT PREPARATION

All reagents are ready for use. Mix each reagent well before use.

STABILITY AND STORAGE

The reagents are stable until expiry date when kept at 2–8 °C. Stability after opening and installation to the instrument is at least 4 weeks if contamination is avoided. Do not freeze.

SPECIMEN COLLECTION

Whole blood collected with EDTA.

Sample preparation:

Hemolyzing solution (R3) 500 µl

Sample 10 µl

Mix and allow to stand for 5 minutes or until complete lysis is apparent.

Sample preparation on-board is not recommended.

Specimen stability:

Whole blood 1 week at 2–8 °C

Hemolysate 10 hours at 15–25 °C

Hemolysate 10 days at 2–8 °C

ASSAY PROCEDURE

Application sheets for automated systems are available on request.

Wavelength 660 nm

Optical path 1 cm

Temperature 37 °C

Measurement Against air

Sample or calibrator	36 µl
Reagent 1	900 µl
Mix, incubate for 5 min., read absorbance 1, then add:	
Reagent 2	300 µl
Mix, read absorbance 2 after exactly 5 min.	

MATERIAL REQUIRED BUT NOT PROVIDED

General laboratory equipment.

CALIBRATION

The concentration of HbA1c in unknown samples is derived from a calibration curve using an appropriate mathematical model such as spline. The calibration curve is obtained with 5 calibrators at different levels, where the first calibrator with zero value is used as blank.

Stability of calibration: 4 weeks

For calibration use the HbA1c 2R CALIBRATOR SET (XSYS0097).

CONTROLS

For internal quality control use the HbA1c 2R CONTROL LOW (XSYS0098) and HbA1c 2R CONTROL HIGH (XSYS0099).

CALCULATION

Results are calculated automatically by the instrument.

EXPECTED VALUES [7]

Reference intervals should be established or verified by the laboratory based on an appropriate non-diabetic patient population.

	% NGSP	% IFCC	mmol/mol
Non-diabetics	4–6	3–4	30–40
Target of therapy	< 7	< 5	< 50
Change of therapy	> 8	> 6	> 60

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of detection: 0.24 % NGSP

Limit of quantification: 0.73 % NGSP

Linearity: 15.6 % NGSP

Measuring range: 0.73–15.6 % NGSP

Accuracy:

	Target (% NGSP)	Result (% NGSP)	bias (%)
CONTROL LOW	5.7	5.42	-4.9
CONTROL HIGH	13.3	14.19	+6.7

Precision:

Intra-assay n = 20	Mean (% NGSP)	SD (% NGSP)	CV (% NGSP)
Sample 1	6.72	0.04	0.65
Sample 2	11.08	0.06	0.58

Inter-assay n = 20	Mean (% NGSP)	SD (% NGSP)	CV (% NGSP)
Sample 1	6.67	0.05	0.81
Sample 2	11.07	0.26	2.39

Method Comparison:

A comparison between XL-Systems HbA1c (y) and commercially available test (x) using 70 samples were carried. Whole blood samples were assayed in parallel and the results compared by linear regression analysis according to Passing-Bablok. The following statistics were obtained:

y = 0.9828 × - 0.048

r = 0.972

Specificity/Interferences:

Due to its antibodies, HbA1c 2R is a specific immunoassay for human HbA1c.

No interference was observed by ascorbic acid up to 50 mg/dl, bilirubin F up to 19.6 mg/dl, bilirubin C up to 20.6 mg/dl, chyle up to 1590 FTU, rheumatoid factor up to 1100 IU/ml and high lipid up to 1400 formazine turbidity units.

Testing of interferences was performed on Hitachi 7180.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1, R2

EUH 208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction

R3

Reagent is not classified as dangerous.

Immediately after HbA1c measurement cleaning of cuvettes is necessary. Use the alkaline cuvette washing solution which is recommended by the analyzer manufacturer. Take necessary precautions for use of laboratory reagent.

WASTE MANAGEMENT

Please refer to local legal requirements.

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ASSAY PARAMETERS (conventional units)

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details						
Test	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R
Test Code	8	8	8	8	8	8
Report Name	HbA1c	HbA1c	HbA1c	HbA1c	HbA1c	HbA1c
Unit	%	%	%	%	%	%
Decimal Places	1	1	1	1	1	1
Wavelength-Primary	660	660	660	660	660	660
Wavelength-Secondary	0	0	0	0	0	0
Assay type	2-Point	2-Point	2-Point	2-Point	2-Point	2-Point
Curve type	Point-To-Point	Point-To-Point	Point-To-Point	Point-To-Point	Point-To-Point	Point-To-Point
M1 Start	20	18	17	27	16	18
M1 End	20	18	17	27	16	18
M2 Start	34	36	50	62	31	34
M2 End	34	36	50	62	31	34
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	0	0	0	0	0	0
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	NA	NA	NA	NA	NA	NA
Technical Maximum	NA	NA	NA	NA	NA	NA
Y=aX+b						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	NA	NA	NA	NA	NA	NA
Reagent Abs Max	0	0	0	0	0	0
Auto Rerun	No	No	No	No	No	No
Total Reagents	2	2	2	2	2	2
Reagent R1	A1C R1	A1C R1	A1C R1	A1C R1	A1C R1	A1C R1
Reagent R2	A1C R2	A1C R2	A1C R2	A1C R2	A1C R2	A1C R2
Reagent R3	NA	NA	NA	NA	NA	NA
Test Volumes						
Test	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R
Sample Type	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD
Sample Volumes						
Normal	8.4	7.2	8.4	7.2	6	8.4
Dilution Ratio	1	1	1	1	1	1
Increase	16.8	14.4	16.8	14.4	12	16.8
Dilution Ratio	1	1	1	1	1	1
Decrease	4.2	3.6	4.2	3.6	3	4.2
Dilution Ratio	1	1	1	1	1	1
Standard volume	8.4	7.2	8.4	7.2	6	8.4
Reagent Volumes and Stirrer speed						
RGT-1 Volume	210	180	210	180	150	210
R1 Stirrer Speed	High	High	High	High	High	High
RGT-2 Volume	70	60	70	60	50	70
R2 Stirrer Speed	High	High	High	High	High	High
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
Reference Ranges						
Test	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R
Sample Type	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD
Reference Range	Default	Default	Default	Default	Default	Default
Category Male						
Normal-Lower Limit	4	4	4	4	4	4
Normal-Upper Limit	6	6	6	6	6	6
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Category Female						
Normal-Lower Limit	4	4	4	4	4	4
Normal-Upper Limit	6	6	6	6	6	6
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Revision Number						
Revision	<A-100- A1C2R-1 03.03.2020>	<A-200- A1C2R-1 03.03.2020>	<A-300/600- A1C2R-1 03.03.2020>	<A-640- A1C2R-1 03.03.2020>	<A-1000- A1C2R-1 03.03.2020>	<A-180- A1C2R-1 03.03.2020>

ASSAY PARAMETERS (SI units)

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details						
Test	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R
Test Code	8	8	8	8	8	8
Report Name	HbA1c	HbA1c	HbA1c	HbA1c	HbA1c	HbA1c
Unit	mmol/mol	mmol/mol	mmol/mol	mmol/mol	mmol/mol	mmol/mol
Decimal Places	1	1	1	1	1	1
Wavelength-Primary	660	660	660	660	660	660
Wavelength-Secondary	0	0	0	0	0	0
Assay type	2-Point	2-Point	2-Point	2-Point	2-Point	2-Point
Curve type	Point-To-Point	Point-To-Point	Point-To-Point	Point-To-Point	Point-To-Point	Point-To-Point
M1 Start	20	18	17	27	16	18
M1 End	20	18	17	27	16	18
M2 Start	34	36	50	62	31	34
M2 End	34	36	50	62	31	34
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	0	0	0	0	0	0
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	NA	NA	NA	NA	NA	NA
Technical Maximum	NA	NA	NA	NA	NA	NA
Y=aX+b						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	NA	NA	NA	NA	NA	NA
Reagent Abs Max	0	0	0	0	0	0
Auto Rerun	No	No	No	No	No	No
Total Reagents	2	2	2	2	2	2
Reagent R1	A1C R1	A1C R1	A1C R1	A1C R1	A1C R1	A1C R1
Reagent R2	A1C R2	A1C R2	A1C R2	A1C R2	A1C R2	A1C R2
Reagent R3	NA	NA	NA	NA	NA	NA
Test Volumes						
Test	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R
Sample Type	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD
Sample Volumes						
Normal	8.4	7.2	8.4	7.2	6	8.4
Dilution Ratio	1	1	1	1	1	1
Increase	16.8	14.4	16.8	14.4	12	16.8
Dilution Ratio	1	1	1	1	1	1
Decrease	4.2	3.6	4.2	3.6	3	4.2
Dilution Ratio	1	1	1	1	1	1
Standard volume	8.4	7.2	8.4	7.2	6	8.4
Reagent Volumes and Stirrer speed						
RGT-1 Volume	210	180	210	180	150	210
R1 Stirrer Speed	High	High	High	High	High	High
RGT-2 Volume	70	60	70	60	50	70
R2 Stirrer Speed	High	High	High	High	High	High
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
Reference Ranges						
Test	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R	A1C2R
Sample Type	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD	BLOOD
Reference Range	Default	Default	Default	Default	Default	Default
Category Male						
Normal-Lower Limit	20	20	20	20	20	20
Normal-Upper Limit	42	42	42	42	42	42
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Category Female						
Normal-Lower Limit	20	20	20	20	20	20
Normal-Upper Limit	42	42	42	42	42	42
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Revision Number						
Revision	<ASI-100- A1C2R-1 03.03.2020>	<ASI-200- A1C2R-1 03.03.2020>	<ASI-300/600- A1C2R-1 03.03.2020>	<ASI-640- A1C2R-1 03.03.2020>	<ASI-1000- A1C2R-1 03.03.2020>	<ASI-180- A1C2R-1 03.03.2020>

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1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 142-48.
2. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B. Saunders Company; 1999.p. 790-6.
3. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002;40:78-89.
4. Hoelzel W, Weykamp C et al. IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardization Schemes in the United States, Japan, and Sweden: A Method-Comparison Study. Clin Chem 2004; 50:1:166-74.
5. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med.1993;329:977-86.
6. Little RR, Rohlfing CL, Wiedmeyer HM, Myers GL et al. The National Glycohemoglobin Standardization Program: A Five-Years Progress Report. Clin Chem 2001;47:1985-92.
7. Miedema K. Standardization of HbA1c and Optimal Range of Monitoring. Scand J Clin Lab Invest 2005;65 (Suppl 240):61-72

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