BICARBONATE

Cat. No.	Pack Name	Packaging (Content)
XSYS0100	CO2	R1: 4 x 34 ml, R2 STD: 1 x 5 ml

EN

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Bicarbonate (CO₂) in human serum and plasma.

CLINICAL SIGNIFICANCE

Measurement of bicarbonate is used in the diagnosis of the acid-base-balance in the blood. Elevated and decreased values indicate disorders associated with disturbances of the metabolic and respiratory systems.

PRINCIPLE

Bicarbonate reacts with phosphoenolpyruvate (PEP), in the presence of phosphoenolpyruvate carboxylase (PEP-C), to form oxaloacetate and phosphate.

The oxaloacetate is then converted to malate by the action of malate dehydrogenase (MDH) and reduced nicotinamide adenine dinucleotide analog (NADH-analog).

The decrease in absorbance at 405 or 415 nm resulting from the oxidation of NADH analog is proportional to the amount of CO₂ in the sample.

REAGENT COMPOSITION

R1

Buffer	0.1 M
PEP	10 mM
PEC-C	> 100 U/l
MDH	> 1000 U/l
NAD-analog	0.5 mM
stabilizer and preservative.	

R2

Standard: bicarbonate solution 30 mmol/l.

REAGENT PREPARATION

Reagents are liquid, ready for use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

After first opening, reagents are stable for 60 days at (+2 to +8) °C if stored at appropriate conditions, closed carefully and without any contamination.

On board stability: min. 28 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum, heparin plasma.

The preferred specimen is from venous blood collected anaerobically in the usual manner for bicarbonate analysis. Bicarbonate content in uncapped tubes decreases approximately 4 mmol/l after one hour. It has been reported that alkalinized serum stored in open cups is stable for up to 4 hours.

Separate from erythrocytes and store tightly stoppered.

It is recommended to follow NCCLS procedures (or similar standardized conditions). Specimens should be protected from direct exposure to light.

Stability of bicarbonate in the dark and closed vials:

3 days at 4–8 °C

1 month at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with the standard included in the kit is recommended.

Calibration frequency: it is recommended to do a calibration:

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This method has been standardized against a primary standard on basis of sodium carbonate.

QUALITY CONTROL

For quality control all control solutions with CO₂ values determined by this method can be used.

CALCULATION

The XL Results are calculated automatically by the instrument.

UNIT CONVERSION

mmol/l = mEq/l

EXPECTED VALUES ¹

Serum, plasma:

Newborn:	13–22 mmol/l – mEq/l
Infant, child:	20–28 mmol/l – mEq/l
Adults:	23–29 mmol/l – mEq/l
Adults > 60yr:	23–31 mmol/l – mEq/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 1.4 mmol/l

Linearity: 50 mmol/l

Measuring range: 1.4–50 mmol/l

PRECISION

Intra-assay (n = 20)	Mean [mmol/l]	SD [mmol/l]	CV (%)
Sample 1	22.9	0.461	2.01
Sample 2	41.2	0.625	1.52

Inter-assay (n = 20)	Mean [mmol/l]	SD [mmol/l]	CV (%)
Sample 1	34.9	0.347	0.99
Sample 2	23.8	0.672	2.83

COMPARISON

A comparison between XL-Systems CO₂ (y) and a commercially available test (x) using 86 samples gave following results:

N = 40

r = 0.986

y = 0.9303 x – 0.0479 mmol/l

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 850 mg/dl

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents of the kit are not classified like dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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