EC cartridge plus iCa/60D

Cat. No.:	Pack Name:	Number of Tests:
REG00070	EC cartridge XS plus iCa/60D	250
REG00068	EC cartridge S plus iCa/60D	500
REG00073	EC cartridge M plus iCa/60D	1000
REG00074	EC cartridge L plus iCa/60D	3000

EN

INTENDED USE

EC cartridge is an IVD kit intended for calibration and quantification of sodium (Na⁺), potassium (K⁺), ionized calcium (iCa²⁺) and chloride (Cl⁻) on the EC 90 analyzer series. Intended only for IVD use by a professionally trained person.

CONTENT

1) EC cartridge

Calibration solution 1:

Na⁺ 150 mmol/l, K⁺ 5 mmol/l, iCa²⁺ 1,25 mmol/l, Cl⁻ 103 mmol/l, buffer, preservative, wetting agent

Calibration solution 2:

Na⁺ 105 mmol/l, K⁺ 2 mmol/l, iCa²⁺ 2,34 mmol/l, Cl⁻ 71,1 mmol/l, buffer, preservative, wetting agent

Reference solution:

K⁺ 1,2 mol/l, buffer, preservative, wetting agent

2) EC Biosensor

PREPARATION

Insert EC Biosensor into the cartridge according to the user manual.

DIRECTIONS FOR USE

1. In the analyzer software, in the maintenance screen, click on CHANGE EC CARTRIDGE button. Follow instructions on the screen.
2. Unpack the new cartridge, remove the cap and seal cap from the top and place it in the analyzer.
3. In the analyzer software click on OK button. The software will detect new cartridge and perform priming automatically.
4. Calibration will be performed automatically.

STABILITY AND STORAGE

Store at 2–30 °C. Note the expiration date on the labels. Protect the cartridge from frost. On board stability: 60 days after opening.

QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat. No. BLT00080) or ERBA NORM 10x5 (Cat. No. XSYS0123) and ERBA PATH 4x5 (Cat. No. BLT00081) or ERBA PATH 10x5 (Cat. No. XSYS0124) are recommended.

PERFORMANCE DATA

Performance is included in the user manual.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Calibration solution 1 and 2:

EUH 208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

Reference solution: The reagent is not classified as dangerous.

WASTE MANAGEMENT

Users should be exceedingly careful when working with this biological material and treat the waste as a potential infectious material. Please refer to local legal requirements.

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