## Opti-sol XL 750 Liquid System

Cat. No.	Product Name	Packaging (Content)
REG00054	Opti-sol XL 750	6 × 50 mL, Instruction for use, RFID tag

### INTENDED USE

Opti-sol XL 750 is the washing solution for hydraulic system of automatic urine analyzer LAURA XL intended for *in vitro* diagnostics. For professional use at clinical laboratories only.

### PRINCIPLE

The performance of analyser LAURA XL depends on stabilization and cleanest purity of optical and flow parts. Opti-sol XL 750 solution contains additives which rinses the hydraulic system between samples in the LAURA XL automatic analyser.

### DESCRIPTION AND COMPOSITION

Purified water	70 %
Detergents, preservatives, and remaining substances	30 %

### WORKING SOLUTION PREPARATION

The content of the bottle of Opti-sol XL 750 dilute in 10 litres of purified water (50 mL of Opti-sol XL 750 + 9950 mL of water) into provided tank and attach to LAURA XL analyser.

### MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

LAURA XL automatic urine analyser, Cat. No. INS00065

### STABILITY AND STORAGE

The washing solution is stable till the expiry date mentioned on the label. The washing solution should be stored in tightly closed original packaging in a dry and dark place at a temperature 2–25 °C.

### WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

#### Liquid System

UFI: D1KX-TJY3-RE7X-NRPC

Warning

**Contains:** 2-Methyl-2H-isothiazol-3-one

#### Hazard statement:

H315 Causes skin irritation.  
H317 May cause an allergic skin reaction.  
H319 Causes serious eye irritation.  
H412 Harmful to aquatic life with long lasting effects.

#### Precautionary statement:

P264 Wash hands and exposed parts of the body thoroughly after handling.  
P280 Wear protective gloves/protective clothing/eye protection.  
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.  
P337+P313 If eye irritation persists: Get medical advice/attention.

#### Working solution:

UFI: 4WWD-K67E-W46N-1ARW

Warning

**Contains:** 2-Methyl-2H-isothiazol-3-one

#### Hazard statement:

H317 May cause an allergic skin reaction.

#### Precautionary statement:

P280 Wear protective gloves.  
P302+P352 IF ON SKIN: Wash with plenty of water and soap.  
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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