Erba H7100 Fluro Ret

Haematology Analyser Staining Agent

Cat. No.	Pack Name	Packaging (Content)
HEM00051	Erba H7100 Fluro Ret (12 mL)	12 mL, Instruction for use, RFID card
HEM00052	Erba H7100 Fluro Ret (24 mL)	24 mL, Instruction for use, RFID card
HEM00068	Erba H7100 Fluro Ret (6 mL)	6 mL, Instruction for use, RFID card

INTENDED USE

Erba H7100 Fluro Ret is an *in vitro* diagnostic reagent used to stain blood cells to observe their morphology and structure for reticulocyte count on automated haematology analysers H7100. It is designed to report on the haematological status of patient samples, supporting the monitoring and diagnosis of pathological conditions. For professional use in clinical laboratory only.

PRINCIPLE

The product is used in conjunction with Erba H7100 Dil Ret to perform reticulocytes differential count by diluting the blood, making the cells spherical, and staining the reticulocytes in sample.

DESCRIPTION AND COMPOSITION

Nucleic acid fluorescent dye; Isopropyl alcohol; Ethylene glycol.

PREPARATION

This reagent is ready for use and can be applied straight from the container; no special reagent preparation is necessary. Avoid reagent contaminations. Before using, leave the reagent at room temperature for at least 24 hours. For reagent connection to the analyser, refer to the H7100 User manual.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Haematology analyser H7100, Cat. No. INS00087
Reagents: Erba H7100 Dil, Cat. No. HEM00044
Erba H7100 Dil Ret, Cat. No. HEM00045, HEM00065
Erba H7100 Lyse1, Cat. No. HEM00046, HEM00047, HEM00066
Erba H7100 Lyse2, Cat. No. HEM00048, HEM00049, HEM00050, HEM00067
Erba H7100 Fluro Diff, Cat. No. HEM00053, HEM00054

SPECIFICATION

The product is a clear liquid without sediment, particles or floccules. Background result: WBC $\leq 0.1 \times 10^9/L$, RBC $\leq 0.02 \times 10^{12}/L$, HGB ≤ 1 g/L, PLT $\leq 5 \times 10^9/L$, HCT ≤ 0.5 %.

STABILITY AND STORAGE

The product is stable till the expiry date mentioned on the label when stored at (+2 to +30) °C. Discard opened container after 90 days.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: S7QQ-PGAT-D35X-Q2N3

Warning

Reagent contains: ethylene glycol

Hazard statements:

H302 Harmful if swallowed.

Precaution statements:

P264 Wash hands and exposed parts of the body thoroughly after handling.

P270 Do not eat, drink or smoke when using this product.

P301 + P312 IF SWALLOWED: Call a doctor if you feel unwell.

WASTE MANAGEMENT

Please refer to local legal requirements.

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