Erba Thrombin Reagent

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Cat. No.	Pack name	Packaging (Content)
EHL00005	Erba Thrombin Reagent	5 × 2 ml
EHL00006	Erba Thrombin Reagent	10 × 2 ml
EHL00025	Erba Thrombin Reagent	10 × 5 ml

POUŽITÍ

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10020367 10020410
10020368 10020411
10020369 10020412

QUALITY SYSTEM CERTIFIED
ISO 13485

Erba Thrombin Reagent For Fibrinogen Determination

Cat. No.	Pack name	Packaging (Content)
EHL00005	Erba Thrombin Reagent	5 × 2 ml
EHL00006	Erba Thrombin Reagent	10 × 2 ml
EHL00025	Erba Thrombin Reagent	10 × 5 ml

INTENDED USE

Erba Thrombin Reagent is intended for use in the quantitative determination of fibrinogen in human plasma using the Clauss Method.

CLINICAL SIGNIFICANCE

Fibrinogen (also called Factor I) is a glycoprotein synthesised in the liver hepatocytes and megakaryocytes. Fibrinogen is composed of two sets of three different chains (α , β , and γ), linked to each other by disulfide bonds. On the alpha and beta chains, there is small peptide sequence (called a fibrinopeptide). Thrombin (Factor IIa) is responsible for converting fibrinogen into fibrin, by splitting out fibrinopeptides monomers. The monomers aggregate to form Fibrin, which is then cross linked by Factor XIII a to form an insoluble fibrin clot.

Fibrinogen is an acute phase reactant protein. Therefore its blood level can be increased as a response of inflammatory states, with infections, during pregnancy, and after trauma. A decrease of the fibrinogen level is found in congenital hypofibrinogenemia and afibrinogenemia. It can also be decreased in cases of Disseminated Intravascular Coagulation (DIC), systemic fibrinolysis, pancreatic and severe hepatic dysfunction.

PRINCIPLE

Clauss 1 developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of dilute plasma after the addition of Thrombin (> 30 NIH units/ml). This clot time is proportional to the fibrinogen concentration.

COMPOSITION

Erba Thrombin Reagent; approximately 100 NIH units/ml bovine thrombin with stabilisers.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

Reagent contains < 30 % bovine thrombin.

Warning

Hazard statement:

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

Precautionary statement:

P302+P352 IF ON SKIN: Wash with plenty of water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

The Erba Thrombin Reagent contains bovine thrombin. Plasma products have been screened and founded negative (unless otherwise stated on the box or vial) for the presence of Hepatitis B Antigen (HBsAg), HIV 1 and 2 antibody and HCV antibody. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

WORKING REAGENT

Reconstitute each vial of Erba Thrombin Reagent with exactly 2 ml (EHL00005/EHL00006) or 5 ml (EHL00025) of purified water. Take care when pipetting to avoid contamination. Swirl gently and allow to stand for 15 minutes. Mix well immediately before use. Do not shake. Avoid reagent contamination

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Once reconstituted, the reagent is stable for:

- 8 hours at 15–30 °C
- 1 week at 2–8 °C
- 1 month at -20 °C.

REAGENTS REQUIRED BUT NOT PROVIDE

Erba Standard Plasma (Cat. No.: EHL00012)

Erba Owrens Veronal Buffer (Cat. No.: EHL00021)

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2 % or 3.8 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept between 2–8 °C or between 18–24 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes.²

PROCEDURE

Manual method

1. Standard Curve Preparation

A new standard curve should be used when reagent lot numbers change or if expected values in quality control testing change. Prepare the following dilutions of Erba Standard Plasma (Cat.No.: EHL00012) in Erba Owrens Veronal Buffer (Cat. No.: EHL00021)

Tube	Dilution	Erba Standard Plasma (ml)	Owren's buffer (ml)
1	1 + 4	0,2	0,8
2	1 + 9	0,1	0,9
3	1 + 19	0,1	1,9
4	1 + 29	0,1	2,9
5	1 + 39	0,1	3,9

2. Patient Sample Preparation:

- Prepare 1 + 9 dilutions of the patient plasma or control plasma in Owrens Buffer.
- Mix without shaking.

3. Testing

- Pipette 0.2 ml of standard, patient or control dilution into a reaction tube and incubate at 37 °C for 2 minutes.
- Add 0.1 ml of Erba Thrombin Reagent (15-30 °C) and start simultaneously a timer.
- Note the time for clot formation to the nearest 0.1 second.
- Plot means Standard Clot Time (Y-axis) versus Fibrinogen Level (X-axis) on the fibrinogen graph paper. A straight line should be obtained. Assign the Erba Standard Plasma (Cat. No.: EHL00012) reference value to the 1 + 9 dilution to allow direct interpolation of patient and control values from the standard curve.

Automated method

Refer to the instrument's operator's manual.

REFERENCES VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should established its own reference ranges.

Usually normal value is interval 170–400 mg/dl.

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015 or Erba Control P Plus, Cat. No.: EHL00017).

LIMITATIONS

Heparin levels > 0.6 units/ml and fibrinolytic degradation product > 100 mg/ml may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, re-assay using an appropriate dilution to bring values into the standard range.

INTERFERENCES

Haemoglobin: No significant interference up to 10 g/l on normal plasma (366.9 mg/dl).

Total Bilirubin: Positive bias from 91 mg/l (155.6 μ mol/l) on normal plasma (373.3 mg/dl).

Turbidity: Positive bias from 3.7 g/l (4.22 mmol/l) equivalent triglycerides on normal plasma (365.1 mg/dl).

PERFORMANCES

These performances have been obtained using an analyzer ECL. Results may vary if a different instrument or manual procedure is used.

Linearity: The Erba Thrombin Reagent assay is designed to give a linear calibration in the range of 100–460 mg/dl.

Precision:

	Intra-assay precision (N = 20)		Inter-assay precision (N = 20)	
Mean (mg/dl)	303.9	144.9	287.7	132.5
CV (%)	4.89	7.19	5.64	5.59

REFERENCES

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur-bestimmung des Fibrinogens, Acta Haematol, 17:237-246
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI:H21-A5
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, Crit.Rev.Clin.Lab.Sci, 8:145-192
- Scully RE et al. (1980) Normal Reference Laboratory Values, New England Journal Medicine, 302(37):37-48
- Okuno T, Selenko V (1972) Plasma fibrinogen determination by automated thrombin time, American Journal of Medical Technology, 38(6):196-201

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