Erba Thrombin Reagent

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LITERATURA

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, Acta Haematol, 17:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, Crit.Rev.Clin.Lab.Sci, 8:145-192
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI:H21-A5
- Hematology of Infancy and Childhood, Nathan and Orkin, Saundersds 2003, P. Monagle: Developmental haemostasis: Impact for clinical haemostasis laboratories, Thromb Haemost 2006; 95: 362-72
- Czech haematological Society CzMA JEP
- Clinical and Laboratory Standards Institute. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline. CLSI document C28-A3

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Erba Thrombin Reagent

For Fibrinogen Determination

Cat. No.	Pack name	Packaging (Content)
EHL00006	Erba Thrombin Reagent 20	4 × 5 mL

EN

INTENDED USE

Erba Thrombin Reagent is intended for use in the quantitative determination of fibrinogen in human plasma using the Clauss Method.¹

CLINICAL SIGNIFICANCE

Fibrinogen (also called Factor I) is a glycoprotein synthesised in the liver hepatocytes and megakaryocytes. Fibrinogen is composed of two sets of three different chains (α, β, and γ), linked to each other by disulfide bonds. On the alpha and beta chains, there is small peptide sequence (called a fibrinopeptide). Thrombin (Factor IIa) is responsible for converting fibrinogen into fibrin, by splitting out fibrinopeptides monomers. The monomers aggregate to form Fibrin, which is then cross linked by Factor XIIIa to form an insoluble fibrin clot.

Fibrinogen is an acute phase reactant protein. Therefore its blood level can be increased as a response of inflammatory states, with infections, during pregnancy, and after trauma.² A decrease of the fibrinogen level is found in congenital hypofibrinogenemia and afibrinogenemia. It can also be decreased in cases of Disseminated Intravascular Coagulation (DIC), systemic fibrinolysis, pancreatic and severe hepatic dysfunction.

PRINCIPLE

Clauss¹ developed a simple method for the quantitative determination of fibrinogen by measuring the clotting time of dilute plasma after the addition of Thrombin (>30 NIH units/mL). This clot time is proportional to the fibrinogen concentration. Determination of biological activity (i.e. the ability to form a coagulum) is more beneficial for diagnosis than determination of the amount of circulating fibrinogen.

COMPOSITION

Erba Thrombin Reagent; Vials of lyophilized bovine thrombin (100 NIH units/mL) with bovine albumin, buffer and preservative.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

Reagent contains <30 % bovine thrombin.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Supplemental information:

EUH208 Contains Thrombin. May produce an allergic reaction.

EUH210 Safety data sheet available on request.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

The Erba Thrombin Reagent contains bovine thrombin. Bovine thrombin contains bovine material. All donor animals were sourced from BSE-free herds. The cattle received ante- and post mortem health inspection by a veterinarian, and they were apparently free from infectious and contagious material. However, the material should be treated as potentially infectious.

WORKING REAGENT

Reconstitute each vial of Erba Thrombin Reagent with exactly 5 mL of purified water. Take care when pipetting to avoid contamination. Swirl gently and allow to stand for 15 minutes. Mix well immediately before use. Do not shake. Avoid reagent contamination

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

The reagent in the original closed vial after reconstitution is stable:

- 3 days on board ECL760 15 °C
- 1 week 2–25 °C in close bottle
- 1 month at -20 °C.

REAGENTS REQUIRED BUT NOT PROVIDE

Erba Standard Plasma (Cat. No.: EHL00012)

Erba Owrens Veronal Buffer (Cat. No.: EHL00021)

SAMPLE COLLECTION AND PREPARATION

It is recommended that specimen collection and storage be carried out in accordance with the CLSI guideline H21-A5.³ Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept between 15–25 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes.³

PROCEDURE

Manual method

1. Standard Curve Preparation

A new calibration should be performed whenever there is a change in the reagent batch or a change in the expected values of the control plasma tested. Prepare the following dilutions of Erba Standard Plasma (Cat.No.: EHL00012) in Erba Owrens Veronal Buffer (Cat. No.: EHL00021)

Tube	ECL 760 Recommended dilution	Dilution (Manual, ECL 105, ECL 412)	Erba Standard Plasma (mL)	Owren's buffer (mL)
1	2/1	1 + 4	0.2	0.8
2	1/1	1 + 9	0.1	0.9
3	1/2	1 + 19	0.1	1.9
4	1/3	1 + 29	0.1	2.9
5	1/4	1 + 39	0.1	3.9

2. Patient Sample Preparation:

- Prepare 1 + 9 dilutions of the patient plasma or control plasma in Owrens Buffer.
- Mix without shaking.

3. Testing

- Pipette 0.2 mL of standard, patient or control dilution into a reaction tube and incubate at 37 °C for 2 minutes.
- Add 0.1 mL of Erba Thrombin Reagent laboratory temperature and start simultaneously a timer.
- Note the time for clot formation to the nearest 0.1 second.
- Plot means Standard Clot Time (Y-axis) versus Fibrinogen Level (X-axis) on the fibrinogen graph paper. A straight line should be obtained. Assign the Erba Standard Plasma (Cat. No.: EHL00012) reference value to the 1 + 9 dilution to allow direct interpolation of patient and control values from the standard curve.

Automated method

Refer to the instrument's operator's manual.

REFERENCES VALUES

Reference limits^{4,5} for Fibrinogen (g/L) in plasma:

0–1 year:	1.50–3.40
1–6 years:	1.70–4.00
6–11 years:	1.55–4.00
11–16 years:	1.55–4.50
16–18 years:	1.60–4.20
≥18 years:	1.80–4.20

Reference intervals vary from laboratory to laboratory depending on the population, the technique and reagent lot. Therefore, each laboratory must establish its own reference intervals or verify them whenever one or more of the aforementioned variables are changed. For more information on establishing reference intervals see CLSI document C28-A3.⁶

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015 or Erba Control P Plus, Cat. No.: EHL00017).

LIMITATIONS

Fibrinolytic degradation product >100 mg/mL may cause falsely low fibrinogen quantitation. If values fall outside the standard curve values for the patient samples, re-assay using an appropriate dilution to bring values into the standard range.

No interference up to:

Heparin UFH: 3.3 IU/mL

Dabigatran: 475 µg/L

PERFORMANCES

These performances have been obtained using an analyzers ECL. Results may vary if a different instrument or manual procedure is used.

Precision

Intra assay	ECL 105/412		ECL 760	
	Clot formation (g/L)	CV (%)	Clot formation (g/L)	CV (%)
Erba Control N	2.56	4.78	3.13	5.19
Erba Control P	2.58	4.04	3.08	5.90
High Control	0.97	6.20	0.98	6.40

Inter assay	ECL 105/412		ECL 760	
	Clot formation (g/L)	CV (%)	Clot formation (g/L)	CV (%)
Erba Control N	2.52	5.11	3.51	10.39
Erba Control P	2.35	5.41	3.11	10.54

Limit of quantification	
ECL 105/412	ECL 760
0.287 g/L	0.406 g/L

Measuring range	
ECL 105	0.59–6.99 (g/L)
ECL 412	0.66–7.22 (g/L)
ECL 760	1.39–5.56 (g/L)

Interfering substance	ECL 105/412	ECL 760
Hemoglobin	min. 2 g/L	min. 3 g/L
Triglycerides	min. 460 mg/dL	min. 320 mg/dL
Bilirubin	min. 32 mg/dL	min. 24 mg/dL

REFERENCES

- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur bestimmung des Fibrinogens, Acta Haematol, 17:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, Crit.Rev.Clin.Lab.Sci, 8:145-192
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI:H21-A5
- Hematology of Infancy and Childhood, Nathan and Orkin, Saunders 2003, P. Monagle: Developmental haemostasis: Impact for clinical haemostasis laboratories, Thromb Haemost 2006; 95: 362-72
- Czech haematological Society CzMA JEP
- Clinical and Laboratory Standards Institute. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline. CLSI document C28-A3

USED SYMBOLS

LOT	Lot Number	IVD	In vitro diagnostic medical device		See Instructions for Use
REF	Catalogue Number		Manufacturer	CONT	Content
	Expiry Date		Storage Temperature		Biological Risks

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- Clauss A (1957) Gerinnungsphysiologische Schnell-methode zur-bestimmung des Fibrinogens, Acta Haema- tol, 17:237-246
- Shaw TS (1977) Assays for Fibrinogen and its Derivatives CRC, Crit.Rev.Clin.Lab.Sci, 8:145-192
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI:H21-A5
- Hematology of Infancy and Childhood, Nathan and Orkin, Saunders 2003, P. Monagle: Developmental haemostasis: Impact for clinical haemostasis laboratories, Thromb Haemost 2006; 95: 362-72
- Czech haematological Society CzMA JEP
- Clinical and Laboratory Standards Institute. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline. CLSI document C28-A3

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