
Erba Chrom Antithrombin III

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LITERATURA
1. Hirsh J et al (1989) Congenital antithrombin III deficiency. Incidence and clinical features, Am J. Med, 87(suppl. 3B): 34S-38S
2. Panicucci F et al (1981) Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition, Haemostasis, 9: 297-302
3. Conlan MG et al (1994) Antithrombin: antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors, Tromb. Haemost, 72(4): 551-556
4. Demers C et al (1993) An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition, Thromb. Haemost, 69:231-235
5. Tølefsen DM, Blank MK (1981) Detection of a new heparin-dependent inhibitor of thrombin in human plasma, J.Clin. Invest, 68:589-596
6. Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5

Číslo šarže

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Datum expirace

In vitro Diagnostikum

Výrobce

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Obsah

Erba Chrom Antithrombin III

Antithrombin III (Chromogenic)

Cat. No.	Pack name	Packaging (Content)
EHL00008	Erba Chrom Antithrombin III	R1: 3 x 2 ml AT III Factor Xa R2: 3 x 2 ml AT III Factor Xa Substrate R3: 3 x 3 ml AT III Diluent

INTENDED USE

Erba Chrom Antithrombin III is intended for the quantitative determination of antithrombin III (AT-III) activity in human citrated plasma by chromogenic assay.

CLINICAL SIGNIFICANCE

AT-III is a major inhibitor of blood coagulation, acting to inhibit plasma serine proteases including factors IXa, Xa, XIa, and thrombin. The rate of the inhibition is significantly increased in the presence of heparin. AT-III deficiency may be congenital or acquired and is associated with an increased risk for thrombosis.¹⁻³ Deficiencies may occur in liver disease disseminated intravascular coagulation (DIC), septicæmia, pulmonary embolism, nephritic syndrome, stroke, and thrombophlebitis.^{2,3}

PRINCIPLE

Most functional assays for AT-III rely on the ability of plasma to inactivate thrombin in the presence of heparin. More recently, it has been shown that a method based on the ability of plasma to inhibit factor Xa in the presence of heparin may be more accurate, since it eliminates interferences due to heparin cofactor II* which does not inhibit factor Xa.⁵

In the present method, factor Xa is added to a plasma dilution containing AT-III in the presence of excess heparin and calcium. After an initial incubation period, residual factor Xa is determined with a factor Xa-specific chromogenic substrate. The residual factor Xa activity is inversely proportional to the AT-III concentration.

COMPOSITION

R1 – AT III Factor Xa: Contains freeze-dried bovine factor Xa.

R2 – AT III Factor Xa Substrate: Contains freeze-dried CH₃OCO-D-CHA-Gly-Arg-pNA-AcOH.

R3 – AT III Diluent: Contains 5x buffer concentrate with sodium azide at <1 % as preservative. When fully diluted, buffer contains 0.05 M Tris-HCl, 0.175 M NaCl, 7.5 mM Na₂EDTA and sodium heparin at pH 8.4.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in this kit contain materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen. and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

WORKING REAGENT.

R1 – AT III Factor Xa: Reconstitute with 2 ml of distilled/deionised water. Allow to stand for 10 minutes and mix before use.

R2 – AT III Factor Xa Substrate: Reconstitute with 2 ml of distilled/deionised water. Allow to stand for 10 minutes and mix before use.

R3 – AT III Diluent: Dilute 1 volume of diluent for 4 volume of distilled/deionised water.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Reconstituted reagents are stable for

R1 – AT III Factor Xa and **R2** – AT III Factor Xa Substrate:

- 2 months at 2–8 °C.
- 1 month at 17 °C
- 2 days at 37 °C.

R3 – AT III Diluent: Store diluted buffer in a tightly stoppered bottle at 2–8 °C and use in one month.

REAGENTS REQUIRED BUT NOT PROVIDED

Erba Standard Plasma (Cat. No.: EHL00012)

Glacial acetic acid

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2 % or 3.8 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at 2–8 °C or 18–24 °C. Testing should be completed within 4 hours of sample collection, or plasma can be storage frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes.⁶

PROCEDURE

Manual Method

Preparation of plasma standard (Erba Standard Plasma), Control and Patient Dilutions:

% AT III	Plasma	Dilution Buffer
100 %	10 µl Standard	990 µl
50 %	500 µl of 100 % Std	500 µl
25 %	500 µl of 50 % Std	500 µl
12.5 %	500 µl of 25 % Std	500 µl
Patient or control	10 µl Plasma	990 µl

A. Semi-Micro Endpoint Method

- Add 200 µl standard / control or patient plasma dilution to a test tube.
- Incubate at 37 °C for 2–4 minutes.
- Add 200 µl Factor Xa Reagent and mix.
- Incubate at 37 °C for EXACTLY 1 minute.
- Add 200 µl Factor Xa Substrate and mix.
- Incubate at 37 °C for EXACTLY 3 minutes.
- Add 200 µl acetic acid (50%) and mix.
- * Add 200 µl water (optional).

The yellow colour of the final reaction product is stable for at least 4 hours. Read absorbance at 405 nm in a 1 cm semi-micro cuvette against a blank prepared in the following order:

- 200 µl acetic acid.
- 200 µl standard dilution.
- 200 µl Factor Xa reagent.
- 200 µl Factor Xa Substrate.
- *200 µl water (optional).

*Some spectrophotometers require a minimum of 1 ml volume in the cuvette.

If the patient plasma is very icteric, a second blank containing patient plasma dilution instead of standard dilution should be prepared and the absorbance subtracted from the absorbance obtained for the patient's AT-III determination

B. Kinetic Method

A kinetic analyser may be used to measure the initial rate of hydrolysis of the chromogenic substrate. The procedure to be used is as follows:

To a reaction cuvette:

- Add 200 µl diluted standard or patient plasma.
- Incubate at 37 °C for 1 minutes.
- Add 200 µl Factor Xa reagent and mix.
- Incubate at 37 °C for 1 minute.
- Add 200 µl Factor Xa Substrate.
- Measure rate of change of absorbance at 405 nm for 1 minute.

Plot the absorbance obtained with each of the AT-III calibration standards on the y-axis against % AT-III on the x-axis using linear graph paper. The line of best fit should be determined by linear regression analysis. The AT-III in plasma samples can be determined by interpolation from the calibration curve. The AT-III concentration in the patient specimen should be adjusted for the AT-III concentration in the standard:

$$\% \text{ AT-III (adjusted) } = \% \text{ AT-III (patient) } \times \% \text{ AT-III (standard) } / 100$$

Automated Method

Refer to the instrument's operator's manual.

REFERENCES VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. The normal range of AT-III is 75–125% in plasma. Plasma activity levels of 30–60% can be observed in patients with hereditary AT-III deficiency.^{1,2} Several clinical conditions associated with acquired AT-III deficiency include: liver disease, DIC, nephrotic syndrome, pulmonary embolism, stroke, and thrombophlebitis. In addition, oral contraceptive use may lead to reduced AT-III levels.²

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P Plus, Cat. No.: EHL00017)

LIMITATIONS

Avoid icteric, lipaemic, and haemolysed samples.

PERFORMANCES

These performances have been obtained using an analyzer ECL. Results may vary if a different instrument or manual procedure is used.

Linearity: The Erba Chrom Antithrombin III reagent is designed to give a linear calibration in the range of 10–140 %.

Precision:

	Intra-assay precision (N = 10)		Inter-assay precision (N = 30)	
% AT III	98.7	55.6	98.4	56.4
CV (%)	0.52	1.38	0.44	2.79

REFERENCES

- Hirsh J et al (1989) Congénital antithrombin III deficiency. Incidence and clinical features, Am J. Med, 87(suppl. 3B): 34S-38S
- Panicucci F et al (1981) Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition, Haemostasis, 9: 297-302
- Conlan MG et al (1994) Antithrombin: antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors, Tromb. Haemost, 72(4): 551-556
- Demers C et al (1993) An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition, Thromb. Haemost, 69:231-235
- Tølefsen DM, Blank MK (1981) Detection of a new heparin-dependent inhibitor of thrombin in human plasma, J.Clin. Invest, 68:589-596
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5

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- Hirsh J et al (1989) Congenital antithrombin III deficiency. Incidence and clinical features, Am J. Med, 87(suppl. 3B): 34S-38S
- Panicucci F et al (1981) Antithrombin III, heparin cofactor and antifactor Xa in relation to age, sex and pathological condition, Haemostasis, 9: 297-302
- Conlan MG et al (1994) Antithrombin: antithrombin III: Associations with age, race, sex and cardiovascular disease risk factors, Tromb. Haemost. 72(4): 551-556
- Demers C et al (1993) An antithrombin III assay based on factor Xa inhibition provides a more reliable test to identify congenital antithrombin III deficiency than an assay based on thrombin inhibition, Tromb. Haemost. 69:231-235
- Tofelsen DM, Blank MK (1981) Detection of a new heparin-dependent inhibitor of thrombin in human plasma, J.Clin. Invest. 68:589-596
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Haemostasis Assays: Approved Guideline, 5th edn. CLSI: H21-A5

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