Erba Chrom Protein C

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LITERATURA

- Walker FJ (1981) Regulation of activated protein C by protein S. The role of phospholipid in factor Va inactivation, J Biol Chem, 256: 11128-31
- Taylor BF, Lockhart MS (1985) A new function for activated protein C: activated protein C prevents inhibition of plasminogen activators by release of mononuclear leukocytes-platelet suspensions stimulated by phorbol diester, Thromb Tes, 37(1): 155-64
- Fulcher CA, Gardiner JE, Griffin JH, Zimmerman TS (1984) Proteolytic inactivation of human factor VIII procoagulant protein by activated protein C and its analogy with factor V, 63(2): 486-9
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of blood specimens for testing plasma-based coagulation assays and molecular haemostasis assays: approved guideline, 5th edn. CLSI: H21-A5
- Young DS et al. Effects of drugs on clinical laboratory tests, 3rd ed., AACC Press Washington, D.C., 1990
- Bertina RM, Broekmans AW, Krommenhoek-Van Es C, Van Wijngaarden A (1984) The use of a functional and immunologic assay for plasma protein C in the study of the heterogeneity of congenital protein C deficiency, Thromb Haemostas, 51:1-5
- Pabinger-Fasching I, Bertina RM, Lechner K, Niesser H, Korninger C (1983) Protein C deficiency in two Austrian families, Thromb haemostas, 50:810-3

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INTENDED USE

Erba Chrom Protein C assay is intended for the quantitative determination of protein C in human plasma using a chromogenic assay method.

CLINICAL SIGNIFICANCE

Protein C is a vitamin K dependent protein which plays an important role in the regulation of anticoagulant mechanisms. It can inhibit coagulation by inactivating factors Va¹ and VIIIa² or, when activated, can stimulate fibrinolysis.³ Protein C circulates as a zymogen, and is converted to an active serine protease by the action of thrombin in the presence of thrombomodulin. Both hereditary and acquired Protein C deficiencies have been shown to be a risk factor for development of venous thrombosis.

PRINCIPLE

Protein C in plasma is activated by a specific fraction from the *Agkistrodon contortrix* snake venom. The amount of activated protein C (APC) is determined by monitoring the rate of hydrolysis of a protein C specific chromogenic substrate. The release of pNA is measured at 405 nm and is proportional to the protein C level.

COMPOSITION

R1 – Protein C substrate: 2,75 µmol lyophilised pyro-Glu-Pro-Arg-pNA.HCl.

R2 – Protein C activator: 0,8 units of activator from snake venom (Protac[®]).

R3 – Protein C diluent: tris buffer with sodium azide as a preservative.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

Reagent 1 are not classified as dangerous.

Reagent 2: EUH 208 Contains Protac[®]. May produce an allergic reaction.

EUH 210 Safety data sheet available on request.

Reagent 3 diluent contains < 1 % imidazole.

Danger

Hazard statements:

H360D May damage the unborn child.

EUH032 Contact with acids liberates very toxic gas.

Precautionary statements:

P308+P313 IF exposed or concerned: Get medical advice/attention.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

The Protein C Activator reagent contains snake venom. Therefore, users of reagent of this type must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of this biological material as if it was infectious.

WORKING REAGENT

R1 – Protein C substrate: reconstitute each vial with 2 ml of Protein C Diluent. If cloudy, warm at 37°C for a few minutes.

R2 – Protein C activator: reconstitute each vial with 2 ml deionised water.

R3 – Protein C diluents: ready for use.

Avoid reagents contamination.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

Reconstituted reagents **R1, R2** are stable:

- 1 week at 2–8°C,
- 1 month at -20°C.

REAGENTS REQUIRED BUT NOT PROVIDED

Erba Standard Plasma (Cat. No.: EHL00012)

NaCl solution (0.9 %)

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2 % or 3.8 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 x g for 15 minutes. Plasma should be kept at 2–8°C or at 18–24°C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for 6 months. Thaw quickly at 37°C prior to testing. Do not keep at 37°C more than 5 minutes.⁴

Erroneous results may be caused by contamination with tissue fluids or stasis. Avoid agitation, air bubbles or foaming. For the effects of commonly administered drugs, refer to Young, et al.⁵

PROCEDURE

Manual method

Calibration

Prepare plasma standards and patients plasma samples as follows: important: use only 0.9 % NaCl solutions for dilutions.

Standard %	Plasma	Buffer
100 %	100 µl Erba Standard Plasma	+ 300 µl saline
50 %	50 µl Erba Standard Plasma	+ 350 µl saline
0 %	—	400 µl saline only
Patient	100 µl plasma	+ 300 µl saline

End Point Method Assay

To a glass or plastic test tube:

- Add 100 µl standard or patients' plasma dilution.
- Incubate at 37°C for 2 minutes.
- Add 200 µl Protein C Activator and mix.
- Incubate at 37°C for 5 minutes.
- Add 200 µl Protein C Substrate and mix.
- Incubate at 37°C for 10 minutes
- Add 200 µl acetic acid and mix.
- * Add 200 µl water (optional).

* Some spectrophotometers require a minimum of 1 ml volume in the cuvette.

Read absorbance at 405 nm in a 1 cm semi-micro cuvette against a blank prepared with deionised water. As many as ten determinations can be performed simultaneously with the same stop watch by staggering pipetting steps at five second intervals.

Automated method

Refer to the instrument's operator's manual.

REFERENCES VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. Protein C values are usually expressed in relative percentage compared to a pooled normal plasma standard. Bertina et al⁶ reported a range of 65–145 % in healthy individuals, with diminished levels found following anticoagulant therapy. There is apparently no difference in protein C between healthy males or females.

QUALITY CONTROL

Each laboratory should establish a quality control program.

Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

It is suggested to use two levels of control, one close to the normal patient values (Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P Plus, Cat. No.: EHL00017).

LIMITATIONS

Icteric, haemolysed or hyperlipemic specimens interfere with absorbance readings, thus requiring the use of plasma blanks for accurate results. Plasma blanks are also needed on patients when contact factor activation is suspended, such as DIC patients, or from individuals on oral contraceptives where cold activation may occur. Blanks are obtained by substituting saline for the Protein C Activator in the test reaction. Subtract the sample blank from the sample test activity.

The naturally-occurring activator of protein C is thrombin in the presence of thrombomodulin. The possible existence of clinical protein C abnormalities not detectable by reactivity with the snake venom activator or with the chromogenic substrate used in this kit should be not precluded. To assure accurate, reproducible results, use accurate pipetting devices and observe recommended procedures with particular emphasis on incubation times and incubation temperature.

INTERFERENCES

Haemoglobin: Negative bias from 2 g/l on normal plasma (106.3 %).

Total bilirubin: Negative bias from 140 mg/l (239.4 µmol/l) on normal plasma (102 %).

Turbidity: Negative bias from 3 g/l (3.42 mmol/l) equivalent triglycerides on normal plasma (106.3 %).

PERFORMANCES

These performances have been obtained using an analyzer ECL. Results may vary if a different instrument or manual procedure is used.

Limit of Detection

The limit of detection for the Erba Chrom Protein C assay was determined at a concentration of 5.3 %.

Linearity

The Erba Chrom Protein C assay is designed to give a linear calibration in the range of 14–110 %.

Precision:

	Intra-assay precision (N = 20)		Inter-assay precision (N = 20)	
Mean (s)	110.9	64.0	113.6	55.2
CV (%)	2.52	2.67	3.60	4.89

REFERENCES

- Walker FJ (1981) Regulation of activated protein C by protein S. The role of phospholipid in factor Va inactivation, J Biol Chem, 256: 11128-31
- Taylor BF, Lockhart MS (1985) A new function for activated protein C: activated protein C prevents inhibition of plasminogen activators by releasate from mononuclear leukocytes-platelet suspensions stimulated by phorbol diester, Thromb Tes, 37(1): 155-64
- Fulcher CA, Gardiner JE, Griffin JH, Zimmerman TS (1984) Proteolytic inactivation of human factor VIII procoagulant protein by activated protein C and its analogy with factor V, 63(2): 486-9
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of blood specimens for testing plasma-based coagulation assays and molecular haemostasis assays: approved guideline, 5th edn. CLSI: H21-A5
- Young DS et al. Effects of drugs on clinical laboratory tests, 3rd ed., AAC Press Washington, D.C., 1990
- Bertina RM, Broekmans AW, Krommenhoek-Van Es C, Van Wijngaarden A (1984) The use of a functional and immunologic assay for plasma protein C in the study of the heterogeneity of congenital protein C deficiency, Thromb Haemostas, 51:1-5
- Pabinger-Fasching I, Bertina RM, Lechner K, Niesser H, Korninger C (1983) Protein C deficiency in two Austrian families, Thromb haemostas, 50:810-3

USED SYMBOLS

LOT	Lot Number	IVD	In vitro Diagnostics	i	See Instruction for Use
REF	Catalogue Number	M	Manufacturer	CONT	Content
E	Expiry Date	T	Storage Temperature		

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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- Walker FJ (1981) Regulation of activated protein C by protein S. The role of phospholipid in factor Va inactivation, J Biol Chem, 256: 11128-31
- Taylor BF, Lockhart MS (1985) A new function for activated protein C: activated protein C prevents inhibition of plasminogen activators by release of factor Va from mononuclear leukocytes-platelet suspensions stimulated by phorbol diester, Thromb Tes, 37(1): 155-64
- Fulcher CA, Gardiner JE, Griffin JH, Zimmerman TS (1984) Proteolytic inactivation of human factor VIII procoagulant protein by activated protein C and its analogy with factor V, 63(2): 486-9
- Clinical and Laboratory Standards Institute (2008) Collection, Transport and Processing of blood specimens for testing plasma-based coagulation assays and molecular haemostasis assays: approved guideline, 5th edn. CLSI: H21-A5
- Young DS et al. Effects of drugs on clinical laboratory tests, 3rd ed., AACC Press Washington, D. C., 1990
- Bertina RM, Broekmans AW, Krommenhoek-Van Es C, Van Wijngaarden A (1984) The use of a functional and immunologic assay for plasma protein C in the study of the heterogeneity of congenital protein C deficiency, Thromb Haemostas, 51:1-5
- Pabinger-Fasching I, Bertina RM, Lechner K, Niesser H, Korninger C (1983) Protein C deficiency in two Austrian families, Thromb haemostas, 50:810-3

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