Erba DDimer Calibrator D-dimer Calibrator

Cat. No.:	Pack name:	Packaging (Content):
EHL00018	Erba DDimer Calibrator	R1: 2 x 1 ml DDimer Calibrator R2: 2 x 7 ml DDimer Diluent
EHL00044	Erba DDimer Calibrator	R1: 1 x 1 ml DDimer Calibrator R2: 1 x 7 ml DDimer Diluent

INTENDED USE

Erba DDimer Calibrator may be used to establish calibration curve in the quantitative determination of D-dimer.

COMPOSITION

R1 – DDimer Calibrator: 1ml of lyophilized human plasma with fibrin D-dimer added.

R2 – DDimer Diluent: containing buffer and preservatives.

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1:
The reagent is not classified as dangerous.

R2:
EUH208 Contains 2-methylisothiazol-3(2H)-one. May produce an allergic reaction.

EUH210 Safety data sheet available on request.

WARNING—POTENTIAL BIOHAZARDOUS MATERIAL

The calibrator provided in this kit contains materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

WORKING REAGENT

R1 – DDimer Calibrator: Reconstitute each vial of DDimer Calibrator with exactly 1 ml of distilled water / saline solution. Swirl gently. Allow to stand for 15 minutes for complete dissolution and mix before use.

R2 – DDimer Diluent: Ready to use.

Avoid calibrator and diluent contamination.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Once reconstituted, the calibrator remains stable for:

- 12 hours at 4–25 °C.

ASSAY VALUE

LOT	XXXXXXX
EXP	MM/RRRR

Product/Test	Unit	Value
DDU*	ng/ml	0000
FEU**	ng/ml	0000

* D-dimer Units (DDU) are assigned against D-dimer in-house reference material. Units of ng/ml = $\mu\text{g/l} = 10^{-3} \mu\text{g/ml}$

** Fibrinogen Equivalent Units (FEU) are calculated from DDU by multiplication with the correlation factor 2.5. Units of ng/ml = $\mu\text{g/l} = 10^{-3} \mu\text{g/ml}$

PROCEDURE

Each calibrator should be treated in the same manner as the patients' sample in accordance with the instructions outlined in each particular assay procedure.

Standard Curve Preparation

A new standard curve should be used when reagent lot numbers change or if expected values in quality control testing change. Prepare the following dilutions of DDimer Calibrator in DDimer Diluent:

Tube	Dilution	DDimer Calibrator (ml)	DDimer Diluent (ml)
1	1 / 1	0.4	0
2	1 / 2	0.2 from tube 1	0.2
3	1 / 4	0.2 from tube 2	0.2
4	1 / 8	0.2 from tube 3	0.2
5	1 / 16	0.2 from tube 4	0.2
6	1 / 32	0.2 from tube 5	0.2

EXPECTED VALUES

The calibrator values for each of the parameters may vary from one lot to another. Refer to assay sheet section of this instruction.

QUALITY CONTROL

To ensure adequate quality, control plasmas are recommended. It is suggested to use two levels of control, one close to the normal patient values and the second representative the pathologic values. Erba DDimer Control kit (Cat. No.: EHL00019) are available for use with this product.

LIMITATIONS

The results obtained with Erba DDimer Calibrator depend on several factors strongly associated with instrumentation, types of reagents, deficient substrates and laboratory to laboratory variations.^{1, 2} Each laboratory should establish an expected range for the particular instrument-reagent system.

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