## Erba Owrens Veronal Buffer

Cat. No.	Product Name	Packaging (Content)
EHL00021	Erba Owrens Veronal Buffer	6 × 25 mL, instruction for use

### INTENDED USE

Erba Owrens Veronal Buffer is a universal buffer used for dilution of reference plasma and human blood plasma samples in quantitative *in vitro* determination of hemostasis parameters. It is intended for semi-automated or automated method. For professional use in clinical laboratory only.

### PRINCIPLE

The diluent is a solution with certain ion strength and conductivity, which can dilute calibrator or human plasma and provide them stable environment.

### DESCRIPTION AND COMPOSITION

Erba Owrens Veronal Buffer: aqueous solution of 28.4 mM barbitol, 125.4 mM sodium chloride, 0.05 % sodium azide as preservative; pH 7.2–7.6.

### PREPARATION

Ready to use. Avoid contamination.

### STABILITY AND STORAGE

The unopened buffer is stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

### SPECIMEN COLLECTION AND HANDLING

It is recommended that specimen collection and storage should be carried out in accordance with the CLSI guideline H21-A5. Plastic or siliconized glass should be used throughout. Blood (9 parts) should be collected into 3.2 % sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 1500 g for 15 minutes. Plasma should be kept between 15–25 °C. Testing should be completed within 4 hours of sample collection, or plasma can be stored frozen at -20 °C for 2 weeks or -70 °C for 6 months. Thaw quickly at 37 °C prior to testing. Do not keep at 37 °C for more than 5 minutes.

### ASSAY PROCEDURE

The Erba Owrens Veronal Buffer should be used in accordance with the instructions for use of the specific reagent kit or user manual of a specific analyser.

### WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

### Hazards identification in accordance with Regulation (EC) No 1272/2008

The reagent is not classified as dangerous.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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