# Erba LA Control High

## Control for the determination of Lupus anticoagulants (LA)

Cat. No.:	Pack name:	Packaging (Content):
EHL00039	Erba LA Control High	6 x 1 ml

EN

### INTENDED USE

Erba LA Control High is prepared with plasma from human donors positive for Lupus anticoagulants. The plasma gives results typical of a Lupus Anticoagulant patient in LA1 Screen, LA2 Confirm and APTT-based tests.

### COMPOSITION

Erba LA Control High: lyophilised, citrated human plasma

### WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and

### WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in these kits contain materials of human and/or animal origin. Whenever human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV-1, HIV-2 and HCV, and for hepatitis B surface antigen. and results are found to be negative. However, no test method can offer complete assurance that infectious agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

### WORKING REAGENT

Reconstitute each vial of Erba LA Control High with exactly 1 ml of distilled/deionised water. Swirl gently and allow to stand for 15 minutes. Mix well before use. Do not shake.

### STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

Reconstituted vials of plasma should be kept on ice and are stable for:

- 4 hours at 2–8 °C.

Keep covered.

### PROCEDURE

Each control should be treated in the same manner as the patients' sample in accordance with the instructions outlined in each particular assay procedure.

### EXPECTED VALUES

The control values for each of the parameters may vary from one lot to another. Refer to the Assay Value insert provided in the box. For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically re-evaluate the mean.

### QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

### LIMITATIONS

The results obtained with Erba LA Control High depends on several factors strongly associated with instrumentation, types of reagents, and laboratory to laboratory variations.<sup>1, 2, 3</sup> Each laboratory should establish an expected range for the particular instrument-reagent system

### REFERENCES

- Kirkwood TBL et al. (1977) Identification of Sources of Variation in Factor VIII Assay, British Journal of Haematology, 37:559-568.
- Goldenfarb MD (1971) Reproducibility in Coagulation Assays, AJCP, 55:561-564.
- Palkuti HA and Longberry JR (1973) A Precision Study of Coagulation Factor Assay Techniques, AJCP, 59:231-235.

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