Erba Protime R

Protrombinový čas

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LITERATURA


- CLSI: Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline- Fifth Edition. CLSI document: H21-A5; 28:5; 2008.
- CLSI: One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline-Second Edition. CLSI document: H47-A2; 28:20; 2008.
- De Caterina R et al: Vitamin K antagonists in heart disease: Current status and perspectives (Section III). Thromb Haemost; 110: 1087-1107; 2013.

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INTENDED USE

Erba Prottime R is a recombinant human thromboplastin reagent, produced by genetic technology in *Escherichia Coli*, with own solvent intended for the determination of the prothrombin time (PT).

CLINICAL SIGNIFICANCE

Erba Prottime R reagent is a recombinant human thromboplastin, which contains recombinant human tissue factor, lipids and calcium ions. The first standardized one-stage prothrombin time test was developed by Dr. Armand Quick in 1935. It has now become the basic coagulation screening test for the diagnosis of congenital and acquired deficiencies of clotting factors from the extrinsic pathway (factors II, V, VII and X). It is also used for the induction and monitoring of oral anticoagulant therapy and can be used to assess the protein synthesis capability of the liver in chronic or acute hepatic disorders.

Therefore, the PT by Erba Prottime R is optimally used for presurgical screening and monitoring for oral anticoagulant therapy (OAT) also. Erba Prottime R reagent with the corresponding deficient plasmas is also suitable for determination of activity of extrinsic coagulation pathway. Furthermore the Erba Prottime R reagent has an increased sensitivity for certain factors, like Factor VII, due to this characteristics in some cases it might cause prolonged coagulation compared with tissue extract thromboplastin.

PRINCIPLE

The one-stage PT measures the clotting time of plasma after adding a source of tissue factor (thromboplastin) and calcium. The recalcification of plasma in the presence of tissue factor generates activated factor Xa. Factor Xa in turn activates prothrombin to thrombin, which converts fibrinogen to an insoluble fibrin clot.

The time of this clotting process is measurable manually or with optical or mechanical coagulation analysers.

COMPOSITION

Erba Prottime R: is a freeze-dried, recombinant human thromboplastin from *Escherichia Coli* with lipids and stabilizers.

Solvent is a buffer, which contains calcium ions and sodium azide (< 0,01 %) as preservative

WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only. These reagents are to be used by certified medical laboratory personnel only.
- Do not ingest.
- Wear gloves when handling all kit components.
- Only use clean or single use laboratory equipment to avoid contaminations.
- The eventual rest of reagents should be disposed of in accordance with the internal regulations and in compliance with local and national regulations relating to the safe handling of waste.
- Reagents containing sodium azide should be discarded with care to prevent the formation of explosive metallic azides. If waste materials are dumped into sinks, use copious quantities of water to flush plumbing thoroughly.

WARNING – POTENTIAL BIOHAZARDOUS MATERIAL

Some reagents provided in these kits contain materials of human and/or animal origin. *Never human plasma is required for the preparation of these reagents, the plasmas are tested for the antibodies to HIV 1, HIV 2 and HCV, and for hepatitis B surface antigen, and results are found to be negative. However, no test method can offer complete assurance that infectious*

agents are absent. Therefore, users of reagents of these types must exercise extreme care in full compliance with regulatory safety precautions in the manipulation of these biological materials as if they were infectious.

WORKING REAGENT

Erba Prottime R reagent is dissolved with the entire contents of one vial Solvent of the same lot. Keep the reagent at room temperature (20–25°C) for at least 30 minutes for proper reconstitution. Swirl the vial gently, horizontally more times (5–10) before using it, but do not shake. Wait until the reagent reaches the working temperature!

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

Once reconstituted, the reagent remains stable for:

- 8 hours at 37°C
- 24 hours at 15–25°C
- 10 days at 2–8°C.

Do not freeze!

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout.

Blood (9 parts) should be collected into 3.2 % sodium citrate anticoagulant (1 part). Avoid hemolysis and contamination by tissue fluids. Samples that have less than 90 % of the expected fill volume should be rejected. Centrifuge blood for 15 minutes at 1500 x g. Test within 2 hours if samples are held at 22–24°C. If testing is not completed within 24 hours, plasma should be frozen at -20°C for up to two weeks or -70°C for up to 6 months.

PROCEDURE

Erba Prottime R is a one-stage PT test which can be used with semi-automated coagulation analysers according to the protocol detailed below.

Manual Method

- Incubate reagent at 37°C for 15 minutes.
- Add 50 µl of the patient plasma or control plasma into a reaction tube and incubate at 37°C for 2 minutes.
- Add 100 µl of reagent and start simultaneously a timer.
- Note the time for clot formation to the nearest 0.1 seconds

Automated method

Refer to the instrument's operator's manual.

INTERPRETATION OF RESULTS

Erba Prottime R test results can be reported in the following units, lot specific sheet in the box will help in the calculation:

- Seconds, which means the observed clotting time.
- Percentage, which means the proportional part of the normal PT activity, which is calculable from the calibration curve. Method dependent master curve in the issued sheet can be used for the calculation.
- International Normalized Ratio (INR), calculated as $INR = (PT/MNPT)^{ISI}$. Where PT means the clotting time of the sample, MNPT means Mean Normal Prothrombin Time and ISI means International Sensitivity Index. The ISI method dependent values in the lot specific values sheet can be used for the calculation.

Notes:

- The INR is the only officially recognized dimension of the result at vitamin K antagonists treated patients. The normal range expressed in INR is 0,8–1,2.
- Every laboratory should determine its own MNPT value and reference range. MNPT value in the issued sheet is only for information, because it depends on the measuring circumstances and population.
- Accurate and general conversion of percentage into INR (or back) is not possible!

REFERENCES VALUES

These values are intended as a guideline only. Each laboratory should establish a Normal Reference Range (NRR) using instrumentation, blood collection methods, and testing techniques used in that laboratory. The NRR should be reestablished or at least verified when changing lot numbers of the same reagent. 3,6 A new NRR should be established with any change in reagents, instrumentation, blood collection techniques, or anticoagulant.

The clotting time of abnormal plasmas will depend on the ISI of the reagent lot in use.

QUALITY CONTROL

Each laboratory should establish a quality control program.

To ensure adequate quality, controls plasma are recommended. It is suggested to use two levels of control, one close to the normal patient values (Erba Control N, Cat. No.: EHL00014 or Erba Control N Plus, Cat. No.: EHL00016) and the second representative the pathologic values (Erba Control P, Cat. No.: EHL00015 or Erba Control P Plus, Cat. No.: EHL00017).

LIMITATIONS

The result of PT test with Erba Prottime R reagent may be influenced by drugs and other pre-analytical interfering agents. The potential limits of these parameters were tested on Diagon analysers (Coag Line) with the following result:

Heparin – 1.25 IU/ml
Hemoglobin – 6.8 g/l
Triglycerides – 9 mmol/l
Bilirubin – 270 µmol/l

PERFORMANCES

These performances have been obtained using an optical coagulation analyzer. Results may vary if a different instrument or manual procedure is used.

Precision:

	Intra-assay precision (N = 10)		Inter-assay precision (N = 10)	
Mean (s)	11.3	16.2	11.8	15.5
CV (%)	0.82	0.49	2.54	1.02

REFERENCE

- CLSI: Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline- Fifth Edition. CLSI document: H21-A5; 28:5; 2008.
- CLSI: One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline-Second Edition. CLSI document: H47-A2; 28:20; 2008.
- De Caterina R et al: Vitamin K antagonists in heart disease: Current status and perspectives (Section III). *Thromb Haemost*; 110: 1087-1107; 2013.

USED SYMBOLS

LOT	Lot Number	IVD	In vitro Diagnostics	i	See Instruction for Use
REF	Catalogue Number	Manufacturer	Manufacturer	CONT	Content
Expiry Date	Expiry Date	Storage Temperature	Storage Temperature		

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com