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Dátum revízie: 20. 12. 2022

**Utilization:**  
The diagnostic strip MicroalbuPHAN® LAURA is intended for semi-quantitative determination of albumin and creatinine in urine. The diagnostic strip MicroalbuPHAN® LAURA can be used for screening of the migroalbuminuria; positive results should be confirmed by quantitative methods.  
The diagnostic strip MicroalbuPHAN® LAURA is intended for *in vitro* diagnostics for professional use only.  
The clinical importance of the positive results of low concentration of albumin in urine is the detection of risk of the beginning kidney damage. The microalbuminuria is the first indicator of the possibility of the incipieny of nephropathy; of other diabetes-related complications. It is a strong predictor of cardiovascular disease and risk factor in essential hypertension and also early marker of complications in diabetic pregnancy.

**Instructions:**  
For objective evaluation of test strips, please, use the user manual for urine analyser LAURA® Smart.  
Use only a freshly voided, well mixed, non centrifuged specimen without preservatives collected in a clear container without traces of detergents and disinfectants. The urine should not be more 4 hours old when tested. Do not divide test strips!  
1. Remove only as many test strip required and reseal the tube immediately after use.  
2. Do not touch test pads of strip.  
3. Immerse completely all reagent pads in specimen (for 1–2 seconds)  
4. Run edge of the strip against rim of urine container to remove excess urine and wipe off the excess urine from the strip to the cellulose cotton wool as well.  
Hold strip in horizontal position.  
5. Evaluate the result using the reader LAURA® Smart, follow enclosed instructions for that instrument.  
6. For the correct interpretation of the results take into account the ratio of albumin / creatinine.

**Notes:** For visual evaluation compare the tests pads to the corresponding colour scale on the label after approx. 60 sec. According to the different spectral sensitivity of the human eye and the optical system it is not possible to guarantee the exact correspondence between visual reading and reading by analyser.

**Calibration:**  
For more accurate results we recommend to calibrate for both diagnostic pads. Six-digit calibration code which is entered into the urine analyser LAURA® Smart is noted on the label always under the batch number.  
A detailed calibration procedure is mentioned in user manual for urine analyser LAURA® Smart.

**Working reagent concentrations:**  
**Creatinine:** 3,5-dinitrobenzoic acid 4.4%, **Albumin:** tetrabromphenol blue 0.17%  
**Principle of test:**  
**Creatinine** – The test is based on reaction of creatinine with 3,5-dinitrobenzoic acid in alkaline medium (Benedict-Behre reaction).  
**Albumin** – The test is based on colour change of acid basic indicator, which is caused by presence of proteins (so-called proteins error of indicator). It is particularly sensitive to albumin, but is much less sensitive to other urinary proteins.

**Limitations:**  
**Creatinine** – In urines with high buffering capacity false negative reading may be obtained. High concentrations of ketone bodies (> 50 mmol/l) false positive reading may be occurred. In case a urine sample is more markedly coloured (i.g. increased bilirubin content), the resulting colour could be affected by a sample coloration.  
**Albumin** – In strongly alkaline urines (pH > 8) with high buffering capacity, high concentrations of creatinine (> 26,5 mmol/l), urines with traces of quinine or quinoline contains drugs false positive reading may be obtained. False positive results may be found when the urine collection vessel contains traces of disinfectants with quarternary ammonium groups. In the presence anionic or non-ionic detergents false negative results may occur. Blood with higher concentration than 2000 Ery/l may cause false positive results.

The determination of both parameters could be limited by extreme dark colours of urine samples which are caused by high concentration of bilirubin, urobilinogen.

**Interpretation of results:**  
**Albumin** – Albumin is normally present in urine at concentrations of less than 0.02 g/l (2 mg/dl). Microalbuminuria is indicated with results of 0.02–0.3 g/l (2–30 mg/dl), high concentrations are indicated clinical albuminuria (proteinuria).  
**Creatinine** – Creatinine is normally present in urine at concentrations 2.2–26.5 mmol/l (0.25–3.0 g/l).  
**Albumin-to-Creatinine Ratio** – Albumin is normally present in urine at concentrations of less than 3.4 mg albumin / mmol creatinine (0.03 g albumin/g creatinine). Microalbuminuria is indicated at ratio results of 3.4–33.9 mg/mmol (30–300 mg/g), high values of ratio are indicated clinical albuminuria (proteinuria).

Albumin / Creatinine (mg/mmol)	=	$\frac{\text{Albumin (g/l)} \times 1000}{\text{Creatinine (mmol/l)}}$	=	$\frac{\text{Albumin (mg/l)}}{\text{Creatinine (mmol/l)}}$
Albumin / Creatinine (mg/g)	=	$\frac{\text{Albumin (g/l)} \times 1000}{\text{Creatinine (g/l)}}$	=	$\frac{\text{Albumin (mg/l)}}{\text{Creatinine (g/l)}}$

<sup>1</sup> Burtis, C.A., Ashwood E.R.: Tietz Textbook of Clinical Chemistry 3<sup>rd</sup> Edition, W.B. Saunders Comp., Philadelphia (1999)

<sup>2</sup> Needleman, S.B., Porvaznik, M., Ander, D.: J. Forensic Sci. 37, 1125 (1992)

<sup>3</sup> Sacks, D.B., Bruns, D.E., Goldstein, D.E., Maclaren, N.K., McDonald, J.M., Parrott, M: Clin. Chem. 48, 436 (2002)

**Albumin / Creatinine ratio is calculated by urine analyser LAURA® Smart automatically for each measured urine sample.**  
**Please note:** Knowledge of the effects drugs or their metabolites upon the individual tests is not yet complete. In doubtful cases, it is advisable to repeat the test after discontinuing a drug. The sensitivity depends upon variability of urines. The semi-quantitative analysis isn't sufficient for the completing of diagnose.

**Quality control:**  
URINORM control urines, Cat. No. REG00053 are designed for verification and confirmation of precision and accuracy of PHAN® diagnostic strips as well as LAURA® and LAURA® Smart. Bio-Rad Liquichek Urinalysis Control can be also used for quality control. It is recommended to perform QC measurements according to your local laboratory guidelines. More information about quality control can be found in LAURA® and LAURA® Smart user manuals.

**Storage conditions:** The recommended storage temperature is between (+2 to +25) °C. After first opening, the tube with remaining urine test strips has to be closed tightly immediately after removing of corresponding number of strips for evaluation and closed tube has to be stored away from direct sunlight at temperature not exceeding +30 °C. Storage above +30 °C will adversely affect the stability and test performance of the product.

**Waste disposal:** Used strip should be treated as potentially infectious and should be liquidated in accordance with local and national regulations relating to the safe handling of such materials. Let waste recycle or put it to municipal waste.

**Performance characteristics:** available at the manufacturer  
**Symbols on packaging:** see the table at the end of the instructions

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