ELite H580 Dil Haematology Analyser Diluent

Cat. No.	Product Name	Packaging (Content)
HEM00019	ELite H580 Dil	20 L, Instruction for use, RFID card

INTENDED USE

ELite H580 Dil is an *in vitro* diagnostic reagent for automated dilution of human blood for counting and sizing of blood cells on haematology analysers ELite 580. It is designed to report on the haematological status of patient samples, supporting the monitoring and diagnosis of pathological conditions. The ELite H580 Dil should be used by trained professionals in the haematology laboratory.

PRINCIPLE

The diluent is a solution with certain ion strength and conductivity, which can dilute the blood and form the sheath flow, providing a stable environment for blood cell counting.

DESCRIPTION AND COMPOSITION

Sodium Chloride; Sodium Sulphate; Buffering Agent; Antifungal and antibacterial agents.

PREPARATION

This reagent is ready for use and can be applied straight from the container; no special reagent preparation is necessary. Avoid reagent contaminations.

Before using, leave the reagent at room temperature for at least 24 hours.

For reagent connection to the analyser, refer to the ELite 580 User manual.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Haematology analyser ELite 580, Cat. No. INS00071

Reagents: ELite H580 Lyse1, Cat. No. HEM00020

ELite H580 Lyse2, Cat. No. HEM00021

ELite H580 Lyse3, Cat. No. HEM00022

SPECIFICATION

The product is a clear liquid without particles, precipitate or floccule. Background result: WBC $\leq 0.2 \times 10^9/L$, RBC $\leq 0.02 \times 10^{12}/L$, HGB ≤ 1 g/L, PLT $\leq 10 \times 10^9/L$.

STABILITY AND STORAGE

The product is stable till the expiry date mentioned on the label when stored at (+2 to +30) °C.

Discard opened container after 60 days.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Supplemental information:

EUH208 Contains reaction mass of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

WASTE MANAGEMENT

Please refer to local legal requirements.

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