## ELite H5 CON

### Haematology Control

Cat. No.	Product Name	Level	Packaging (Content)
HEM00071	ELite H5 CON	Low, Normal, High	3 × 3 mL, Instruction for use
HEM00024	ELite H5 CON L	Low	1 × 3 mL, Instruction for use
HEM00025	ELite H5 CON N	Normal	1 × 3 mL, Instruction for use
HEM00026	ELite H5 CON H	High	1 × 3 mL, Instruction for use

#### INTENDED USE

ELite H5 CON is an *in vitro* diagnostic reagent used for control of the accuracy and precision of automated haematology analysers ELite 580 and H560. For professional use in clinical laboratory only.

#### PRINCIPLE

It is an established laboratory practice to use a stable control to monitor the performance of diagnostic tests. This control is composed of stable materials that provide a means of monitoring the performance of haematology blood cell counters. It is sampled in the same manner as a patient specimen.

#### DESCRIPTION AND COMPOSITION

ELite H5 CON is an *in vitro* diagnostic material composed of human erythrocytes and mammalian leukocytes and platelets suspended in a plasma-like fluid with preservatives.

#### MEASURED PARAMETERS

Parameter	Abbreviation	Default Unit
White Blood Cell count	WBC count	10 <sup>9</sup> /L
Number of Neutrophils	Neu#	10 <sup>9</sup> /L
Number of Lymphocytes	Lym#	10 <sup>9</sup> /L
Number of Monocytes	Mon#	10 <sup>9</sup> /L
Number of Eosinophils	Eos#	10 <sup>9</sup> /L
Number of Basophils	Bas#	10 <sup>9</sup> /L
Percentage of Neutrophils	Neu%	%
Percentage of Lymphocytes	Lym%	%
Percentage of Monocytes	Mon%	%
Percentage of Eosinophils	Eos%	%
Percentage of Basophils	Bas%	%
Red Blood Cell count	RBC count	10 <sup>12</sup> /L
Hemoglobin Concentration	HGB concentration	g/L
Mean Platelet Volume	MPV	fL
Mean Corpuscular Volume	MCV	fL
Platelet count	PLT count	10 <sup>9</sup> /L

#### CALCULATED PARAMETERS

Parameter	Abbreviation	Default Unit
Mean Corpuscular Hemoglobin	MCH	pg
Mean Corpuscular Hemoglobin Concentration	MCHC	g/L
Red Blood Cell Distribution Width – Standard Deviation	RDW-SD	fL
Red Blood Cell Distribution Width – Coefficient of Variation	RDW-CV	%
Hematocrit	HCT	%
Platelet Distribution Width – Coefficient of Variation	PDW-CV	%
Platelet Distribution Width – Standard Deviation	PDW-SD	fL
Plateletcrit	PCT	%
Platelet-large cell ratio	P-LCR	%
Platelet-large cell count	P-LCC	10 <sup>9</sup> /L

#### PREPARATION

The haematology controls are ready to use.

#### STABILITY AND STORAGE

Store ELite H5 CON upright at 2–8 °C when not in use. Protect tubes from overheating and freezing. Unopened tubes are stable through the expiration date. Opened tubes are stable for 16 days at 2–8 °C, if they are handled properly.

#### MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Haematology analysers: ELite 580, Cat. No. INS00071  
H560, Cat. No. INS00078

Haematology calibrator: ELite H CAL, Cat. No. HEM00027

Assay values sheet, available electronically on the website <https://www.erbalachema.com/en/hematology-files/>

#### USED SYMBOLS

**REF** Catalogue number

**LOT** Lot number

**IVD** *In vitro* diagnostic medical device

 Manufacturer

 Consult instructions for use

**CONT** Content

#### ASSAY PROCEDURE

- Remove tubes from the refrigerator and allow them to warm to room temperature (15–25 °C) for 15 minutes before mixing.
- To mix, hold a tube horizontally between the palms of the hands. Do not premix on a mechanical mixer.
  - Roll the tube back and forth for 20–30 seconds; occasionally invert the tube. Mix vigorously, but do not shake.
  - Continue to mix in this manner until the red cells are completely suspended. Tubes stored for a long time may require extra mixing.
  - Gently invert the tube 8–10 times immediately before sampling.
- Analyze the sample as instructed in the Quality Control section of the User Manual for your instrument.
- After sampling:
  - If tube has been open for sampling, clean residual material from the cap and tube rim with a lint-free tissue. Replace the cap tightly.
  - Return tubes to refrigerator within 30 minutes of use.

Note: After mixing, product should be similar in appearance to fresh whole blood. In unmixed tubes, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration. Other discoloration, very dark red supernatant or unacceptable results may indicate deterioration. Do not use the product if deterioration is suspected.

#### EVALUATION

The ELite 580 or H560 analysers perform the evaluation automatically.

#### EXPECTED RESULTS

You can find the assay values sheet on the website <https://www.erbalachema.com/en/hematology-files/>. Search the assay values according the name of the lot. Verify that the lot number on the tube matches the lot number on the table of assay values. Assay values are determined on well-maintained, properly calibrated instruments using the instrument manufacturer's recommended reagents. Reagent differences, maintenance, operating technique, and calibration may contribute to inter-laboratory variation.

Assigned values are presented as a Mean and Range. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material.

Assay values of a new lot of control should be confirmed before the new lot is put into routine use. Test the new lot when the instrument is in good working order and quality control results on the old lot are acceptable. The laboratory's recovered mean should be within the assay range.

For greater control sensitivity each laboratory should establish its own mean and acceptable range and periodically reevaluate the mean. The laboratory range may include values outside of the assay range. The user may establish assay values not listed on the assay sheet, if the control is suitable for the method.

Erba Lachema recommends running the QC program on a daily basis with low, normal and high level controls. Perform the calibration, if the results of control bloods are out of the limitation.

#### TRACEABILITY

The controls are traceable to the internationally agreed reference method.

- WBC:** A series of 1:500 dilutions are made with calibrated glassware. Counting is performed on a Coulter Counter Z series instrument. All counts are corrected for coincidence.
- RBC:** A series of 1:50,000 dilutions are made with calibrated glassware. Counting is performed on a Coulter Counter Z series instrument. All counts are corrected for coincidence.
- HGB:** Haemoglobin value is determined by spectrophotometric procedure according to CLSI Standard H15-A3 and is traceable to ICSS/WHO International Haemoglobin-cyanide Standard.
- HCT:** Packed cell volume (PCV) is measured by the microhaematocrit procedure according to CLSI Standard H07-A3. No correction is made for trapped plasma.
- MCV:** MCV is the calibrated parameter instead of the HCT. The MCV is calculated from the HCT and RBC using the formula:  $MCV = HCT \times 10/RBC$ .
- PLT:** A series of 1:126 dilutions are made using calibrated glassware in 1 % ammonium oxalate. Platelets are counted using a haemocytometer and phase contrast microscopy.

#### PERFORMANCE CHARACTERISTICS

All performance characteristics of the ELite 580 and H560 haematology systems are listed in the user manual of a specific analyser in the specification chapter.

#### LIMITATIONS

The performance of this product is assured only if it is properly stored and used as described in this insert. Incomplete mixing of a tube prior to use invalidates both the sample withdrawn and any remaining material in the tube.

#### WARNING AND PRECAUTION

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

#### POTENTIAL BIOHAZARDOUS MATERIAL

Each human donor/unit used in the preparation of this product has been tested and found to be negative or non-reactive for the presence of HBsAg, Anti-HCV, NAT testing for HIV-1, HCV (RNA) and HIV-1/2. Each unit is also negative by a serological test for Syphilis (RPR or STS). Because no test method can offer complete assurance that infectious agents are absent, this material should be handled as potentially infectious. When handling or disposing of tubes follow precautions for patient specimens as specified in biosafety procedures.

#### Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents of the kit are not classified as dangerous.

#### WASTE MANAGEMENT

ELite H5 CON should be disposed as potentially infectious material. Please refer to local legal requirements.

#### REFERENCES

- Clinical Laboratory Standards Institute. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood: Approved Standard-Third Edition. CLSI document H15-A3.
- Clinical Laboratory Standards Institute. Procedure for Determining Packed Cell Volume by the Microhaematocrit Method: Approved Standard-Third Edition. CLSI document H07-A3.

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### REFERENCE

- Clinical Laboratory Standards Institute. Reference and Selected Procedures for the Quantitative Determination of Hemoglobin in Blood: Approved Standard-Third Edition. CLSI document H15-A3.
- Clinical Laboratory Standards Institute. Procedure for Determining Packed Cell Volume by the Microhematocrit Method: Approved Standard-Third Edition. CLSI document H07-A3.