Erba H560 Lyse2

Haematology Analyser Lyse

Cat. No.	Product Name	Packaging (Content)
HEM00032	Erba H560 Lyse2	500 mL, Instruction for use, RFID card

INTENDED USE

Erba H560 Lyse2 is an *in vitro* diagnostic reagent used to lyse red blood cells for quantitative 5-part differentiation of white blood cells in human blood on automated haematology analysers H560. It is designed to report on the haematological status of patient samples, supporting the monitoring and diagnosis of pathological conditions. The Erba H560 Lyse2 should be used by trained professionals in the haematology laboratory.

PRINCIPLE

The product is added to lyse the red blood cells (RBC) and maintain the morphology of cells for white blood cell (WBC) differential.

DESCRIPTION AND COMPOSITION

Surfactant; Buffering Agent; Antifungal and antibacterial agents.

PREPARATION

This reagent is ready for use and can be applied straight from the container; no special reagent preparation is necessary. Avoid reagent contaminations. Before using, leave the reagent at room temperature for at least 24 hours. For reagent connection to the analyser, refer to the H560 User manual.

MATERIAL REQUIRED BUT IS NOT PROVIDED WITH THE DEVICE

Haematology analyser H560, Cat. No. INS00078
Reagents: Erba H560 Dil, Cat. No. HEM00030
Erba H560 Lyse1, Cat. No. HEM00031

SPECIFICATION

The product is a clear liquid without particles, precipitate or floccule. Background result: WBC $\leq 0.2 \times 10^9/L$, RBC $\leq 0.02 \times 10^{12}/L$, HGB ≤ 1 g/L, PLT $\leq 10 \times 10^9/L$.

STABILITY AND STORAGE

The product will be stable till the expiry date mentioned on the label when stored at (+2 to +30) °C. Discard opened container after 60 days.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and / or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008
UFI: WWCF-C3YA-3N4R-SY41

Warning

Hazard statement:

H315 Causes skin irritation.
H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.
P302+P352 IF ON SKIN: Wash with plenty of water and soap.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Supplemental information:

EUH208 Contains reaction mass of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

WASTE MANAGEMENT

Please refer to local legal requirements.

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