FERRITIN

Cat. No.:	Pack name:	Packaging (Content):
BLT20015	FERR	1 x 24 ml Pufr 1 x 8 ml Latex

*Important note: R1 and R2 are lot dependent and not interchangeable.

INTENDED USE

Quantitative determination of Ferritin (FER) in human serum by turbidimetric immunoassay.

DIAGNOSTIC IMPLICATIONS

The plasma Ferritin concentration declines very early in the development of iron deficiency. On the other hand, a large number of chronic diseases result in increased serum Ferritin concentrations. These diseases include chronic infections, chronic inflammatory disorders such as rheumatoid arthritis or renal disease, Gaucher's disease, and numerous types of malignancies, especially lymphomas, leukaemias, breast cancer and neuroblastoma. Increase in plasma Ferritin concentration also occurs in viral hepatitis or following toxic liver injury because of release of Ferritin from damaged liver cells. Plasma Ferritin concentration is also increased with increases of iron stores, as seen in patients with hemosiderosis or hemochromatosis. Besides the use of Ferritin as an iron metabolism parameter, Ferritin as also gained importance as a tumour marker for therapeutic drug monitoring and follow-up.

METHOD

Measurement of antigen/latex-antibody reaction by the end-point method.

REAGENTS PROVIDED

Buffer
Hepes buffer
Sodium azide

Latex

Solution of suspended latex microparticles sensitized with rabbit IgG anti-human ferritin

Preparation and Stability of Reagents

Reagent Preparation

Liquid reagents, ready to use.

Stability and Storage

The reagents are stable until expiry date when kept at 2–8 °C. Stability in the instrument is at least 4 weeks if contamination is avoided. Do not freeze.

REAGENTS REQUIRED BUT NOT SUPPLIED

- Saline (9 g/l)
- Calibrators and Controls

Cat.No.	Product name	Pack name	Content
BLT20016	Ferritin Calibrator Set	FERR CAL SET	5 x 1 ml
BLT20017	Ferritin Control High	FERR CON H	1 x 1 ml

Pooled human serum, liquid and stabilized. Value is stated in the insert.

SAMPLE COLLECTION

Use fresh serum. If the test can not be carried out on the same day, the serum may be stored at 2–8 °C for 48 hours. If stored for a longer period, the sample should be frozen.

GENERAL ASSAY PROCEDURE

Application sheets for automated systems on clinical chemistry analyzers or manual procedures are available upon request.

Wavelength 600 nm

Sample, Control or calibrator	20 µl
Add Reagent 1	150 µl
Mix, incubate for 2 minutes, then add	
Reagent 2	50 µl
Mix and read absorbance, incubate and read absorbance after 5 minutes.	

REFERENCE VALUES

Men: 20–300 ng/ml

Woman: 15–200 ng/ml

This range is given for orientation only. Each Laboratory should establish his own reference values.

PERFORMANCES

The performance characteristics for the Ferritin reagents were measured on a clinical chemistry analyzer (ERBA XL) according to the procedure R1: 150 µl and R2: 50 µl.

Measuring Range: 0–500 ng/ml

LDD: 5 ng/ml

Hook effect: No risk

Precision:

[%CV]

	Low	Medium	High
Intra-Run	0.76	0.63	0.61
Inter-Run	4.14	3.80	4.11

Accuracy:

[ng/ml]

Control	Assigned	Measured
ERBA	85.2 (72.4–98.0)	92.1
ERBA	224 (190–258)	221.1
SIEMENS	87.7 (70.2–105.2)	82.5

Specificity: Monospecific

Interferences: No interference for Heparin (50 mg/dl), Na Ci (1000 mg/dl), Triglyceride (2500 mg/dl), Bilirubin (30 mg/dl), Heamoglobin (250 mg/dl), EDTA (5 mg/ml)

Limitations: None

Comparison with turbidimetry, external ferritin reagents:

$$y = 0.9261x + 15.596$$

$$r = 0.9868$$





PRECAUTIONS AND WARNINGS

- In vitro* diagnostic use only.
- Reagents of the kit are not classified as dangerous.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing, which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.

REFERENCES

- Lipzchitz DA, Cook JD, Finch CA. A clinical evaluation of serum ferritin as an index of iron stores. N Engl J Med. 1974; 290(22): 1213-1216
- Worwood M. Ferritin in human tissues and serum. Clin Heamatol. 1982; 11(2): 275-307
- Worwood M. Serum ferritin. Clin Sci(Lond)1986; 70(3): 215-220
- Warr GW, Magor KE, Higgins DA. IgY: clues to the origins of modern antibodies. Immunology Today 1995; 16: 92-8.

USED SYMBOLS

 REF	Catalogue Number	 IVD	In vitro Diagnostics	 i	See Instruction for Use
 LOT	Lot Number		Manufacturer	 CONT	Content
	Expiry Date		Storage Temperature		

 Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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