# CRP-HS CALIBRATOR

Lot: 0000000

Cat. No.	Pack Name	Packaging (Content)
XSYS0088	CRP-HS CAL	1 x 1 ml

## INTENDED USE

Preparation of reference curves for quantitative immunochemical high sensitive determination of CRP.

## COMPOSITION

A dilution of purified CRP with phosphate buffered saline. Ready to use.

## STORAGE AND STABILITY

The expiry date of the product at 2–8 °C is listed on the label. Do not freeze. After first opening the container, the serum can be used for at least 6 weeks if stored tightly closed at 2–8 °C after use. Do not freeze.

## PRECAUTIONS AND WARNING

- For *in Vitro* Diagnostic use.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.

## ASSIGNED VALUES

The value for CRP has been assigned to BCR – ERM-DA474/IFCC:

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