# APOLIPOPROTEIN B

Cat. No.	Pack Name	Packaging (Content)
XSYS0103	APO B	R1: 2 × 20.5 ml (Buffer), R2: 2 × 5.3 ml (Antiserum)

## INTENDED USE

Reagent kit for immunoturbidimetric determination of apolipoprotein B (Apo B) in human serum and plasma.

## CLINICAL SIGNIFICANCE

Apo B is the main protein component of LDL (Low Density Lipoprotein). Apo B is necessary for the reaction with LDL receptors in the liver and on cell walls and is thus involved in transporting cholesterol from the liver to the vessel cells.

Elevated levels of Apo B are frequently found in atherosclerotic vascular changes and are a risk factor for atherosclerosis.

## METHODOLOGY

Measurement of antigen-antibody reaction by the end-point method.

## REAGENT COMPOSITION

<b>R1</b> (Buffer)	
TRIS Buffer (pH 8)	200 mmol/l
PEG	56 g/l
Sodium chloride	50 mmol/l
Sodium azide	0.9 g/l

<b>R2</b> (Antiserum)	
Sterile delipidated goat serum against human Apo B	
HEPES buffer (pH 7.4)	50 mmol/l
EDTA	9 mmol/l
Sodium azide	0.9 g/l

## REAGENT PREPARATION

Liquid reagents, ready to use.

## STABILITY AND STORAGE

The reagents are stable until expiry date when kept at 2–8 °C. Stability in the instrument is at least 8 weeks if contamination is avoided. Do not freeze.

## SAMPLE COLLECTION

Serum or plasma (heparin, NaF, EDTA, sodium citrate), samples should be fresh and not hemolysed. Prior the analysis, the serum can be stored a week at 2–8 °C or 3 months at -20 °C. Avoid repeated freezing and thawing.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance accurately at 340 nm may be used
- Analyser specific consumables such as sample cups
- Calibrators
- Controls
- Saline (9 g/l NaCl)

## ASSAY PROCEDURE

Refer to the assay parameters for details.

## CALIBRATION

Blank: saline

Cat. No.	Product name	Pack name	Content
<b>XSYS0105</b>	APOLIPOPROTEIN CALIBRATOR SET	APO CAL SET	5 × 0.5 ml

Calibration curve: generate a 5 point calibration curve by using of set of calibrators APO CAL SET.

**Calibration frequency** it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

## QUALITY CONTROL

For quality control use

Cat. No.	Product name	Pack name	Content
<b>BLT00074</b>	ERBA LIPID NORM	ERBA LIPID NORM	3 × 5 ml
<b>BLT00075</b>	ERBA LIPID PATH	ERBA LIPID PATH	3 × 5 ml

## CALCULATION

Results are calculated automatically by the instrument.

## EXPECTED VALUES

### Reference Values

4–5 y	Male	58–103 mg/dl	resp. 0.58–1.03 g/l
	Female	58–104 mg/dl	resp. 0.58–1.04 g/l
6–11 y	Male	56–105 mg/dl	resp. 0.56–1.05 g/l
	Female	57–113 mg/dl	resp. 0.57–1.13 g/l
12–19 y	Male	55–110 mg/dl	resp. 0.55–1.10 g/l
	Female	53–119 mg/dl	resp. 0.53–1.19 g/l
20–29 y	Male	59–130 mg/dl	resp. 0.59–1.30 g/l
	Female	59–132 mg/dl	resp. 0.59–1.32 g/l
30–39 y	Male	63–143 mg/dl	resp. 0.63–1.43 g/l
	Female	70–132 mg/dl	resp. 0.70–1.32 g/l
40–49 y	Male	71–152 mg/dl	resp. 0.71–1.52 g/l
	Female	75–136 mg/dl	resp. 0.75–1.36 g/l
50–59 y	Male	75–160 mg/dl	resp. 0.75–1.60 g/l
	Female	75–168 mg/dl	resp. 0.75–1.68 g/l
60–69 y	Male	81–156 mg/dl	resp. 0.81–1.56 g/l
	Female	75–173 mg/dl	resp. 0.75–1.73 g/l
>69 y	Male	73–152 mg/dl	resp. 0.73–1.52 g/l
	Female	79–168 mg/dl	resp. 0.79–1.68 g/l

**Reference values are given for orientation only. Each laboratory should establish its own reference values.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Measuring Range:** 0–280 mg/dl (resp. 0–2.80 g/l)

**Limit of quantification:** 6.45 mg/dl (resp. 0.064 g/l)

**Hook Effect:** from 1700 mg/dl (resp. 17.0 g/l)

## PRECISION

Intra-assay precision Within run (n = 20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
<b>Sample 1</b>	82.2	3.85	4.69
<b>Sample 2</b>	121.3	3.44	2.84

Inter-assay precision Run to run (n = 20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
<b>Sample 1</b>	83.8	2.47	2.95
<b>Sample 2</b>	123.1	2.83	2.30

## COMPARISON

A comparison between XL-Systems APO B (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.955 x + 5.893 mg/dl

r = 0.959

## INTERFERENCES

Following substances do not interfere:

haemoglobin up to 3 g/l, bilirubin up to 30 mg/dl, triglycerides up to 1000 mg/dl

## WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use. To be handled by entitled and professionals educated person.
2. Reagents of the kit are not classified like dangerous.
3. Products from human source were tested and founded free from HbsAg and antibodies to HCV and HIV but this material should be treated just as carefully as potentially infective.
4. Product contains sodium azide <0.1 %.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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## REFERENCES / LITERATURA

1. Marcovina, S. M., Albers, International Federation of Clinical Chemistry: Standardization Project for Measurements of Apolipoproteins A1 and B. Clin. Chem. 37/10, (1991) 1676-1682
2. Marcovina, S. M., Albers, International Federation of Clinical Chemistry: Standardization Project for Measurements of Apolipoproteins A1 and B.III. Comparability of apolipoproteins A1 values by use of international reference material Clin. Chem. 39, (1993) 773-781
3. Tietz, N. W. Fundamentals of Clinical Chemistry. Saunders, Philadelphia 1987.
4. Burtis, C.A., Ashwood, E.R., Bruns, D.E. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 5th edition: W.B. Saunders Company Ltd., 2012.

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