Hb-Vario Kit

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EN

## INTENDED USE

The Hb-Vario Kit is intended for the *in vitro* quantitative determination of glycated human hemoglobin, referred to as Hemoglobin A1c or HbA1c. (IFCC mmol/mol and NGSP %) in human blood using ion-exchange high performance liquid chromatography (HPLC) on the Hb-Vario™ automated hemoglobin testing system (Hb-Vario™ Analyser). Hemoglobin A1c measurements are used as an aid in diagnosis of diabetes mellitus, as an aid to identify patients who may be at risk for developing diabetes mellitus, and for the monitoring of long-term blood glucose control in individuals with diabetes mellitus.

The Hb-Vario Kit and Hb-Vario™ automated hemoglobin testing system (Hb-Vario™ Analyser) are intended for Professional Use Only.

## Calibrators:

Hb-Vario HbA1c Calibrator Set CM is for the calibration of the Hb-Vario™ automated hemoglobin testing system for the quantitative determination of as Hemoglobin A1c in human whole blood.

## Controls:

Hb-Vario HbA1c Control Set is for the estimation of test precision and the detection of systematic analytical deviations of Hemoglobin A1c assay on the Hb-Vario™ automated hemoglobin testing system (Hb-Vario™ Analyser).

## SPECIAL CONDITIONS FOR USE STATEMENT(S)

The Hb-Vario HbA1c test should not be used to replace glucose testing in pediatric patients, pregnant women, or patients with Type 1 diabetes.

HbA1c should not be used to diagnose diabetes during pregnancy or to diagnose gestational diabetes. HbA1c reflects the average blood glucose levels over the preceding 3 months (the average life of a red blood cell), and therefore may be falsely low during pregnancy or any other condition associated with recent onset of hyperglycemia and/or decreased red cell survival.

The Hb-Vario HbA1c test should not be used to diagnose diabetes in patients with the following conditions:

- Any condition that alters the life span of the red blood cells, including recent blood loss, transfusion, significant iron deficiency, hemolytic anemia (including hereditary spherocytosis) or other hemolytic diseases, hemoglobinopathies and thalassemias, as the altered red blood cell turnover interferes with the relationship between mean blood glucose and HbA1c values.

- Malignancies or severe chronic hepatic and renal disease.

In cases of rapidly evolving Type 1 diabetes, the increase of HbA1c values might be delayed compared to the acute increase in glucose concentrations. In these conditions, diabetes mellitus must be diagnosed based on plasma glucose concentration and/or the typical clinical symptoms.

## CLINICAL SIGNIFICANCE

Human hemoglobin inside the erythrocytes found in blood undergoes a non-enzymatic chemical reaction with glucose and its metabolites.<sup>1</sup> The rate and extent of this reaction are thought to be dependent on the average blood glucose concentration during the lifetime of the erythrocyte<sup>2,3</sup>.

Glycated hemoglobin (hemoglobin A1c, HbA1c, A1c, or Hb1c) is a form of hemoglobin that is measured primarily to identify the average plasma glucose concentration over prolonged periods of time. It is formed in a non-enzymatic glycation pathway by hemoglobin's exposure to plasma glucose. Normal levels of glucose produce a normal amount of glycated hemoglobin. As the average amount of plasma glucose increases, the fraction of glycated hemoglobin increases in a predictable way. Because red blood cells in the human body survive for 8–12 weeks before renewal, measuring glycated hemoglobin (or HbA1c) can be used to reflect average blood glucose levels over that duration, providing a useful longer-term gauge of blood glucose control.

In diabetes mellitus, higher amounts of glycated hemoglobin, indicating poorer control of blood glucose levels, have been associated with cardiovascular disease, nephropathy, and retinopathy. Monitoring HbA1c in type 1 diabetic patients may improve outcomes.

Foetal hemoglobin (HbF) is the major form of hemoglobin in the foetus and neonate. Following birth, the level of HbF falls and is replaced by HbA in normal individuals. Within 12 months the level of HbF will have normally fallen <1 % of the total hemoglobin. Higher levels of HbF are often associated with a range of hereditary haematological disorders including sickle cell anemia, and thalassemia.

## PRINCIPLE

The Hb-Vario Analyser uses pressure cation exchange high performance liquid chromatography (HPLC) in conjunction with gradient elution to separate human hemoglobin subtypes and variants from hemolysed whole blood<sup>6,7</sup>. The separated hemoglobin fractions are monitored by absorption of light at 415 nm.

The chromatogram obtained is recorded and stored by the internal computer. A software program analyses the chromatogram and generates a results report on the monitor that can be printed.

## REAGENT COMPOSITION

## Reagent A

Phosphate Buffer (pH 7.00)

Sodium Azide

Preservatives

&lt;0.05 %

## Reagent B

Phosphate Buffer (pH 6.57)

Sodium Chloride

Sodium Azide

Preservatives

&lt;0.05 %

## Hemolysers Reagent

Borate buffer

&lt;3 %

## Cleaner

Tris Buffer (pH 8.00)

Subtilisin

Detergent

2.5 mmol/l

&lt;0.1 %

&lt;2 %

## Rinse

Sodium azide

0.05 %

## REAGENT PREPARATION

Reagent A and B, Hemolysers Reagent and Cleaner and Rinse are liquid, ready to use.

## STABILITY / STORAGE

Hb-Vario Kit is stored between 2 and 25 °C. Product is stable until the expiration date stated on the product label.

## SAMPLES

Blood samples should be collected into EDTA-coated tubes.

Whole blood samples can be directly positioned on the instrument rack. The analyzer performs automatically the hemolysate at the correct dilution ratio – refer to Hb Vario manual for instruction of hemolysate preparation. Samples are stable for seven days between 2 and 8°C.

## CALIBRATION

Calibration with Hb-Vario HbA1c Calibrator Set CM, REG00039, or Hb-Vario HbA1c Calibrator Set, REG00040 is recommended.

**Note: for Hb-Vario Kit it is necessary to order separately pack of Hb-Vario HbA1c Calibrator Set CM with the same LOT.**

## QUALITY CONTROL

For quality control Hb-Vario HbA1c Control Set, REG00041 is recommended.

**Note: for Hb-Vario Kit it is necessary to order separately pack of Hb-Vario HbA1c Control Set.**

## UNIT CONVERSION

IFCC = (NGSP - 2.15) / 0.915

NGSP = 0.915 × IFCC + 2.15

mmol/mol = 10 × IFCC

## REFERENCE VALUES

HbA1c (mmol/mol)

for non-diabetic

20–42

for glycemic control of a person with diabetes

43–60

According to IFCC (%)

for non-diabetic

2.0–4.2

for glycemic control of a person with diabetes

&lt;6.0

According to NGSP (%)

for non-diabetic

4.0–6.0

for glycemic control of a person with diabetes

&lt;7.0

## PERFORMANCE DATA

## Linearity:

Measuring range:

Linear regression:

y = 0.9985x + 0.0538 %      R<sup>2</sup> = 0.9999

17.7 %

3.8–17.7 %

## Precision:

Intra-assay (N=20)	Normal	Threshold	Pathological
Mean (%)	5.4	6.30	8.30
CV (%)	1.35	0.82	1.05

Inter-assay (N=80)	Normal	Threshold	Pathological
Mean (%)	5.05	5.97	7.89
CV (%)	3.12	2.76	2.2

## COMPARISON WITH COMMERCIAL METHOD (Tosoh G8)

Linear regression:

N = 77 samples; range: 5.15–15.23 %; Bias: +/-6 %

y = 0.94x - 0.2284%      R<sup>2</sup> = 0.9884

## INTERFERENCES

The data are summarized below:

Labile HbA1c – No interference up to 1500 mg/dl.

Triglycerides – No interference up to 1090 mg/dl.

Bilirubine – No interference up to 14.5 mg/dl.

Acetyl Salicylic acid – No interference up to 40 mg/dl.

Carb-Hb – interference in presence of Carb-Hb at the tested level of 20 mg/dl.

## HAEMOGLOBIN VARIANT INTERFERENCES

The data are summarized below:

HbF – No interference up to 10% of F

HbAS – No interference up to 34.6 % of S

HbAC – No interference up to 36.5 % of C

HbAD – No interference up to 37.7 % of D

HbAE – No interference up to 26.3 % of E

## CARRYOVER

The system demonstrates a carryover equal or less than 0.2 % HbA1c.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

**Hazards identification in accordance with Regulation (EC) No 1272/2008**

Hb-Vario Hemolysers Reagent

## Danger

Contains: boric acid

**Hazard statement:**

H360FD May damage fertility. May damage the unborn child.

**Precautionary statement:**

P280 Wear protective gloves / protective clothing.

P308 + P313 IF exposed or concerned: Get medical advice / attention.

UFI: WETA-Y3ND-WN4U-EM50

Hb-Vario Cleaner

## Warning

**Hazard statement:**

H319 Causes serious eye irritation.

H412 Harmful to aquatic life with long lasting effects.

**Precautionary statement:**

P273 Avoid release to the environment.

P280 Wear eye protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove

contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

UFI: 8HTA-G3AT-7N4A-3XR2

**Hb-Vario Reagent A, B**

**Supplemental information:**

EUH208 Contains reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

## Hb-Vario Rinse

Reagent of the kit is not classified as dangerous.

## FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water.

On eye contact rinse the eye quickly and thoroughly with the jet of tap water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## ASSAY PROCEDURE

1. USE EDTA-anticoagulated whole blood

2. Place the samples on the instrument rack

3. Set up a batch of samples for analysis

Refer to the Hb-Vario user Manual for detailed instructions on using the Hb-Vario Analyser.

## REFERENCES

1. Thorpe SR & Baynes JW(1982) in The Glycoconjugates (Horowitz M.ed) Vol.3, p113-132, Academic Press, New York, USA.
2. Rahbar S. (1968) Clin Chim Acta, 22, 296.
3. Koenig R.J. et al (1976) N Engl J Med 295, 417.
4. Garlick RL. et al (1983) J Clin Invest 71, 1062.
5. Ersser RS. et al (1986) Biomed. Chromatogr 1, 183
6. Ersser RS. et al (1987) J Auto Chem 9, 92.
7. Ersser R.S. et al (1991) Biomed. Chromatogr 5, 226.

## USED SYMBOLS

REF	Catalogue Number	IVD	<i>In vitro</i> diagnostic medical device	i	See Instructions for Use
LOT	Lot Number		Manufacturer	CONT	Content
	Expiry Date		Storage Temperature		

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