ALPHA AMYLASE

Cat. No.	Pack Name	Packaging (Content)
BLT00007	AMY 250	R1: 4 x 50 ml, R2: 1 x 50 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of alpha-Amylase in human serum, plasma and urine.

CLINICAL SIGNIFICANCE

α -Amylases are hydrolytic enzymes, which break down starch into maltose. In the human body α -amylases originate from various organs: the pancreatic amylase is produced by the pancreas and released into the intestinal tract, the salivary amylase is synthesized in the salivary glands and secreted into saliva. The amylase present in the blood is eliminated through the kidney and excreted into the urine. Therefore, elevation of serum activity is reflected in a rise of urinary amylase activity.

Measurement of α -amylase in serum and urine is mainly used for the diagnosis of pancreatic disorders as well as for detecting the development of complications. In acute pancreatitis the blood amylase activity increases within few hours after onset of abdominal pain, peaks after approx. 12 hours and returns to values within the reference range at the latest after 5 days. The specificity of α -amylase for pancreatic disorders is not very high as elevated levels are measured also in various non-pancreatic diseases, e.g. parotitis and renal insufficiency. Therefore, for confirmation of an acute pancreatitis measurement of lipase should be additionally performed.

PRINCIPLE

4,6 –Ethyliden – 4-nitrophenyl- α -D-maltoheptaoside serves as the substrate for determinativ of the catalytic concentration of alpha-amylase. After hydrolysis of the internal part of the chain of AMS and the rest of the chain by glukosidase, p-nitrophenol is gradually released. Its absorbance is measured at 405 nm. The absorbance increase is proportional to the catalytic concentration of alpha-amylase in the sample.

REAGENT COMPOSITION

R1

Good's buffer, pH 7.10	125 mmol/l
NaCl	62.5 mmol/l
MgCl ₂	12.5 mmol/l
α -glukosidase	> 42 μ kat/l

R2

Good's buffer, pH 7.10	500 mmol/l
4,6 –Ethyliden – 4-nitrophenyl-maltoheptaoside	8 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready for use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

Two reagents method – substrate start

Reagents R1 and R2 are liquid, ready to use.

After opening, reagents are stable until expiry date at 2–8°C if stored at appropriate conditions, closed carefully and without any contamination.

Monoreagent method – sample start

Mix 4 portion of reagent R1 with 1 portion of reagent R2.

Stability:	4 weeks	at 15–25°C	in the dark
	6 months	at 2–8°C	in the dark

SPECIMEN COLLECTION & HANDLING

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Serum, plasma (heparin, EDTA), urine.

Stability

in serum / plasma:	8 weeks	at 4–8°C
in urine:	26 weeks	at 4–8°C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator Lyonorm Calibrator is recommended.

QUALITY CONTROL

For quality control Lyonorm HUM N and Lyonorm HUM P are recommended.

UNIT CONVERSION

U/l x 0.017 = μ kat/l

EXPECTED VALUES³

fS alpha-amylase (μ kat/l) 37 °C up to 1.67

fU alpha-amylase (μ kat/l) 37 °C up to 8.35

The range of reference values is only approximate; it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

PERFORMANCES

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0,16 μ kat/l

Linearity: 36 μ kat/l

Measuring range: 0,16 – 36 μ kat/l

PRECISION

Intra-assay (n=20)	Mean (μ kat/l)	SD (μ kat/l)	CV (%)
Sample 1	1.31	0.013	0.96
Sample 2	3.06	0.040	1.30

Inter-assay (n=20)	Mean (μ kat/l)	SD (μ kat/l)	CV (%)
Sample 1	1.35	0.039	2.89
Sample 2	2.96	0.046	1.57

COMPARISON

A comparison between XL-Systems Amylase (y) and a commercially available test (x) using 40 samples gave following results:

N = 40

r = 0.989

y = 1.012 x + 0.030 μ kat/l

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 1000 mg/dl.

HEALTH PROTECTION

For in vitro diagnostic use. To be handled by entitled and professionally educated person. Reagents of the kit are not classified like dangerous.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials. Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

ASSAY PROCEDURE

Wavelength: 405 (400 - 420) nm

Cuvette: 1 cm

Temperature: 37°C

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

Two reagents method – substrate start

	Reagent blank	Calibrator	Sample
Reagent 1	0.800 ml	0.800 ml	0.800 ml
Sample	-	-	0.020 ml
Calibrator	-	0.020 ml	-
Distilled water	0.020 ml	-	-
Mix and after 1 min. incubation (at 37 °C) and			
Reagent 2	0.200 ml	0.200 ml	0.200 ml

Mix and measure the initial absorbance after 2 minutes from the addition of the reagent 2 (at 37 °C). Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA).

Monoreagent method – sample start

	Reagent blank	Calibrator	Sample
Reagent 1	1.00 ml	1.00 ml	1.00 ml
Sample	-	-	0.020 ml
Calibrator	-	0.020 ml	-
Distilled water	0.020 ml	-	-

Mix and measure the initial absorbance after 2 minutes from the addition of the reagent 2 (at 37 °C). Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA).

CALCULATION

$$1. \text{ alpha-amylase } (\mu\text{kat/l}) = \frac{\Delta A_{\text{sam}}}{\Delta A_{\text{cal}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

2. Using factor:

alpha-amylase (μ kat/l) = f x ΔA /min.

f = factor:

405 nm	Serum/plasma	Urine
Serum/plasma	94.5	187.5
Urine	75.9	150.3

NOTE

Saliva and skin contain alpha-amylase therefore never pipette reagents by mouth and avoid contamination of samples and reagents. Even trace contamination can affected results.

Applications for automatic analysers will be supplied on request.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erbamannheim.com

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