BILIRUBIN TOTAL & DIRECT

| Cat. No. | Pack Name | Packaging (Content) |
|----------|-------------|---|
| BLT00011 | BIL T&D 200 | R1: 2 × 40 mL (BIT), R2: 2 × 40 mL (BID), R3: 1 × 40 mL, instruction for use |

INTENDED USE

The kit is intended for *in vitro* photometric quantitative determination of bilirubin total and direct in human serum and plasma on various automatic systems. In combination with other parameters it is intended for screening, monitoring and diagnosis of liver diseases. For professional use in clinical laboratory only.

CLINICAL SIGNIFICANCE

Red blood cells at the end of their circulating lives are broken down in the reticuloendothelial system, mainly the spleen. The resulting heme is converted to bilirubin upon removal of iron. This process accounts for about 80% of the 500 μmol (292 mg) of bilirubin formed daily. Other sources of bilirubin include the breakdown of myoglobin and cytochromes and the catabolism of immature red blood cells in the bone marrow.

Once formed, bilirubin is transported to the liver bound to albumin. This fraction of bilirubin is referred to as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid (mono- and diglucuronides) by the enzyme uridyl diphosphate glucuronyl transferase to form conjugated bilirubin. Conjugated bilirubin or direct bilirubin is excreted via the biliary system into the intestine, where it is metabolized by bacteria to a group of products known collectively as stercobilinogen. Elimination is almost complete and serum levels are normally negligible.


Total bilirubin is the sum of the unconjugated and conjugated fractions. Total bilirubin is elevated in hepatitis, cirrhosis, hemolytic disorders, several inherited enzyme deficiencies, and conditions causing hepatic obstruction.

TOTAL BILIRUBIN = INDIRECT BILIRUBIN + DIRECT BILIRUBIN

PRINCIPLE

Total Bilirubin: bilirubin and bilirubin glucuronate react with diazotized sulphanilic acid in a strongly acidic medium and form an intensely coloured diazo dye – azobilirubin. Bilirubin glucuronate is soluble in water and reacts directly, bilirubin associated with albumin must be previously hydrolysed in presence of Cetrimide.

Direct Bilirubin: bilirubin reacts with diazotized sulphanilic acid in a strongly acidic medium and form an intensely coloured diazo dye – azobilirubin^{1,2,3,4,5}.

Absorbance of the Azobilirubin measured at 546 nm is proportional to the concentration of the total or direct bilirubin in the sample.

REAGENT DESCRIPTION AND COMPOSITION

| R1 | | R2 | |
|------------------|-------------|------------------|-------------|
| Sulphanilic acid | 28.9 mmol/L | Sulphanilic acid | 28.9 mmol/L |
| HCl | 58.8 mmol/L | HCl | 23.0 mmol/L |
| Cetrimide | 68.6 mmol/L | | |
| R3 | | | |
| Sodium nitrite | 2.90 mmol/L | | |

COMPOSITION OF REACTION MIXTURE

| Bilirubin Total | | Bilirubin Direct | |
|------------------|-------------|------------------|-------------|
| Sulphanilic acid | 22.0 mmol/L | Sulphanilic acid | 22.0 mmol/L |
| HCl | 44.8 mmol/L | HCl | 17.5 mmol/L |
| Cetrimide | 52.2 mmol/L | Sodium nitrite | 0.55 mmol/L |
| Sodium nitrite | 0.55 mmol/L | | |

REAGENT PREPARATION

Reagents are liquid, ready to use. Prepare working reagent by mixing of 4 portion of reagent R1 (BIT) or reagent R2 (BID) with 1 portion of reagent R3.

MATERIAL REQUIRED BUT NOT PROVIDED WITH THE DEVICE

Any instrument with temperature control of 37 ±0.5 °C that is capable of reading absorbance at 546 nm may be used, general laboratory equipment.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. The working reagent is stable for 7 days at 2–8 °C, when protected from contamination and light. It is recommended to prepare fresh working solution before assay is performed.

SPECIMEN COLLECTION AND HANDLING

It is recommended to follow ISO 15189 and laboratory instruction. For specimen collection and preparation only use suitable tubes or collection containers. Only the specimens listed below were tested and found acceptable.

Serum.
Plasma: Li-heparin and K₂-EDTA plasma.
The sample types listed were tested with a selection of sample collection tubes that were commercially available at the time of testing, i.e. not all available tubes of all manufacturers were tested. Sample collection systems from various manufacturers may contain differing materials which could affect the test results in some cases. When processing samples in primary tubes (sample collection systems), follow the instructions of the tube manufacturer.
Centrifuge samples containing precipitates before performing the assay.
See the limitations and interferences section for details about possible sample interferences.

| Stability in serum / plasma ^{6,7} : | | Bilirubin Direct: | |
|--|---|--|--|
| Bilirubin Total: | 1 day at 15–25 °C 7 days at 4–8 °C 6 months at -20 °C | 2 days at 15–25 °C 7 days at 4–8 °C 6 months at -20 °C | |

Protect sample from light.
Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL 4x3, Cat. No. XSYS0034 or XL MULTICAL 10x3, Cat. No. XSYS0122 is recommended.

2-point calibration (blank and calibrator): distilled water is recommended as blank

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5, Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

The control intervals and limits should be adapted according to each individual laboratory's requirements. Values obtained should fall within the defined intervals. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

TRACEABILITY

This method, calibrator XL MULTICAL and controls ERBA NORM and ERBA PATH have been standardized against the reference material NIST SRM 916 (Total Bilirubin) or Doumas method¹¹ (Direct Bilirubin).

ASSAY PROCEDURE

Wavelength: 546 / 670 nm
Cuvette: 1 cm

| | Reagent blank | Calibrator | Sample |
|-----------------|---------------|------------|---------|
| Working reagent | 1.00 mL | 1.00 mL | 1.00 mL |
| Sample | – | – | 0.05 mL |
| Calibrator | – | 0.05 mL | – |
| Distilled water | 0.05 mL | – | – |

Mix and incubate 5 min. at 37 °C. Measure absorbance of the sample A_{sam} and calibrator A_{cal} against reagent blank.

CALCULATION

$$1. \text{ Direct or Total Bilirubin (mg/dL)} = \frac{A_{\text{sam}}}{A_{\text{cal}}} \times C_{\text{cal}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

$$2. \text{ Using factor (f):} \quad \text{Direct or Total Bilirubin (mg/dL)} = A_{\text{sam}} \times f \quad \begin{matrix} \text{Total Bilirubin: } f = 24 \\ \text{Direct Bilirubin: } f = 16 \end{matrix}$$

ASSAY PARAMETERS FOR PHOTOMETERS

| Parameter | Total Bilirubin | Direct Bilirubin |
|-----------------------------|-----------------|------------------|
| Mode | End point | End point |
| Wavelength 1 (nm) | 546 | 546 |
| Wavelength 2 (nm) | 670 | 670 |
| Sample Volume (μL) | 25/50 | 25/50 |
| Working Reagent Volume (μL) | 500/1000 | 500/1000 |
| Incubation time (min.) | 5 | 5 |
| Reaction temperature (°C) | 37 | 37 |
| Reaction direction | Increasing | Increasing |
| Normal Low (mg/dL) | 0 | 0 |

| | | |
|-------------------------|---------|---------|
| Normal High (mg/dL) | 2 | 0.2 |
| Linearity Low (mg/dL) | 0.02 | 0.04 |
| Linearity High (mg/dL) | 30 | 20.5 |
| Factor | 24 | 16 |
| Blank with | Reagent | Reagent |
| Absorbance limit (max.) | 0.15 | 0.15 |
| Units | mg/dL | mg/dL |

UNIT CONVERSION

mg/dL × 17.1 = μmol/L

EXPECTED VALUES¹²

In serum:

| Total Bilirubin: | Premature | Full term |
|------------------|-----------|----------------|
| Cord | <2.0 | <2.0 mg/dL |
| 0–1 d | 1.0–8.0 | 2.0–6.0 mg/dL |
| 1–2 d | 6.0–12.0 | 6.0–10.0 mg/dL |
| 3–5 d | 10.0–14.0 | 4.0–8.0 mg/dL |
| Adult | | 0–2.0 mg/dL |

Direct Bilirubin:

Adults and infants: 0–0.2 mg/dL

It is recommended that each laboratory verifies this range or derives reference interval for the population it serves.

ANALYTICAL PERFORMANCE

Data contained within this section is representative for performance on ERBA XL-640 automatic system. Data obtained in your laboratory may differ from these values.

Limit of quantification:

Bilirubin Total: 0.02 mg/dL
Bilirubin Direct: 0.04 mg/dL

Limit of quantification represents the lowest measurable analyte level. It is calculated as the determined activity of diluted sample to have CV <20 % (n = 30)

Linearity:

Bilirubin Total: 30 mg/dL
Bilirubin Direct: 20.5 mg/dL

Linearity is the highest measured activity with recovery within ±10 % from theoretical value.

Precision:

Precision was determined by using controls in an internal protocol with repeatability (n = 20) and intermediate precision (2 aliquots per run, 2 run per day, 20 days). The following results were obtained:

| Bilirubin Total / Repeatability | Mean (mg/dL) | SD (mg/dL) | CV (%) |
|---------------------------------|--------------|------------|--------|
| Sample 1 | 1.56 | 0.007 | 0.44 |
| Sample 2 | 5.25 | 0.028 | 0.53 |

| Bilirubin Direct / Repeatability | Mean (mg/dL) | SD (mg/dL) | CV (%) |
|----------------------------------|--------------|------------|--------|
| Sample 1 | 0.59 | 0.005 | 0.80 |
| Sample 2 | 1.56 | 0.013 | 0.82 |

| Bilirubin Total Intermediate precision | Mean (mg/dL) | SD (mg/dL) | CV (%) |
|--|--------------|------------|--------|
| Sample 1 | 1.37 | 0.027 | 1.96 |
| Sample 2 | 4.76 | 0.094 | 1.98 |

| Bilirubin Direct Intermediate precision | Mean (mg/dL) | SD (mg/dL) | CV (%) |
|---|--------------|------------|--------|
| Sample 1 | 0.55 | 0.012 | 2.27 |
| Sample 2 | 1.62 | 0.033 | 2.01 |

Accuracy

Bilirubin Total: two different validated control materials were used. Determined bias is and -0.8 % at the target value 1.68 mg/dL and 6.7 % at the target value 4.55 mg/dL.

Bilirubin Direct: two different validated control materials were used. Determined bias is and 6.1 % at the target value 0.35 mg/dL and 7.9 % at the target value 0.99 mg/dL.

Comparison

A comparison between XL-640 automatic system BILIRUBIN TOTAL (y) and a commercially available test (x) using 150 samples gave following results:

Linear regression: $y = 1.014x - 0.056 \text{ mg/dL}$ $r = 0.998$

Passing-Bablok¹³: $y = 0.920x + 0.046 \text{ mg/dL}$ $r = 0.916$

A comparison between XL-640 automatic system BILIRUBIN DIRECT (y) and a commercially available test (x) using 120 samples gave following results:

Linear regression: $y = 0.922x + 0.154 \text{ mg/dL}$ $r = 0.992$

Passing-Bablok¹³: $y = 0.933x + 0.134 \text{ mg/dL}$ $r = 0.882$

Interferences

Criterion: Recovery within ±10 % of initial value of total or direct bilirubin concentration in the sample without interfering substance.

Following substances do not interfere:

Total Bilirubin: haemoglobin up to 12.5 g/L, triglycerides up to 850 mg/dL

Direct Bilirubin: haemoglobin up to 3 g/L, triglycerides up to 850 mg/dL.

Limitations:

- Deteriorated reagents (e.g. exceeding the storage temperature) may give incorrect results. Maximum allowable absorbance of the reagent blank measured at 546 nm against the distilled water is 0.15.

- High concentration of haemoglobin and triglycerides in sample can interfere with determination of direct bilirubin. See paragraph Interferences.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

Hazards identification in accordance with Regulation (EC) No 1272/2008

R1

UFI: 0V0F-WE2R-2H7R-9P9K

Danger

Hazard statements:

H314 Causes severe skin burns and eye damage.

Precautionary statements:

P260 Do not breathe vapours.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Supplemental information

EUH 208 Contains sulphanilic acid. May produce an allergic reaction.

R2

UFI: U61F-WEUA-9H7R-X1MU

Danger

Hazard statements:

H314 Causes severe skin burns and eye damage.

Precautionary statements:

P260 Do not breathe vapours.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301 + P330 + P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Supplemental information

EUH208 Contains sulphanilic acid. May produce an allergic reaction.

R3

Reagent is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Pearlman FC and Lee RTY. Detection and measurement of total bilirubin in serum, with use of surfactants as solubilizing agents. Clin Chem 20: 447-453, 1974.
2. Zoppi F, Peracino A, Fenili D, Marcovina S and Ramella C. Metodo per la determinazione della bilirubina totale e coniugata. Uso di un tensioattivo cationico come agente solubilizzante. Giorn It Chim Cl 1: 343-359, 1976.
3. Malloy HT, Evelyn KA: J. Biol. Chem. 119: 481-490, 1937.
4. Walters MI and Gerarde HW, An ultramicromethod for the determination of conjugated and total bilirubin in serum or plasma, Microchemical Journal. 15: 231-243, 1970.
5. Henry, RJ, Editor, Clinical Chemistry, Principles and Technics, p.1058, Harper and Row, Publishers, Hagerstown, Maryland, 1974.
6. Quality of Diagnostic Samples, Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory medicine, 3rd completely revised ed. 2010.
7. Use of Anticoagulants in Diagnostic Laboratory Investigations. WHO Publication WHO/DIL/LAB/99.1 Rev. 2: Jan 2002.
8. Balistreri WF, Shaw LM. Liver function. In: Tietz NW, ed. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia: WB Saunders: 729-761, 1987.
9. Thomas L ed. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft: 192-202, 1998.
10. Tolman KG, Rej R. Liver function. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1125-77, 1999.
11. Dumas BT, Perry BW, Jendrzyszczak B, et al. Pitfalls in the American Monitor Kit Methods for Determination of Total and Direct Bilirubin. Clin Chem 28 (11), 2305-2308, 1982.
12. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
13. Bablok W, Passing H, Bender R, et al. A general regression procedure for method transformation. Application of linear regression procedures for method comparison studies in clinical chemistry, Part III. J Clin Chem Clin Biochem, Nov;26(11): 783-790, 1988.

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