# BILIRUBIN TOTAL JENDRASSIK GROF

Cat. No.	Pack Name	Packing (Content)
BLT00012	BIL T JG 350	R1: 1 x 250 ml, R2: 1 x 90 ml, R3: 1 x 3 ml

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Total Bilirubin in human serum and plasma.

## CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of haemoglobin. Bilirubin formed in the reticulo endothelial system is transported bound by albumin to the liver. This bilirubin is water insoluble and is known as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid to form direct bilirubin. Conjugated bilirubin is excreted via the biliary system into the intestine. Here it is metabolised by bacteria to urobilinogen and stercobilinogen.

### TOTAL BILIRUBIN = INDIRECT BILIRUBIN + DIRECT BILIRUBIN

Total Bilirubin is elevated in obstructive conditions of the bile duct, hepatitis, cirrhosis, in haemolytic disorders and several inherited enzyme deficiencies.

Indirect Bilirubin is elevated by pre-hepatic causes such as haemolytic disorders or liver diseases resulting in impaired entry transport or conjugation within the liver.

Monitoring of indirect bilirubin in neonates is of special importance as it is the indirect (or free) bilirubin bound to albumin that is able to cross the blood brain barrier more easily increasing the danger of cerebral damage.

## PRINCIPLE

Reagent kit for determination of Total bilirubin in human serum and plasma according to Jendrassik-Grof modif. method, convenient for use on automatic analysers.

Nonconjugated (free) bilirubin is released from the bound to albumin in the presence of accelerator. Bilirubin (free and direct) then react in strongly acid medium with diazotized sulfanilic acid forming red coloured azobilirubin, useful for photometric determination.

## REAGENT COMPOSITION

### R1 (ACCELERATOR)

Sodium benzoate	387 mmol/l
EDTA	2.7 mmol/l
Caffeine	192 mmol/l
Sodium acetate	680 mmol/l

### R2

Hydrochloric acid	≥ 170 mmol/l
Sulfanilic acid	29 mmol/l

### R3

Sodium nitrite	72.5 mmol/l
----------------	-------------

### REACTION MIXTURE

Sodium benzoate	285.1 mmol/l
EDTA	2.0 mmol/l
Caffeine	141.8 mmol/l
Sodium acetate	501.4 mmol/l
Hydrochloric acid	≥ 34.7 mmol/l
Sulfanilic acid	5.9 mmol/l
Sodium nitrite	0.48 mmol/l

## REAGENT PREPARATION

Reagents are liquid, ready to use.

Reagent 1 is liquid ready to use.

Preparation of working solution:  
mix reagents R2 and R3 in ratio 31+1.  
Stability: 1 week at 2–8 °C in the dark  
5 hours at 15–25 °C in the dark

## STABILITY AND STORAGE

If stored at 2–25 °C, reagents are stable until expiry date, that is stated on the package.

After opening, all reagents are stable until expiry date at 2–25 °C if stored at appropriate temperature conditions, closed carefully and without any contamination.

## SPECIMEN COLLECTION & HANDLING

Use serum, plasma (heparin, EDTA). Protected from the light!

It is recommended to follow NCCLS procedures (or similar standardized conditions).

### Stability

2 days at 5–25 °C

7 days at 2–8 °C

3 months at -20 °C

**FREEZE ONLY ONCE!**

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator Lyonorm Calibrator is recommended.

## QUALITY CONTROL

For quality control, it is recommended to use following materials:

Lyonorm HUM N

Lyonorm HUM P

## UNIT CONVERSION

mg/dl x 17.1 = μmol/l

## EXPECTED VALUES <sup>4</sup>

newborn ≤ 171

adults and children (older than 1 month) 3.4 – 17.1

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 1.37 μmol/l

**Linearity:** 500 μmol/l

**Measuring range:** 1.37 – 500 μmol/l

Intra-assay precision Within run (n=20)	Mean (μmol/l)	SD (μmol/l)	CV (%)
Sample 1	24.63	0.34	1.38
Sample 2	89.1	0.96	1.08

Inter-assay precision Run to run (n=20)	Mean (μmol/l)	SD (μmol/l)	CV (%)
Sample 1	24.5	0.86	3.52
Sample 2	86.6	1.59	1.83

## COMPARISON

A comparison between BIL T JG 350 (y) and a commercially available test (x) using 40 samples

gave following results:

y = 1.022 x + 0.360 μmol/l

r = 0.997

## INTERFERENCES

Following substances do not interfere:

triglycerides up to 750 mg/dl, haemoglobin up to 10 g/l.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous.

Reagent 2

EUH210 Safety data sheet available on request.

EUH208 Contains sulphanilic acid. May produce an allergic reaction.

## FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

## WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

## PROCEDURE

Wavelength 546 (540–560) nm

Cuvette 1 cm

Temperature 37 °C

Serum/reaction mixture ratio 1/19

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

	Sample Blank	Sample	Standard (calibrator) Blank	Standard (calibrator)
Reagent R1	0.7 ml	0.7 ml	0.7 ml	0.7 ml
Sample	0.05 ml	0.05 ml	-	-
Standard (calibrator)	-	-	0.05 ml	0.05 ml

Mix immediately and incubate for 5 minutes in the dark, then add

Working solution	-	0.2 ml	-	0.2 ml
Reagent R2	0.2 ml	-	0.2 ml	-

Mix and after 5 minutes incubation in the dark read the absorbance of the sample blank  $A_1$ , the absorbance of the sample  $A_2$ , the absorbance of the standard blank  $A_3$  and the absorbance of the standard  $A_4$  against reagent blank. The coloration is stable 30 minutes protected from the light.

## CALCULATION

$$\text{Bilirubin } (\mu\text{mol/l}) = \frac{(A_2 - A_1)}{(A_4 - A_3)} \times C_{\text{st}}$$

$C_{\text{st}}$  = standard (calibrator) concentration

## NOTE

1. It is recommended to use 0.9% saline for reagent blank for automatic procedure, in case of manual measurement can be used distilled water.
2. The analysis are to be carried out in a dark place. In the light, the coloration of azobilirubin drops of about 10% during 15 minutes.
3. In case of manual procedure, measure at 15–25 °C.
4. Without measurement of sample and standard (calibrator) blank, false higher concentration of bilirubin might be found, by about 2–14 μmol/l in dependence of the nature of the sera.


**Applications for automatic analysers will be supplies on request.**

## REFERENCES

1. Jendrassik, L., Gróf, P.: Biochem. Z. 297, 81, 1938
2. Doumas, B. T., Kwok-Cheung, P. P., Perry, B. W. a kol.: Clin. Chem. 3, 1779, 1985
3. Garber, C. C.: Clin. Chem. 27, 1410, 1981
4. Tietz, N.W.: Textbook Of Clin. Chem., 1169–71, W.B.Saunders Co., Philadelphia, 1999

## USED SYMBOLS

 Lot Number  In vitro Diagnostics  See Instruction for Use

 Catalogue Number  Manufacturer  Content

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