BILIRUBIN DIRECT JENDRASSIK GROF

Cat. No.	Pack Name	Packing (Content)
BLT00014	BIL D JG 300	R1: 1 x 250 ml, R2: 1 x 45 ml, R3: 1 x 3 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Bilirubin Direct in human serum and plasma.

CLINICAL SIGNIFICANCE

Bilirubin is a breakdown product of haemoglobin. Bilirubin formed in the reticulo endothelial system is transported bound by albumin to the liver. This bilirubin is water insoluble and is known as indirect or unconjugated bilirubin. In the liver, bilirubin is conjugated to glucuronic acid to form direct bilirubin. Conjugated bilirubin is excreted via the biliary system into the intestine. Here it is metabolised by bacteria to urobilinogen and stercobilinogen.

TOTAL BILIRUBIN = INDIRECT BILIRUBIN + DIRECT BILIRUBIN

Total Bilirubin is elevated in obstructive conditions of the bile duct, hepatitis, cirrhosis, in haemolytic disorders and several inherited enzyme deficiencies.

Indirect Bilirubin is elevated by pre-hepatic causes such as haemolytic disorders or liver diseases resulting in impaired entry transport or conjugation within the liver.

Monitoring of indirect bilirubin in neonates is of special importance as it is the indirect (or free) bilirubin bound to albumin that is able to cross the blood brain barrier more easily increasing the danger of cerebral damage.

PRINCIPLE

Reagent kit for determination of direct bilirubin in human serum and plasma according to Jendrassik-Grof modif. method, convenient for use on automatic analysers.

Direct bilirubin in strongly acid medium reacts with diazotized sulfanilic acid forming red coloured azobilirubin, useful for photometric determination.

REAGENT COMPOSITION

R1
Hydrochloric acid 50.6 mmol/l

R2
Hydrochloric acid ≥ 170 mmol/l
Sulfanilic acid 29 mmol/l

R3
Sodium nitrite 72.5 mmol/l

REACTION MIXTURE

Hydrochloric acid ≥ 62.3 mmol/l
Sulfanilic acid 3.7 mmol/l
Sodium nitrite 0.36 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

Reagent 1 is liquid ready to use.

Preparation of working solution:

mix reagents R2 and R3 in ratio 26+1.

Stability: 1 week at 2–8°C in the dark
5 hours at 15–25°C in the dark

STABILITY AND STORAGE

If stored at 2–25°C, reagents are stable until expiry date, that is stated on the package.

After opening, all reagents are stable until expiry date at 2–25°C if stored at appropriate temperature conditions, closed carefully and without any contamination.

SPECIMEN COLLECTION & HANDLING

Use serum, plasma (heparin, EDTA). Protected from the light!

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

2 days at 15–25°C

7 days at 2–8°C

3 months at -20°C

FREEZE ONLY ONCE!

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator Lyonorm Calibrator.

QUALITY CONTROL

For quality control, it is recommended to use following materials:

Lyonorm HUM N

Lyonorm HUM P

UNIT CONVERSION

mg/dl x 17.1 = µmol/l

EXPECTED VALUES ⁴

Adults and infants : ≤ 3.4 µmol/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves..

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 1.46 µmol/l

Linearity: 180 µmol/l

Measuring range: 1.46 – 180 µmol/l

Intra-assay precision Within run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	11.50	0.394	3.43
Sample 2	44.9	1.67	3.72

Inter-assay precision Run to run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	11.56	0.378	3.28
Sample 2	42.5	1.40	3.30

COMPARISON

A comparison between BIL D JG 300 (y) and a commercially available test (x) using 40 samples

gave following results:

y = 0.997 x + 0.220 µmol/l

r = 0.999

INTERFERENCES

Following substances do not interfere:

triglycerides up to 750 mg/dl, haemoglobin up to 10 g/l.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous.

Reagent 2:

EUH210 Safety data sheet available on request.

EUH208 Contains sulphanilic acid. May produce an allergic reaction.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

PROCEDURE

Wavelength 546 (540–560) nm

Cuvette 1 cm

Temperature 37 °C

Serum/reaction mixture ratio 1/15

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

	Sample Blank	Sample	Standard (calibrator) Blank	Standard (calibrator)
Reagent R1	1.2 ml	1.2 ml	1.2 ml	1.2 ml
Sample	0.1 ml	0.1 ml	-	-
Standard (calibrator)	-	-	0.1 ml	0.1 ml
Mix immediately and incubate for 5 minutes in the dark, then add				
Working solution	-	0.2 ml	-	0.2 ml
Reagent R2	0.2 ml	-	0.2 ml	-
Mix and after 5 minutes incubation in the dark read the absorbance of the sample blank A ₁ , the absorbance of the sample A ₂ , the absorbance of the standard blank A ₃ and the absorbance of the standard A ₄ against reagent blank. The coloration is stable 30 mins protected from the light.				

CALCULATION

$$\text{Bilirubin } (\mu\text{mol/l}) = \frac{(A_2 - A_1)}{(A_4 - A_3)} \times C_{\text{st}}$$

C_{st} = standard (calibrator) concentration

NOTE

- It is recommended to use 0.9% saline for reagent blank for automatic procedure, in case of manual measurement can be used distilled water.
- The analysis are to be carried out in a dark place. In the light, the coloration of azobilirubin drops of about 10% during 15 minutes.
- In case of manual procedure, measure at 15–25°C.
- Without measurement of sample and standard (calibrator) blank, false higher concentration of bilirubin might be found, by about 2-14 µmol/l in dependence of the nature of the sera.
- Procedure in case of calibration with Bilirubin Standard or with using calibration curve obtained by BLT Bilirubin Total Liquid 350 S is described in Russian version of this instruction. In other languages it is available on request.

Applications for automatic analysers will be supplies on request.

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