CREATININE

Cat. No.	Pack Name	Packaging (Content)
BLT00020	CREA 200	R1: 2 x 50 ml, R2: 2 x 50 ml, R3 STD: 1 x 5 ml
BLT00019	CREA 1000	R1: 2 x 250 ml, R2: 2 x 250 ml, R3 STD: 2 x 5 ml

EN

IVD

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Creatinine in human serum, plasma and urine.

CLINICAL SIGNIFICANCE

Creatinine is a waste product formed in muscle from the high energy storage compound, creatine phosphate. The amount of creatinine produced is fairly constant (unlike Urea) and is primarily a function of muscle mass. It is not greatly affected by diet, age, sex or exercise. Creatinine is removed from plasma by glomerular filtration and then excreted in urine without any appreciable resorption by the tubules.

Creatinine is used to assess renal function, however, serum creatinine levels do not start to rise until renal function has decreased by at least 50 %.

PRINCIPLE

Creatinine reacts with alkaline picrate to produce an orange-red complex (Jaffe reaction). This is a non-specific reaction and is given by many other substances. Specificity of the assay has been improved by the introduction of a kinetic method, however, the cephalosporin antibiotics are still major interferents.

The rate of complex formation or change of absorbance, measured at 490–510 nm, is proportional to the creatinine concentration in the sample.

REAGENT COMPOSITION

R1

Sodium Hydroxide 394 mmol/l

R2

Picric Acid 11 mmol/l

R3 standard see bottle label

WORKING REAGENT PREPARATION

Reagents are liquid, ready to use.

Mix equal volumes of Reagent 1 and Reagent 2. Wait for 15 minutes before use.

The Working Reagent is stable for 1 month at 2–8 °C. The absorbance of the reagent blank should be approximately <0.15 at 505 nm when read against distilled water.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–25 °C.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum, plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum / plasma: 7 days at 4–25 °C
at least 3 months at -20 °C

in urine: 2 days at 20–25 °C
6 days at 4–8 °C
6 months at -20 °C

For the determination in urine use 24 hours specimen. It is important to exactly measure the volume of collected urine. Dilute urine samples in 1+19 ratio with distilled water and multiply results by 20.

Discard contaminated specimens.

CALIBRATION

Calibration with the standard included in the kit or with the calibrator XL MULTICAL 4x3 (Cat. No. XSYS0034) or XL MULTICAL 10x3 (Cat. No. XSYS0122) is recommended.

QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat. No. BLT00080) or ERBA NORM 10x5 (Cat. No. XSYS0123) and ERBA PATH 4x5 (Cat. No. BLT00081) or ERBA PATH 10x5 (Cat. No. XSYS0124) are recommended.

UNIT CONVERSION

mg/dl x 88.4 = µmol/l

EXPECTED VALUES ³

	Serum	Urine
Males:	0.7–1.3 mg/dl	14–26 mg/kg/day
Females:	0.6–1.1 mg/dl	11–20 mg/kg/day
Newborn:	0.3–1.0 mg/dl	
Infant:	0.2–0.4 mg/dl	8–20 mg/kg/day
Child:	0.3–0.7 mg/dl	8–22 mg/kg/day
Adolescent:	0.5–1.0 mg/dl	8–30 mg/kg/day

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA CHEM 7 instrument. Data obtained in your laboratory may differ from these values.

Limit of quantification: 0.14 mg/dl

Linearity: 40 mg/dl

Measuring range: 0.14–40 mg/dl

PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	1.45	0.02	1.56
Sample 2	2.37	0.04	1.74

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	1.48	0.03	1.77
Sample 2	2.39	0.04	1.58

COMPARISON

A comparison between XL-Systems Creatinine (y) and a commercially available test (x) using 70 samples gave following results:

y = 1.0515 x - 0.0636 mg/dl

r = 0.9986

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 10 mg/dl, triglycerides up to 1700 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 127 2/2008

R1: UF! JCCF-NE95-FH74-0ADF

Reagent 1 contains < 1.6 % sodium hydroxide.

Warning

Hazard statement:

H315 Causes skin irritation

H319 Causes serious eye irritation

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Reagents R2 and R3 of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength: 492 (490–510) nm

Cuvette: 1 cm

	Reagent blank	Standard (Calibr.)	Sample
Working Reagent	1.00 ml	1.00 ml	1.00 ml
Sample	–	–	0.100 ml
Standard (Calibr.)	–	0.100 ml	–
Distilled water	0.100 ml	–	–

Mix well and read initial absorbance (A_1) 30 seconds after mixing and final absorbance (A_2) 60 seconds after mixing.

CALCULATION

$\Delta A = A_2 - A_1$

$$\text{Creatinine (mg/dl)} = \frac{\Delta A_{\text{sam}}/\text{min.} - \Delta A_{\text{bl}}/\text{min.}}{\Delta A_{\text{st}}/\text{min.} - \Delta A_{\text{bl}}/\text{min.}} \times C_{\text{st}}$$

C_{st} = standard (calibrator) concentration

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	2-point
Wavelength 1 (nm)	505
Sample Volume (µl)	50/100
Working Reagent Volume (µl)	500/1000
Lag time (sec.)	30
Kinetic interval (sec.)	60
No. of readings	1
Reaction temperature (°C)	37
Reaction direction	Increasing
Normal Low (mg/dl)	0.7
Normal High (mg/dl)	1.3
Linearity Low (mg/dl)	0.14
Linearity High (mg/dl)	40
Blank with	Reagent
Absorbance limit (max.)	0.15
Concentration of Standard	See bottle label
Units	mg/dl

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erbamannheim.com

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1. Fabiny DL and Ertinghausen G. Clin Chem, 1971 : 71 : 391.
2. Kroll MH and Elin RJ. Clin Chem, 1983 : 29 : 2044.
3. Tietz Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.
4. Young DS, Effects of Drugs on Clinical Laboratory Tests. Third Edition 1990; 3 : 122-32.
5. Henry RJ (Ed), Clinical Chemistry : Principles and Technics (2nd Ed), Clinical Chemistry : Principles and Technics (2nd Ed), Harper and Row, 1974; 548-551.
6. Wachtel et al, Creation and Verification of Reference Intervals. Laboratory Medicine 1995; 26-593-7.
7. National Committee for Clinical Laboratory Standards. User Evaluation of Precision Performance of Clinical Chemistry Devices. NCCLS; 1984, NCCLS Publication EP5-T.

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