CREATININE

Cat. No.	Pack Name	Packaging (Content)
BLT00022	CREA 500	R1: 4 x 100 ml, R2: 1 x 100 ml, R3: 1 x 10 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Creatinine in human serum, plasma and urine by Jaffé method.

CLINICAL SIGNIFICANCE

Creatinine is a waste product formed in muscle from the high energy storage compound, creatine phosphate. The amount of creatinine produced is fairly constant (unlike Urea) and is primarily a function of muscle mass. It is not greatly affected by diet, age, sex or exercise. Creatinine is removed from plasma by glomerular filtration and then excreted in urine without any appreciable resorption by the tubules.

Creatinine is used to assess renal function, however, serum creatinine levels do not start to rise until renal function has decreased by at least 50 %.

PRINCIPLE

Creatinine reacts with alkaline picrate to produce a reddish colour (Jaffe reaction). This is a non-specific reaction and is given by many other substances. Specificity of the assay has been improved by the introduction of a kinetic method¹, however, the cephalosporin antibiotics are still major interferents.

REAGENT COMPOSITION

R1	
Sodium Hydroxide	0,24 mol/l
R2	
Picric Acid	26 mmol/l
R3	
Zero calibrator	

COMPOSITION OF REACTION MIXTURE

Sodium hydroxide	0,183 mol/l
Picric acid	5 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

Two – reagent method

Reagents R1, R2 and R3 are liquid, ready to use.

If stored at 2–8°C, reagents are stable until expiry date, that is stated on the package.

After opening, reagents are stable until expiry date at 2–8°C if stored at appropriate conditions, in the dark place, closed carefully and without any contamination.

Mono – reagent method

Mix 4 portions of reagent R1 with 1 portion of reagent R2.

Stability: 1 week at 2–25°C in the dark.

SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability in serum / plasma:

7 days	at 4–25°C
at least 3 months	at -20°C

Stability in urine:

2 days	at 20–25°C
6 days	at 4–8°C
6 months	at -20°C

For the determination in urine use 24 hours specimen. It is important to exactly measure the volume of collected urine. Dilute urine samples in 1+19 ratio with distilled water and multiply results by 20.

Discard contaminated specimens.

CALIBRATION

For determination in serum and plasma, it is recommended two-point linear calibration with Lyonorm Calibrator and reagent R3 from the kit as 0,0 (zero) calibrator (see note 1).

QUALITY CONTROL

For quality control it is recommended to use LYONORM HUM N, Cat. No. BLT00070 and LYONORM HUM P, Cat. No. BLT00071.

UNIT CONVERSION

mg/dl x 88.4 = µmol/l

EXPECTED VALUES ⁴

fS Creatinine (µmol/l)	
male	55–110
female	44–95
dU Creatinine (mmol/24 hrs)	5–18

The range of reference values is only approximate; it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 7,07 µmol/l

Linearity: 1591 µmol/l

Measuring range: 7,07–1591 µmol/l

PRECISION

Intra-assay precision Within run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	262.5	3.80	1.45
Sample 2	396.9	4.60	1.16

Inter-assay precision Run to run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	110.5	1.86	1.71
Sample 2	292.6	2.74	0.95

COMPARISON

A comparison between CREAT 500 (y) and a commercially available test (x) using 40 samples gave following results:

r = 0.997

y = 0.998 x + 2.564 µmol/l

INTERFERENCES

Following substances do not interfere:

bilirubin up to 15 mg/dl, hemoglobin up to 10 g/l, triglycerides up to 1000 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

The reagent R1 contains irritating 1.0 % sodium hydroxide.

Warning

Hazard statement:

H315 Causes skin irritation

H319 Causes serious eye irritation

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

PROCEDURE

Wavelength: 492 (490–510) nm

Cuvette: 1 cm

Temperature: 37 °C

Serum/reaction mixture ratio 1/21

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

Two-reagent method

	Reagent blank	Calibrator (Standard)	Sample
Reagent R1	0,80 ml	0,80 ml	0,80 ml
Sample	–	–	0,05 ml
Calibrator (Standard)	–	0,05 ml	–
Reagent R3* (Distilled water)	0,05 ml	–	–

Mix and incubate 1–5 min. Then add:

Reagent R2	0,20 ml	0,20 ml	0,20 ml
------------	---------	---------	---------

Mix and after exactly 1 minute incubation read the initial absorbance for blank A_{bl} , sample A_{sam} and calibrator (standard) A_{st} . Exactly after 2 minutes read the final absorbance of blank A_{bl} , sample A_{sam} and calibrator (standard) A_{st} . Calculate resulting absorbance as the difference between the final and initial absorbance ($\Delta A/min$).

Mono-reagent method

	Reagent blank	Calibrator (Standard)	Sample
Working reagent	1,00 ml	1,00 ml	1,00 ml
Sample	–	–	0,05 ml
Calibrator (Standard)	–	0,05 ml	–
Reagent R3* (Distilled water)	0,05 ml	–	–

Mix and after exactly 1 minute incubation read the initial absorbance for blank A_{bl} , sample A_{sam} and calibrator (standard) A_{st} . Exactly after 2 minutes read the final absorbance of blank A_{bl} , sample A_{sam} and calibrator (standard) A_{st} . Calculate resulting absorbance as the difference between the final and initial absorbance ($\Delta A/min$).

* see note 1

CALCULATION

$$\text{Creatinine } (\mu\text{mol/l}) = \frac{\Delta A_{sam}/min. - \Delta A_{bl}/min.}{\Delta A_{st}/min. - \Delta A_{bl}/min.} \times C_{st}$$

C_{st} = standard (calibrator) concentration

NOTE

The determination of serum/plasma creatinine passed on Jaffe method is not fully specific. Alkaline picrate reacts also with some others substances present in serum matrix. Therefore, to correct the matrix effects, it is recommended to perform two-point calibration with Lyonorm calibrator and zero calibrator (reagent R3). Zero calibrator is used as STD1 (blank) instead of water or saline.

Applications for automatic analysers are available on request.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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