# GLUCOSE

Cat. No.	Pack Name	Packaging (Content)
BLT00025	GLU 500	R1: 2 × 250 ml, R2 standard: 2 × 5 ml
BLT00026	GLU 4×250	R1: 4 × 250 ml, R2 standard: 2 × 5 ml
BLT00027	GLU 1000	R1: 1 × 1000 ml

## INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Glucose in human serum, plasma and urine.

## CLINICAL SIGNIFICANCE



Accurate measurement of glucose in body fluid is important in diagnosis and management of diabetes, hypoglycemia, adrenal dysfunction and various other conditions. High levels of serum glucose may be seen in case of Diabetes mellitus, in patients receiving glucose containing fluids intravenously, during severe stress and in cerebrovascular accidents.

Decreased levels of glucose can be due to insulin administration, as a result of insulinoma, inborn errors of carbohydrate metabolism or fasting.

## PRINCIPLE

Trinder's method

Glucose in the sample is oxidised to yield gluconic acid and hydrogen peroxide in the presence of Glucose oxidase. The enzyme peroxidase catalyses the oxidative coupling of 4-aminoantipyrine with phenol to yield a coloured quinonemine complex, with absorbance proportional to the concentration of glucose in sample.

## REAGENT COMPOSITION

### R1

Phosphate buffer	250 mmol/l
Glucose oxidase	> 25 U/ml
Peroxidase	> 2 U/ml

Phenol	5 mmol/l
4-aminoantipyrine	0.5 mmol/l
<b>R2 standard</b>	See bottle label

## REAGENT PREPARATION

Reagents are liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

## SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum, plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

**Stability after addition of a glycolytic inhibitor (Fluoride, monoiodoacetate, mannose):** <sup>4</sup>

2 days	at 20–25 °C
7 days	at 4–8 °C

**Stability in serum (separated from cellular contents, hemolysis free) without adding a glycolytic inhibitor:** <sup>2,5</sup>

8 hours	at 25 °C
72 hours	at 4 °C

**Stability of glucose in urine:** 24 hours at 4–8 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with the standard included in the kit or with the calibrator XL MULTICAL 4x3 (Cat. No. XSYS0034) or XL MULTICAL 10x3 (Cat. No. XSYS0122) is recommended.

## QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat. No. BLT00080) or ERBA NORM 10x5 (Cat. No. XSYS0123) and ERBA PATH 4x5 (Cat. No. BLT00081) or ERBA PATH 10x5 (Cat. No. XSYS0124) are recommended.

## UNIT CONVERSION

mg/dl × 0.056 = mmol/l

## EXPECTED VALUES <sup>2</sup>

### Serum:

#### Glucose Fasting:

Cord:	45–96 mg/dl
Newborn, 1 d:	40–60 mg/dl
Newborn, >1 d:	50–80 mg/dl
Child:	60–100 mg/dl
Adult:	74–100 mg/dl
> 60 y:	82–115 mg/dl
> 90 y:	75–121 mg/dl
<b>Glucose 2 h Postprandial:</b>	<120 mg/dl
<b>WB (Hep) Adult:</b>	65–95 mg/dl
<b>Urine:</b>	1–15 mg/dl

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 2.34 mg/dl

**Linearity:** 450 mg/dl

**Measuring range:** 2.34–450 mg/dl

## PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
<b>Sample 1</b>	153.75	1.61	1.05
<b>Sample 2</b>	239.64	3.93	1.66

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
<b>Sample 1</b>	97.14	1.61	1.64
<b>Sample 2</b>	258.39	3.04	1.19

## COMPARISON

A comparison between XL-Systems Glucose (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.000 x - 0.714 mg/dl

r = 0.999

## INTERFERENCES

Following substances do not interfere:

haemoglobin up to 7.5 g/l, bilirubin up to 30 mg/dl, triglycerides up to 750 mg/dl. Interference by N-acetylcysteine (NAC), acetaminophen and metamizole causes falsely low results. To carry out the test, blood withdrawal should be performed prior to administration of drugs.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

**Hazards identification in accordance with Regulation (EC) No 1272/2008**

UFI: WY5U-AJ46-VE76-0WED

Warning

## Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

## Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

Reagent R2 is not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## ASSAY PROCEDURE

**Wavelength** (490–550) maximum at 500 nm

**Cuvette** 1 cm

	Reagent blank	Standard (Cal.)	Sample
Reagent 1	1000 µl	1000 µl	1000 µl
Sample	-	-	10 µl
Standard (Cal.)	-	10 µl	-
Distilled water	10 µl	-	-

Mix and incubate 5–10 min. at 37 °C. Measure absorbance of the sample A<sub>sam</sub> and standard A<sub>st</sub> against reagent blank.

## CALCULATION

$$\text{Glucose (mg/dl)} = \frac{\Delta A_{\text{sam}}}{\Delta A_{\text{st}}} \times C_{\text{st}} \quad C_{\text{st}} = \text{standard (calibrator) concentration}$$

**Applications for automatic analysers are available on request.**

## ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point
Wavelength 1 (nm)	505
Wavelength 2 (nm)	670
Sample Volume (µl)	5/10
Reagent Volume (µl)	500/1000
Incubation time (min.)	5
Incubation temp. (°C)	37
Normal Low (mg/dl)	74
Normal High (mg/dl)	100
Linearity Low (mg/dl)	2.34
Linearity High (mg/dl)	450
Concentration of Standard	See bottle label
Blank with	Reagent
Absorbance limit (max.)	0.2
Units	mg/dl

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: diagnostics@erba.com, www.erbamannheim.com

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# REFERENCES / ЛИТЕРАТУРА / ЛІТЕРАТУРА / LITERATURA / LITERATÚRA / REFERENCIAS

1. Thomas L.: Clinical Laboratory Diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998, p. 131 - 7.
2. Tietz N. W., (Ed.),Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.
3. Barham, D., Trinder, P.: An improved color reagent for the determination of blood glukose by the oxidase system. Analyst, 1972, 97; 142 - 5.
4. Guder WG, Zawta B et al. The quality of Diagnostic Samples. 1st ed. Darmstadt: GIT verlag; 2001;p.30-1.
5. Snacks DB, Bruns DE, Goldstein DE, Mac Laren NK, Mc Donald JM, Parrott M. Guidelines and recommendations for laboratory analysis in the diagnosis and management of diabetes mellitus. Clin Chemi 2002; 48:436-72.

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