CHLORIDE

Cat. No.	Pack Name	Packaging (Content)
BLT00033	CL 250	R1: 1 × 250 ml, R2 standard: 1 × 5 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Chloride in human serum, plasma and urine.

CLINICAL SIGNIFICANCE

Chloride is estimated in conjunction with the other electrolytes and used to support interpretation of the other electrolytes. The anion gap (Na + K) – (Cl + HCO₃) between these electrolytes is frequently affected in disease states. Increased levels can occur in nephritis, prostatic obstruction, eclampsia and dehydration. Decreased levels can occur with impaired gastrointestinal or renal function.

PRINCIPLE

When chloride is mixed with a solution of undissociated mercuric thiocyanate, the chloride preferentially combines with mercury forming mercuric chloride. The thiocyanate that is released then combines with ferric ions present in the solution forming strongly coloured ferric thiocyanate with an absorption maxima at 480 nm.

REAGENT COMPOSITION

R1	
Mercuric Thiocyanate	1.01 mmol/l
Ferric Nitrate	37.13 mmol/l
Nitric Acid	62.43 mmol/l
Mercuric Nitrate	105 µmol/l
R2 standard	See bottle label

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

SPECIMEN COLLECTION AND HANDLING

Use serum or plasma (heparin) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability:	7 days	at 20–25 °C
	7 days	at 4–8 °C
	at least 1 year	at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with the standard included in the kit or with the calibrator XL MULTI-CAL 4x3 (Cat. No. XSYS0034) or XL MULTICAL 10x3 (Cat. No. XSYS0122) is recommended.

QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat. No. BLT00080) or ERBA NORM 10x5 (Cat. No. XSYS0123) and ERBA PATH 4x5 (Cat. No. BLT00081) or ERBA PATH 10x5 (Cat. No. XSYS0124) are recommended.

UNIT CONVERSION

mmol/l = mEq/l

EXPECTED VALUES¹

Serum:

Adult	98–107 mmol/l
Cord	96–104 mmol/l
Premature	95–110 mmol/l
0–30 d	98–113 mmol/l

Urine:

Adult	110–250 mmol/24 hours
-------	-----------------------

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	3.7 mmol/l
Linearity:	160 mmol/l
Measuring range:	3.7–160 mmol/l

Intra-assay precision Within run (n=20)	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	116.1	1.25	1.08
Sample 2	108.3	1.64	1.52

Inter-assay precision Run to run (n=20)	Mean (mmol/l)	SD (mmol/l)	CV (%)
Sample 1	117	1.47	2.75
Sample 2	94.6	1.43	1.52

COMPARISON

A comparison with a commercially available test (x) using 40 samples gave following results:

$$y = 1.000 x + 1.0 \text{ mmol/l}$$

$$r = 0.927$$

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 500 mg/dl.

Note:

Grossly haemolysed, lipaemic and jaundiced samples give falsely elevated results and should not be used. High levels of immunoglobulins, as found in case of multiple myeloma will interfere due to the development of turbidity.

Do not use normal saline for diluting samples.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

UFI: 058F-WERG-3H7T-HPU2

Danger

Contains: methanol, Mercury dithiocyanate, Mercury dinitrate.

Hazard statement:

H302 Harmful if swallowed.

H314 Causes severe skin burns and eye damage.

H371 May cause damage to organs.

Precautionary statement:

P260 Do not breathe vapours.

P280 Wear protective gloves/protective clothing/eye protection.

P301+P330+P331 IF SWALLOWED: Rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P308+P313 IF exposed or concerned: Get medical advice/attention.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength	492 (470–630) nm
Cuvette	1 cm

	Reagent blank	Standard (Cal.)	Sample
Reagent 1	1000 µl	1000 µl	1000 µl
Sample	–	–	10 µl
Standard (Cal.)	–	10 µl	–
Distilled water	10 µl	–	–

Mix and incubate 5 min. at 37 °C. Measure absorbance at of the sample A_{sam} and standard A_{st} against reagent blank.

CALCULATION

$$\text{Chloride (mmol/l)} = \frac{\Delta A_{\text{sam}}}{\Delta A_{\text{st}}} \times C_{\text{st}} \quad C_{\text{st}} = \text{standard (calibrator) concentration}$$

Applications for automatic analysers are available on request.

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