LDL DIRECT

Cat. No.	Pack Name	Packaging (Content)
BLT00041	LDL 80	R1: 2 x 30 ml, R2: 2 x 10 ml

EN

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of LDL Cholesterol in human serum and plasma.

CLINICAL SIGNIFICANCE

Low Density Lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on triglyceride-rich Very Low Density Lipoproteins (VLDLs). Specific LDL receptors exist to facilitate the elimination of LDL from plasma by liver parenchymal cells. It has been shown that most of the cholesterol stored in atherosclerotic plaques originates from LDL. For this reason the LDL Cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis. ^{2,4}

Accurate measurement of LDL Cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture. Can be applied on automated analyzers.

PRINCIPLE

The assay is based on a modified polyvinyl sulfonic acid (PVS) and polyethylene-glycol methyl ether (PEGME) coupled classic precipitation method with the improvements in using optimized quantities of PVS/PEGME and selected detergents. LDL, VLDL, and chylomicron (CM) react with PVS and PEGME and the reaction results in inaccessibility of LDL, VLDL and CM by cholesterol oxidase (CHOD) and cholesterol esterase (CHER), whereas HDL reacts with the enzymes. Addition of R2 containing a specific detergent releases LDL from the PVS/PEGME complex. The released LDL reacts with the enzymes to produce H₂O₂ which is quantified by the Trinder reaction.

HDL + LDL + VLDL + CM $\xrightarrow{\text{PVS / PEGME}}$ HDL + (LDL + VLDL + CM) • PVS / PEGME

HDL $\xrightarrow{\text{CHOD, CHER}}$ Fatty Acid + H₂O₂

(LDL + VLDL + CM) • PVS / PEGME $\xrightarrow{\text{Detergent}}$ LDL + (VLDL + CM) • PVS / PEGME

LDL $\xrightarrow{\text{CHOD, CHER}}$ Fatty Acid + H₂O₂

2H₂O₂ + 4-AA+TODB $\xrightarrow{\text{Peroxidase}}$ Quinone + 5 H₂O

REAGENT COMPOSITION

R1

MES buffer (pH 6.5)	50 mmol/l
Polyvinylsulfonic acid	50 mg/l
Polyethyleneglycolmethylester	30 ml/l
4-aminoantipyrine	0.9 g/l
Cholesterol esterase	5 kU/l
Cholesterol oxidase	20 kU/l
Peroxidase	5 kU/l
Detergent	

R2

MES buffer (pH 6.5)	50 mmol/l
Detergent	
TODB N,N-Bis(4-sulfobutyl)-3-methylaniline	3 mmol/l

REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. Once opened both reagents R1 & R2 are stable for 60 days at 2–8 °C, when protected from contamination.

Reagents are light-sensitive. Do not let bottles remain open. Keep containers tightly closed.

SPECIMEN COLLECTION AND HANDLING

Use serum or heparin plasma.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability in serum/plasma:	12 hours	at 20–25 °C
	10 days	at 4–8 °C
	12 weeks	at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with HDL/LDL CAL, Cat. No. XSYS0061 is recommended.

QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5, Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

UNIT CONVERSION

mg/dl x 0.026 = mmol/l

EXPECTED VALUES ¹¹

Less than 100 mg/dl	– optimal
100–129 mg/dl	– near/above optimal
130–159 mg/dl	– borderline high
160–189 mg/dl	– high
≥ 190 mg/dl	– very high

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Limit of quantification:	2.60 mg/dl
Linearity:	263 mg/dl
Measuring range:	2.60–263 mg/dl

PRECISION

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	50.346	0.885	1.72
Sample 2	82.308	1.808	2.21

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	47.00	0.885	1.91
Sample 2	92.69	1.500	1.61

COMPARISON

A comparison between LDL DIRECT (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.964x -1.615 mg/dl
r = 0.995

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl. Interference by N-acetylcysteine (NAC), acetoaminophen and metamizole causes falsely low results. To carry out the test, blood withdrawal should be performed prior to administration of drugs.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength:	600/700 nm
Cuvette:	1 cm

Pipette in Tube	Reagent blank	Sample / Calibrator
Reagent 1	375 µl	375 µl
D.D water	3 µl	-
Sample / Calibrator	-	3 µl
Mix and incubate at 37°C for 5 min.		
Add Reagent 2	125 µl	125 µl
Mix and incubate at 37 °C for 5 min.		

Read final absorbances at the specified wavelength against reagent blank.

CALCULATION

LDL C = $\frac{(\text{Abs. of Sample} - \text{Abs. of Sample Blank})}{(\text{Abs. of Cal.} - \text{Abs. of Cal. Blank})} \times \text{Concentration of Calibrator}$

ASSAY PARAMETERS

Mode	1-Point End
Wavelength (Primary)	600 nm
Wavelength (Secondary)	700 nm
Sample Volume	3 µl
Reagent 1 Volume	375 µl
Reagent 2 Volume	125 µl
Incubation time	5 min.
Incubation temperature	37 °C
Normal low	0
Normal high	130
Linearity low	2.6
Linearity high	263
Blank with	Reagent
Absorbance limit (max.)	0.3
Units	mg/dl

Program parameters for specific clinical analyzers are available on request.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
e-mail: diagnostics@erba.com, www.erbamannheim.com

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1. "Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult treatment Panel III)", JAMA, 285:2486 (2001).
2. Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26; 566 (1985).
3. Castelli WP et al., Cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
4. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med., 11; 480 (1951).
5. Badimon JJ, Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85:1234-41(1990).
6. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am.J. Med., 62; 707 (1977).
6. Kannel WB, Castelli WP, Gordon T., Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Am. Intern. Med., 90; 85 (1979).
7. William P., Robinson D., Baily A., High density lipoprotein and coronary risk factor, Lancet, 1; 72 (1979).
8. Castelli, W. P., et al, Cholesterol and other lipids in coronary heart disease.
9. Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of the High Blood Cholesterol in Adults (Adult Treatment Panel III).
10. Pisani T, GebSKI CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 1995; 119:1127)
11. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.

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