MAGNESIUM

Cat. No.	Pack Name	Packaging (Content)
BLT00049	MG 250	R1: 1 x 250 ml, R2 STD: 1 x 5 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Magnesium in human serum, plasma, cerebrospinal fluid and urine.

CLINICAL SIGNIFICANCE

Magnesium is an essential nutrient which is involved in many biochemical functions. It has a structural role in nucleic acids and ribosomal particles, required as an activator for many enzymes and has a role in energy producing oxidative phosphorylation.

Hypomagnesaemia results in the impairment of neuromuscular functions and may develop in severe prolonged diarrhoea, malabsorption syndromes, primary aldosteronism and diuretic therapy. Hypermagnesaemia is seen in renal glomerular failure and diabetic coma.

PRINCIPLE

Magnesium reacts with Xylidyl Blue to form a colored compound in alkaline solution. The intensity of the colour formed is proportional to the magnesium concentration in the sample. Interference with calcium is prevented by the use of GEDTA.

REAGENT COMPOSITION

R1	
Xylidyl Blue (I)	110 µmol/l
Ethanolamine	1 mol/l
GEDTA	60 µmol/l
R2 standard	See bottle label

REAGENT PREPARATION

Reagent is liquid, ready for use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8°C.

SPECIMEN COLLECTION & HANDLING

Use serum, plasma, cerebrospinal fluid (CSF) or urine. Do not use EDTA plasma. It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum/plasma:	7 days	at 20–25°C
	7 days	at 4–8°C
	1 year	at -20°C
in urine:	3 days	at 20–25°C
	3 days	at 4–8°C
	1 year	at -20°C

Discard contaminated specimens.

Acidify urine with some drops of conc. HCl to pH 3–4, then dilute 1+4 with dist. water; multiply the result by 5.

CALIBRATION

Calibration with the standard included in the kit or calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

UNIT CONVERSION

mg/dl x 0.4114 = mmol/l

EXPECTED VALUES⁶

Serum or plasma	
Newborn, 2 - 4 d	1.5 – 2.2 mg/dl
5 mo – 6 y	1.7 – 2.3 mg/dl
6 – 12 y	1.7 – 2.1 mg/dl
> 12 y	1.6 – 2.6 mg/dl
CSF	2.1 – 3.3 mg/dl
Urine	12 – 291 mg / 24 h

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	0.16 mg/dl
Linearity:	5.85 mg/dl
Measuring range:	0.16 - 5.85 mg/dl

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	3.305	0.122	3.66
Sample 2	5.854	0.120	2.04

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	2.341	0.039	1.65
Sample 2	4.634	0.095	2.05

COMPARISON

A comparison between XL-Systems Magnesium (y) and a commercially available test (x) using 40 samples gave following results:

$$y = 1.014 x - 0.005 \text{ mg/dl}$$

$$r = 0.995$$

INTERFERENCES

Following substances do not interfere:
 bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl.
 Haemoglobin interferes because magnesium is released by erythrocytes.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 contains 2-aminoethanol.

Danger

Hazard statement(s):

H315 Causes skin irritation.

H318 Causes serious eye damage.

Precautionary statement(s):

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water and soap.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength: 520 (500 – 550) nm

Cuvette: 1 cm

Temperature: 37°C

Serum / reaction mixture ratio 1/101

	Reagent blank	Standard (Calibrator)	Sample
Reagent 1	1.00 ml	1.00 ml	1.00 ml
Sample	-	-	0.01 ml
Standard (calibrator)	-	0.01 ml	-
Distilled water	0.01 ml	-	-

Mix and incubate 5 min. at 37 °C. Measure absorbance of the sample A_{sam} and standard A_{st} against reagent blank. The coloration is stable during one hour.

CALCULATION

$$\text{Magnesium (mg/dl)} = \frac{\Delta A_{\text{sam}}}{\Delta A_{\text{st}}} \times C_{\text{st}}$$

C_{st} = standard (calibrator) concentration

Applications for automatic analysers are available on request.

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	End Point
Wavelength 1 (nm)	505
Sample Volume (µl)	5/10
Reagent Volume (µl)	500/1000
Incubation time (min.)	5
Incubation temp. (°C)	37
Normal Low (mg/dl)	1.6
Normal High (mg/dl)	2.6
Linearity Low (mg/dl)	0.16
Linearity High (mg/dl)	5.85
Concentration of standard	See bottle label
Blank with	Reagent blank
Absorbance limit (max.)	2.0
Units	mg/dl

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ
 e-mail: diagnostics@erbamannheim.com, www.erbamannheim.com

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5. Sitzmann FC. Normalwerte. München: Hans Marseille Verlag GmbH; 1986. p. 166.