
AST/GOT

Cat. No.	Pack Name	Packaging (Content)
BLT00050	AST/GOT 250	R1: 4 x 50 ml, R2: 1 x 50 ml
BLT00051	AST/GOT 500	R1: 4 x 100 ml, R2: 1 x 100 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of AST/GOT (Aspartate Amino-transferase) in human serum and plasma.

CLINICAL SIGNIFICANCE

AST/GOT occurs in all human tissues and is present in large amounts in liver, renal, cardiac and skeletal muscle tissue.

Increased levels are associated with liver diseases or damage myocardial infarction, muscular dystrophy and cholecystitis.

Decreased levels are observed in patients undergoing renal dialysis and those with B6 deficiency. Monitoring the change in levels over a period of time is beneficial to the physician evaluating myocardial infarction or following chronic or resolving hepatitis.

PRINCIPLE

International Federation of Clinical Chemistry (IFCC), without Pyridoxal Phosphate.

AST: Aspartate aminotransferase

LDH: Lactate dehydrogenase

MDH: Malate dehydrogenase

The rate of absorbance change at 340 nm is directly proportional to AST/GOT activity in the specimen.

REAGENT COMPOSITION

R1

Tris Buffer (pH 7.8) 110 mmol/l

L-Aspartate 340 mmol/l

LDH ≥ 4000 U/l

MDH ≥ 750 U/l

R2

CAPSO 20 mmol/l

2-Oxoglutarate 85 mmol/l

NADH 1.05 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

Two reagents method – substrate start

Reagents R1 and R2 are liquid, ready to use. After the first opening the vials, reagents are stable for 30 days at 2–8 °C in the dark.

Monoreagent method – sample start

Mix 4 portion of reagent R1 with 1 portion of reagent R2.

Stability: 5 days at 20–25 °C in the dark

4 weeks at 2–8 °C in the dark

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (EDTA, heparin).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Loss of activity: at 2–8 °C < 8 % within 3 days
at 15–25 °C < 10 % within 3 days

Stability at -20 °C at least 3 months.

Discard contaminated specimens.

CALIBRATION

Calibration with the calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

UNIT CONVERSION

U/l x 0.017 = μ kat/l

EXPECTED VALUES ⁴

At 37 °C Men up to 35 U/l

Women up to 31 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 3.84 U/l

Linearity: 390 U/l

Measuring range: 3.84 – 390 U/l

PRECISION

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	103.2	0.60	0.54
Sample 2	313.2	1.68	0.54

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	43.8	0.60	1.37
Sample 2	115.2	1.08	0.92

COMPARISON

A comparison between XL-Systems AST/GOT (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.941 x - 3.96 U/l

r = 0.996

INTERFERENCES

Following substances do not interfere:

bilirubin up to 30 mg/dl, triglycerides up to 2000 mg/dl, haemolysis interferes due to ASAT activity from erythrocytes.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent 1 of the kit contains Asparatic acid 4,5 %, Tris(hydroxymethyl)aminomethane 1,3 %, hidróxido de sodio < 0,5 %.

Warning

Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water and soap.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+P313 If eye irritation persists: Get medical advice/attention.

Reagent 2 is not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

ASSAY PROCEDURE

Wavelength 340 nm, Hg 334 nm, Hg 365 nm

Cuvette 1 cm

Two reagents method – substrate start

Reagent 1 (buffer)	0.800 ml
Sample	0.100 ml

Mix and incubate for 5 min. at 37 °C. Then add:

Reagent 2 (substrate)	0.200 ml
-----------------------	----------

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA /min).

Monoreagent method – sample start

Working solution	1.000 ml
Sample	0.100 ml

Mix, incubate 1 min. at 37 °C and then measure the initial absorbance of calibrator and sample against reagent blank. Measure the absorbance change exactly after 1, 2 and 3 min. Calculate 1 minute absorbance change (ΔA /min).

CALCULATION

$$1. \quad \text{AST/GOT (U/l)} = C_{\text{cal}} \times \frac{\Delta A_{\text{sam}}/\text{min}}{\Delta A_{\text{cal}}/\text{min}} \quad C_{\text{cal}} = \text{calibrator concentration}$$

$$2. \quad \text{Using factor:} \quad \text{AST/GOT} = f \times \Delta A/\text{min} \quad f = \text{factor}$$

Factors:	Substrate Start:	25° or 30°C	37°C
	Factor at 340 nm	1151	2143
	Factor at 334 nm	1173	2184
	Factor at 365 nm	2132	3971
	Sample Start:	25° or 30°C	37°C
	Factor at 340 nm	952	1745
	Factor at 334 nm	971	1780
	Factor at 365 nm	1765	3235

ASSAY PARAMETERS FOR PHOTOMETERS

Mode	Kinetic
Wavelength (nm)	340
Sample Volume (μl)	50 / 100
Reagent Volume (μl)	500 / 1000
Lag time (sec.)	60
Kinetic interval (sec.)	60
No. of readings	3
Kinetic factor	1745
Reaction temperature (°C)	37
Reaction direction	Decreasing
Normal low	0
Normal high (U/l)	31
Linearity low (U/l)	3.84
Linearity high (U/l)	390
Absorbance limit (min.)	1.1
Blank with	Water
Units	U/l

Program parameters for specific clinical analyzers are available on request.

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N/35/22/H/INT

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REFERENCES / ЛИТЕРАТУРА / ЛІТЕРАТУРА / LITERATURA / LITERATÚRA / REFERENCIAS

1. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 55-65.
2. Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.
3. Schumann G, Bonora R, Ceriotti F, Férard G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002;40:725-33.
4. Tietz Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.

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