URIC ACID AOD

Cat. No.	Pack Name	Packaging (Content)
BLT00063	UA 500	R1: 4 x 100 ml, R2: 1 x 100 ml
BLT00064	UA 500 S	R1: 4 x 100 ml, R2: 1 x 100 ml, R3 STD: 1 x10 ml Only for CIS

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Uric Acid in human serum, plasma and urine.

CLINICAL SIGNIFICANCE

Uric acid is a metabolite of purines, nucleic acids and nucleoproteins, consequently, abnormal levels may be indicative of a disorder in the metabolism of these substances. Hyperuricaemia may be observed in renal dysfunction, gout, leukemia, polycythaemia, atherosclerosis, diabetes and hypothyroidism. Decreased levels are present in patients with Wilson's Disease.

PRINCIPLE

Uric acid is oxidized by oxygen to hydrogen peroxide and allantoin. This oxidation is catalysed by uricase. The produced hydrogen peroxide is determined by oxidative coupling with a natrium salt of N-ethyl-N-(2-hydroxy-3-sulphopropyl)-m-toluidine and 4-aminoanti-pyrene. The reaction is catalysed by peroxidase.

REAGENT COMPOSITION

R1
Phosphate buffer, pH 7.8 (25 °C) 0.1 mol/l
TOOS 0.626 mmol/l
AOD ≥ 20.0 µkat/l

R2
Phosphate buffer, pH 7.8 (25 °C) 0.1 mol/l
POD ≥ 100.0 µkat/l
Uricase ≥ 10.0 µkat/l
4-AAP 2.5 mmol/l

R3 standard
Uric acid 357 µmol/l

COMPOSITION OF REACTION MIXTURE

Phosphate buffer, pH 7.8 (25 °C) 0.098 mol/l
N-ethyl-N-(2-hydroxy-3-sulphopropyl)-
-m-toluidine, natrium salt (TOOS) 0.49 mmol/l
4-aminoantipyrene (4-AAP) 0.49 mmol/l
Uricase ≥ 1.96 µkat/l
Peroxidase (POD) ≥ 19.60 µkat/l
Ascorbate oxidase (AOD) ≥ 15.80 µkat/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

Two - reagent method

Reagents R1, R2 and R3 are liquid, ready to use.

If stored at 2–8°C, reagents are stable until expiry date, that is stated on the package.

After first opening, reagents are stable for 90 days at 2–8°C if stored at appropriate conditions, closed carefully and without any contamination.

Mono - reagent method

Mix 4 portions of reagent R1 with 1 portion of reagent R2.

Stability: 3 weeks at 2–8 °C in the dark
5 days at 15–25 °C in the dark

SPECIMEN COLLECTION & HANDLING

Use serum or plasma (heparin, EDTA) or urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum / plasma:

3 days at 20–25°C
7 days at 4–8 °C
6 months at -20°C

in urine:

4 days at 20–25°C

For the determination in urine use 24 hours specimen. To prevent the precipitation of uric acid add 15 ml 5 mol/l NaOH into the urine collector to ensure urine pH > 8. Dilute urine samples in 1+9 ratio with distilled water and multiply results by 10.

Discard contaminated specimens.

CALIBRATION

For the calibration it is recommended to use Lyonorm Calibrator or the standard included in the set.

QUALITY CONTROL

For quality control, it is recommended to use Lyonorm HUM N, Lyonorm HUM P.

UNIT CONVERSION

mg/dl x 59.48 = µmol/l

EXPECTED VALUES ²

fS uric acid (µmol/l)

men 220–420

female 140–340

dU uric acid (mmol/24 hours) 1.5–4.5

The range of reference values is only approximate; it is recommended that each laboratory verify the extent of the reference interval for their particular examined population.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 4.5 µmol/l

Linearity: 1500 µmol/l

Measuring range: 4.5–1500 µmol/l

Intra-assay precision Within run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	269	7.7	2.86
Sample 2	658	3.8	1.39

Inter-assay precision Run to run (n=20)	Mean (µmol/l)	SD (µmol/l)	CV (%)
Sample 1	282	11.5	4.07
Sample 2	669	17.2	2.58

COMPARISON

A comparison between UA 500/UA 500 S (y) and a commercially available test (x) using 40 samples gave following results:

r = 0.999

y = 1.015x + 2.4000 µmol/l

INTERFERENCES

Following substances do not interfere:

bilirubin up to 15 mg/dl, hemoglobin up to 5 g/l, triglycerides up to 1500 mg/dl. Interference by N-acetylcysteine (NAC), acetoaminophen and metamizole causes falsely low results. To carry out the test, blood withdrawal should be performed prior to administration of drugs.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the kit are not classified like dangerous.

EUH 210 Safety data sheet available on request.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

PROCEDURE

Wavelength: 550 nm

Cuvette: 1 cm

Temperature: +37 °C

Serum/reaction mixture ratio 1/51

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

Two-reagent method

	Reagent blank	Calibrator (Standard)	Sample
Reagent R1	0.80 ml	0.80 ml	0.80 ml
Sample	–	–	0.02 ml
Calibrator (Standard)	–	0.02 ml	–
Distilled water	0.02 ml	–	–

Mix and after 1–5 min. incubation read the initial absorbance for blank A_{bi} , sample A_{sam} and standard (calibrator) A_{st} . Then add:

Reagent R2	0.2 ml	0.2 ml	0.2 ml
------------	--------	--------	--------

Mix and after 1–2 min. incubation read the absorbance for blank A_{bi} , sample A_{sam} and standard (calibrator) A_{st} .

Calculate resulting absorbance like the difference between the final and initial absorbance $A = (A_{FINAL} - A_{INITIAL})$.

Mono-reagent method

	Reagent blank	Calibrator (Standard)	Sample
Working reagent	1.0 ml	1.0 ml	1.0 ml
Sample	–	–	0.02 ml
Calibrator (Standard)	–	0.02 ml	–
Distilled water	0.02 ml	–	–

Mix and after 1–2 min. incubation read the absorbance for blank A_{bi} , sample A_{sam} and standard (calibrator) A_{st} .

CALCULATION

$$\text{Uric acid } (\mu\text{mol/l}) = \frac{\Delta A_{sam} - \Delta A_{bi}}{\Delta A_{st} - \Delta A_{bi}} \times C_{st}$$

C_{st} = standard (calibrator) concentration)

NOTE

For the uric acid determination it is recommended to use two reagents procedure (with sample blank) due to an interference-suppression of ascorbic acid and visual sample defects.

Applications for automatic analysers will be supplies on request.

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