CREATININE ENZYMATIC

Cat. No.	Pack Name	Packaging (Content)
BLT00065	CREA ENZ 204	R1: 3 x 50 ml, R2: 3 x 18 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Creatinine in human serum, plasma and urine.

CLINICAL SIGNIFICANCE

Creatinine is a waste product formed in muscle from the high energy storage compound, creatine phosphate. The amount of creatinine produced is fairly constant (unlike Urea) and is primarily a function of muscle mass. It is not greatly affected by diet, age, sex or exercise. Creatinine is removed from plasma by glomerular filtration and then excreted in urine without any appreciable resorption by the tubules. Creatinine is used to assess renal function, however, serum creatinine levels do not start to rise until renal function has decreased by at least 50%.

PRINCIPLE

In the first reaction, creatinase and sarcosine oxidase are used in the enzymatic hydrolysis of endogenous creatine to produce hydrogen peroxide, that is eliminated by catalase. Creatininase and 4-aminoantipyrine are added, and only the creatine generated from creatinine by creatininase is hydrolysed sequentially by creatinase and sarcosine oxidase to produce hydrogen peroxide. This newly formed hydrogen peroxide is measured in a coupled reaction catalysed by peroxidase, with N-ethyl-N-sulphopropyl-m-toluidine (ESPM) as a chromogen.

The absorbance of the produced complex at 546 nm is proportional to the creatinine concentration in the sample.

REAGENT COMPOSITION

R1	
MOPS pH 7.5	25 mmol/l
TOPS	0.5 mmol/l
Creatinase	10 kU/l
Sarcosine Oxidase	5 kU/l
Catalase	3 kU/l
EDTA	1 mmol/l
R2	
MOPS pH 7.5	90 mmol/l
Creatininase	30 kU/l
Peroxidase	10 kU/l
Sodium azide	0.5 g/l

REACTION MIXTURE

MOPS pH 7.5	41 mmol/l
TOPS	0.37 mmol/l
Creatinase	7.4 kU/l
Sarcosine oxidase	3.7 kU/l
Catalase	2.2 kU/l
EDTA	0.7 mmol/l
Creatininase	7.4 kU/l
Peroxidase	2.5 kU/l

REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

After first opening, reagents are stable for 30 days at 2–8 °C if stored at appropriate conditions, closed carefully and without any contamination.

SPECIMEN COLLECTION & HANDLING

Use serum, plasma (heparin, EDTA), urine.

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum / plasma:	7 days	at 4–25 °C
	at least 3 months	at -20 °C
in urine:	2 days	at 20–25 °C
	6 days	at 4–8 °C
	6 months	at -20 °C

For the determination in urine use 24 hours specimen. It is important to exactly measure the volume of collected urine. Dilute urine samples in 1+19 ratio with distilled water and multiply results by 20.

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator Lyonorm Calibrator is recommended.

QUALITY CONTROL

For quality control, it is recommended to use following materials: Lyonorm HUM N, Lyonorm HUM P.

UNIT CONVERSION

mg/dl x 88.4 = µmol/l

EXPECTED VALUES

fS Creatinine (µmol/l)	
male	44–97
female	44–80
dU Creatinine (mmol/24 hrs)	5–18

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification:	0.064 mg/dl
Linearity:	85.8 mg/dl
Measuring range:	0.064–85.8 mg/dl

Intra-assay (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	1.01	0.006	0.63
Sample 2	3.81	0.016	0.42

Inter-assay (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	1.02	0.013	1.26
Sample 2	3.76	0.065	1.72

COMPARISON

A comparison between XL-Systems CREA ENZ (y) and a commercially available test (x) using 50 samples gave following results:

N = 50
r = 0.999
y = 0.9595 x – 0.016 mg/dl

INTERFERENCES

Following substances do not interfere: haemoglobin up to 12.5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 850 mg/dl. Interference by N-acetylcysteine (NAC), acetoaminophen and metamizole causes falsely low results. To carry out the test, blood withdrawal should be performed prior to administration of drugs.

N-acetyl-p-benzoquinone imine (metabolite of paracetamol) could generate erroneously low results in samples for patients that have taken toxic doses of paracetamol.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Reagents are not classified as hazardous. R2 contains < 0,1% sodium azide.

WASTE MANAGEMENT

Please refer to local legal requirements.

FIRST AID

In case of an accidental ingestion, wash up the mouth and drink about 0.5 l of water. On eye contact rinse the eye quickly and thoroughly with the jet of tap of water. Contaminated skin should be washed with warm water and soap. In all serious cases of health damage consult a physician.

WASTE DISPOSAL

All tested samples should be treated as potentially infectious and with an eventual rest of reagents should be disposed in accordance with the internal regulations for dangerous waste, in compliance with local and national regulations relating to the safe handling of dangerous materials.

Paper packing and others should be handed over for recycling or discarded as sorted waste (paper, glass, plastic).

PROCEDURE

Wavelength:	546 (530–560) nm
Cuvette:	1 cm
Temperature:	37 °C
Serum/reaction mixture ratio	1/61

Reagents and sample volume can be modified, by respecting reagents/sample volume ratio.

Two-reagent method

	Reagent blank	Calibrator (Standard)	Sample
Reagent R1	0.90 ml	0.90 ml	0.90 ml
Sample	-	-	0.02 ml
Calibrator (Standard)	-	0.02 ml	-
Distilled water	0.02 ml	-	-

Mix and incubate 3-5 min. at 37 °C. Measure absorbance A1 of the sample (A_{sam}), calibrator (A_{cal}) and blank (A_{bl}). Then add:

Reagent R2	0.30 ml	0.30 ml	0.30 ml
------------	---------	---------	---------

Mix and incubate 5–10 min. at 37 °C. Measure absorbance A2 of the sample (A_{sam}), calibrator (A_{cal}) and blank (A_{bl}). Calculate absorbance change ΔA = A2 – A1 for sample, calibrator and blank.

CALCULATION

$$\text{Creatinine } (\mu\text{mol/l}) = \frac{\Delta A_{\text{sam}} - \Delta A_{\text{bl}}}{\Delta A_{\text{cal}} - \Delta A_{\text{bl}}} \times C_{\text{cal}}$$

C_{cal} = calibrator (standard) concentration

Applications for automatic analysers will be supplied on request.

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