LYONORM HUM P

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Characteristics

LYONORM HUM P is a lyophilized control serum based on human matrix with most constituent concentrations and activities in the pathological range. The added enzymes of animal origin have the properties as these enzymes in human serum.

Use

Control serum LYONORM HUM P is intended for the control of accuracy and precision of laboratory testing procedures for the analytes listed in the package insert.

Reconstitution of lyophilizate

Open a bottle very carefully avoiding any loss of the lyophilized material. Add exactly 5,00 ml of redistilled water. Close the bottle carefully, swirl gently and let stand for 30 min at (+15 to +25)°C. In this time invert the bottle several times until the solution is homogenous (do not shake, avoid the formation of the foam). Control serum can be used 30 min after its reconstitution. Only for Alkaline Phosphatase determination wait for another 30 min.

The reconstituted control serum is to be treated like the serum or blood specimen.

Stability and storage

When stored at a dark place, the LYONORM HUM P is stable in the lyophilized form at temperature of (+2 to +8) °C till the date of expiration stated on the vial. The expiry date for the lyophilised control serum is given on each vial.

The constituents when protected from light and contamination are stable:

- 1 day at 18–25°C
- 1 week at 2–8°C (bilirubin: 1 day)
- 1 month at –20 °C (bilirubin: 2 weeks)

Health protection

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagents of the set contain biological material, which is potentially infectious. This control serum has been prepared exclusively from the blood of donors tested individually and shown to be free from antibodies to HIV (human immunodeficiency virus also known as HTLV III/LAV) and HB_s Ag. Testing for HIV antibodies in serum or in plasma as well as in the products of them is less reliable. Because no known method can eliminate the risk of infection the reconstituted control serum is to be treated like potentially infectious specimen. During the work all principles of personal hygiene are to be maintained. Do not eat, drink, smoke. Use the personal protection aids. After finishing the work wash your hands with soap and warm water.

When the bottle is broken the contaminated place must be disinfected with Persteril or Chloramine.

Waste disposal

The set should be treated as potentially infectious and the contingent rest of the reagent should be liquidated in accordance with any other local and national regulations relating to the safe handling of such materials. Put packaging paper waste to recycling or in municipal waste.

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