HDL CHOLESTEROL CALIBRATOR

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Characteristics

HDL CHOLESTEROL CALIBRATOR is lyophilized human serum. It is a monospecific calibrator. The value of the calibrator was determined by the procedure recommended by the National Reference System for cholesterol (NRS/CHOL) and complies with significant medical values.

Intended use

The HDL CHOLESTEROL CALIBRATOR is intended for the quantitative determination of HDL cholesterol by the HDL C DIRECT 240 (Cat. No. BLT00030) test in serum.

Batch	0000000
Expiration	MM/RRRR

Certified values of calibrator

Parameter	Concentration	Unit
HDL cholesterol	0.00	mmol/l

Unit conversion:

mg/dl x 0.026 = mmol/l

Characteristics

HDL CHOLESTEROL CALIBRATOR is lyophilized human serum. It is a monospecific calibrator. The value of the calibrator was determined by the procedure recommended by the National Reference System for cholesterol (NRS/CHOL) and complies with significant medical values.

Reagents

HDL CHOLESTEROL CALIBRATOR (1 bottle for 3 ml)
It contains lyophilized human serum.

Reagent preparation

Open the bottle very carefully and add exactly 3.00 ml of redistilled or deionized water.

Reconstitute the calibrator within 30 min at room temperature of (+15 to +25) °C in a dark place. Swirl the bottle carefully, do not shake!

Stability

When stored at a dark place, the HDL CHOLESTEROL CALIBRATOR is stable in the lyophilized form:

at temperature of (+2 to +8) °C till the date of expiration stated on the bottle.

The reconstituted calibrator is stable for at least:

32 h at (+15 to +25) °C in dark

15 days at (+2 to +8) °C in dark

1 month at (-20) °C in dark, when frozen once

Data concerning calibrator stability are valid only for the calibrator not contaminated microbially.

Do not use the calibrator with a turbidity developed pointing at a possible bacterial growth.

Warning

The HDL CHOLESTEROL CALIBRATOR is a diagnosticum for *in vitro* use.

Not to be used internally in humans or animals.

Do not use this product past the expiration date stated on the bottle.

Do not use this product for any purpose other than described herein.

Health protection

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

The control sera have been prepared from normal human serum of donors tested to be free of antibodies of HIV (human immunodeficiency virus), HCV (hepatitis C virus) and HBsAg by specific procedures (ELISA and RIA).

The negative results of tests for HIV antibodies in the case of lyophilized product were confirmed by immunofluorescent method.

However, since the risk of infection can never be excluded, the calibrator is to be handled like a patient specimen.

During the work all principles of general safety as well as the principles of work with biological materials are to be maintained. Above all observe the principles of personal hygiene, use personal protection aids (clothes, shoes, etc.), do not eat, drink and smoke. After finishing the work and before a meal wash your hands with soap and warm water and apply a suitable regenerating cream.

When the bottle is broken, the contaminated surface must be disinfected with Persteril or Chloramin.

Reagent contains < 0.025 % reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one.

Warning

Hazard statement:

H317 May cause an allergic skin reaction.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

First aid

In case of an accidental ingestion wash up the mouth and drink about 0.5 litre of water.

On eye contact rinse quickly and thoroughly with a jet of tap water.

Contaminated skin wash with warm water and soap.

In all serious cases of health damage consult a physician.

Waste disposal

The kit should be held as a potentially infectious material. The eventual rest of the reagent must be disposed according to inner regulations for dangerous waste, in compliance with local and national regulations. Paper packing are handed over for recycling or discarded as communal waste.

Date of revision: 8. 7. 2020

Used Symbols

	Catalogue number		Lot number		In vitro diagnostics		Expiry date		Manufacturer		Storage temperature		See instruction for use
---	------------------	---	------------	---	----------------------	---	-------------	---	--------------	---	---------------------	---	-------------------------

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