ERBA LIPID NORM

CE

IVD

EN

CONTROL SERUM FOR BIOCHEMISTRY

Assayed Values for ERBA Reagents in NORMAL range

Cat. No.:	Pack name:	Packaging (Content):
BLT00074	ERBA LIPID NORM	3 × 5 mL

INTENDED USE
Control serum (assayed) is intended especially for control of accuracy and precision of the determination of HDL and LDL cholesterol, cholesterol, triglycerides and apolipoprotein A1 and apolipoprotein B in normal or slightly increased range.

PRINCIPLE OF THE PROCEDURE
This human serum control is to be used in the same way as an unknown specimen that would be used in accordance with the direction of the assay procedure. The results obtained for the control are to be compared with the assigned values given in the "assay data" section of the insert and an evaluation made by standard statistical techniques to determine if the procedure is within the control limits.

COMPOSITION
The control serum is prepared from human serum with chemical additives and tissue extracts of human and biological origin. Bacteriostatic agents have been added. The control is provided in lyophilized form for increase stability.

WARNINGS
For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. The control sera have been prepared from normal human serum, which was tested for HIV and HBsAg antibodies by ELISA or RIA methods with negative results. The negative results of tests for HIV antibodies in the case of lyophilized products were confirmed by immunofluorescent method performed by the authorized institution. However, since the risk of infection can never be excluded, the product is to be handled like potentially infectious human material, similar to sera of patients.

During the work all principles of personal hygiene are to be maintained. Do not eat, drink, smoke. Use the personal protection aids. After finishing the work wash your hands with soap and warm water. When the bottle is broken the contaminated place must be disinfected.

RECONSTITUTION
Open the bottle very carefully avoiding any loss of lyophilized material. Add exactly 5.00 ml of redistilled water. Close the bottle tightly, swirl gently and allow stand for 30 min in the dark. Let the content to dissolve completely avoiding any formation of the foam. Do not shake the bottle.

STORAGE & STABILITY
Prior to reconstitution
The control should be stored at 2–8 °C and is stable till the expiry date printed on the label. Protect from light.
After reconstitution
The constituents when protected from light and contamination are stable:
1 day at 15–20 °C
1 week at 2–8 °C
1 month at -20 °C

LIMITATIONS
The results obtained using the control are dependent upon several factors. Erroneous results can occur from reconstitution inaccuracy and the technique errors associated with the assay procedure. Improper storage or handling of the control can also affect the results. If there is a visible evidence of microbial growth in a vial, do not use that vial.

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