CRP CONTROL

Lot: 0000000

| Cat. No. | Pack Name | Packaging (Content) |
|----------|-----------|---------------------|
| BLT20013 | CRP CON H | 1 x 1 ml |

INTENDED USE

Accuracy control for the determination of C-Reactive Protein (CRP) in serum by turbidimetry and nephelometry.

COMPOSITION

A dilution of human sera containing high levels of CRP with phosphate buffered saline. Liquid stabilised.

Preservative: 0.095 % sodium azide.

Ready for use.

STORAGE AND STABILITY

The expiry date of the product at 2–8 °C is listed on the label.

After first opening the container, the serum can be used for 6 weeks if stored tightly closed at 2–8 °C after use.

Do not freeze.

PRECAUTIONS AND WARNING

- For *in Vitro* Diagnostic use.
- Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV-1 and HIV-2 by FDA required test. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent is not classified as dangerous.

ASSIGNED VALUES

The value for CRP has been assigned to ERM-DA474/IFCC.

CRP:

| Value | Interval | 1 SD | Unit |
|-------|-----------|------|-------|
| x.xx | x.xx–x.xx | x.xx | mg/dl |
| xx.x | xx.x–xx.x | x.xx | mg/l |

Assigned values for XL-Erba analysers with calibration using automatic dilution – see the end of the document.

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