FERRITIN CALIBRATOR SET

Lot: xxxxxxxx

Expiry Date: YYYY-MM

Cat. No.	Pack name	Packaging
BLT20016	FERR CAL SET	5 x 1 ml

INTENDED USE

Preparation of reference curves for quantitative determination of Ferritin in human serum by turbidimetric immunoassay. The Calibrator Set has to be used with the FERRITIN kit, Cat. No. BLT20015.

COMPOSITION

Ferritin Calibrator Set are dilutions of human plasma, liquid stabilised and filtered through 0.2 µm.

Preservative: Sodium azide (< 0.1 %).

Ready to use.

STORAGE AND STABILITY

After first opening the bottle the calibrator can be used for at least 6 weeks if stored tightly closed at 2–8 °C.

WARNING AND PRECAUTIONS

1. For *in Vitro* Diagnostic use only
2. Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV-1 and HIV-2 by CE certified tests. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.
3. Reagents containing sodium azide must be handled with caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent is not classified as dangerous.

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