MAL CONTROL

Lot: xxxxxxxx

Expiry Date: YYYY-MM

Cat. No.	Pack Name	Packaging (Content)
BLT20033	MAL CON	1 × 1 ml

INTENDED USE

Accuracy control for quantitative immunochemical determination of Microalbumin in urine by turbidimetry and nephelometry.

COMPOSITION

A dilution of defibrinated human plasma with phosphate buffered saline, liquid stabilised and filtered through 0.2 µm.
Preservative: 0.095 % sodium azide.
Ready for use.

STORAGE AND STABILITY

After first opening the container, the control can be used for 6 weeks if stored tightly closed at 2–8 °C after use.
Do not freeze.

PRECAUTIONS AND WARNING

- For *in Vitro* Diagnostic use.
- Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV-1 and HIV-2 by FDA required test. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent is not classified as dangerous.

ASSIGNED VALUES

The value for microalbumin has been assigned to ERM-DA470k/IFCC.

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