ALKALINE PHOSPHATASE

Cat. No.	Pack Name	Packaging (Content)
XSYS0002	ALP 110	R1: 2 x 44 ml, R2: 2 x 11 ml

INTENDED USE

This reagent is intended for *in vitro* quantitative determination of ALP in human serum or plasma.


CLINICAL SIGNIFICANCE

Human ALP consists of a group of enzymes which hydrolyse phosphates at an alkaline pH. ALP is found in practically all tissues of the body but in high concentrations in the osteoblasts of bone, liver, placenta, kidney, intestinal wall and lactating mammary glands. In adults the ALP normally found circulating in the serum is largely derived from the liver. In children or in adolescents going through pubertal growth spurts, there is an additional contribution from bone and this accounts for the higher reference interval for these groups. Pregnancy also raises the normal values of ALP.

Raised ALP levels are often observed in bone disease or liver disease involving the biliary tract. If the source of the isoenzyme is not apparent then estimation of GGT may help differentiate between the two. A raised GGT in the presence of a raised ALP would suggest the liver is the primary source. Increased ALP (usually normal GGT) is seen in Osteomalacia and Rickets, primary hyperparathyroidism with bone involvement, Pagets disease, secondary carcinoma in bone and some cases of osteogenic sarcoma. Increased levels of ALP (usually with a raised GGT) is seen in cholestasis, hepatitis, cirrhosis, space occupying lesions and malignancy with bone or liver involvement or direct production. Low levels of ALP may be observed in conditions which cause arrested bone growth or in hypophosphatasia.

PRINCIPLE

The method according to IFCC recommendation. This method utilises 4-nitrophenyl phosphate as the substrate. Under optimised conditions ALP present in the sample catalyses the following reaction.

At the pH of the reaction, 4-nitrophenol has an intense yellow colour. The reagent also contains a metal ion buffer system to ensure that optimal concentrations of Zinc and Magnesium are maintained. The metal ion buffer can also chelate other potentially inhibitory ions which may be present. The reaction is monitored by measuring the rate of increase in absorbance at 405 or 415 nm which is proportional to the activity of ALP in the serum.

REAGENT COMPOSITION

R1	
AMP buffer, pH 10.4	434 mmol/l
Magnesium acetate	2.48 mmol/l
Zinc sulfate	1.24 mmol/l
HEDTA	2.48 mmol/l
R2	
p-nitrophenyl phosphate	81.6 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 12 days if refrigerated (2–10 °C) and not contaminated and if chimneys are used.

Absorption of atmospheric CO₂ in open reagent bottle R1 leads to the reduced stability of reagents. To minimize the rate of absorption of CO₂ the use of "chimneys" provided with the kit is recommended. Use of chimneys has shown improvement in calibration stability, provided the storage and use of reagents is as per instructions. Chimneys provided with the kit are to be placed directly into the reagent bottle immediately after opening the cap.

Chimneys can be re-used again for the bottle of the same kit.

Before using the chimney in a new bottle, it is recommended to wash chimneys well with distilled water and dry them thoroughly.

SPECIMEN COLLECTION AND HANDLING

Use serum, plasma (heparin).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability

in serum / plasma:	4 hours	at 20–25 °C
	3 days	at 4–8 °C
	2 months	at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL, Cat. No. XSYS0034 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to original method IFCC.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

Results are calculated automatically by the instrument.

UNIT CONVERSION

U/l x 0.017 = µkat/l

EXPECTED VALUES ⁴

At 37 °C

Females:	4–15 years:	54–369 U/l
	20–50 years:	42–98 U/l
	≥ 60 years:	53–141 U/l
Males:	1–12 years:	54–369 U/l
	20–50 years:	53–128 U/l
	≥ 60 years:	56–119 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 4.5 U/l
Linearity: 1300 U/l
Measuring range: 4.5–1300 U/l

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	84.2	0.58	0.69
Sample 2	217.8	2.34	1.07

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	84.6	1.98	2.34
Sample 2	205.2	2.80	1.37

COMPARISON

A comparison between XL-Systems ALP (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.965 x - 1.68 U/l

r = 0.998

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Reagent R2 also contains <0.7% KOH and is classified as irritant.

Warning

Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302+P352 IF ON SKIN: Wash with plenty of water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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Dátum revízie: 14. 6. 2022

ASSAY PARAMETERS (conventional units)

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details						
Test	ALP	ALP	ALP	ALP	ALP	ALP
Test Code	4	4	4	4	4	4
Report Name	Alkaline Phosphatase	Alkaline Phosphatase	Alkaline Phosphatase	Alkaline Phosphatase	Alkaline Phosphatase	Alkaline Phosphatase
Unit	U/l	U/l	U/l	U/l	U/l	U/l
Decimal Places	0	0	0	0	0	0
Wavelength-Primary	405	405	415	415	405	405
Wavelength-Secondary	700	700	700	700	700	700
Assay type	Rate A	Rate A	Rate A	Rate A	Rate A	Rate A
Curve type	Linear	Linear	Linear	Linear	Linear	Linear
M1 Start	0	0	0	0	0	0
M1 End	0	0	0	0	0	0
M2 Start	21	21	21	32	16	21
M2 End	31	31	39	52	26	25
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	2.4	2.4	2.4	2.4	2.4	2.4
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	4.5	4.5	4.5	4.5	4.5	4.5
Technical Maximum	1300	1300	1300	1300	1300	1300
Y=aX+b						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	0	0	0	0	0	0
Reagent Abs Max	1.4	1.4	1.4	1.4	1.4	1.4
Auto Rerun	No	No	No	No	No	No
Total Reagents	2	2	2	2	2	2
Reagent R1	ALP R1	ALP R1	ALP R1	ALP R1	ALP R1	ALP R1
Reagent R2	ALP R2	ALP R2	ALP R2	ALP R2	ALP R2	ALP R2
Reagent R3	NA	NA	NA	NA	NA	NA
Test Volumes						
Test	ALP	ALP	ALP	ALP	ALP	ALP
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Sample Volumes						
Normal	4	4	4	4	3	4
Dilution Ratio	1	1	1	1	1	1
Increase	8	8	8	8	6	8
Dilution Ratio	1	1	1	1	1	1
Decrease	2	2	2	2	3	2
Dilution Ratio	5	5	5	5	5	5
Standard volume	4	4	4	4	3	4
Reagent Volumes and Stirrer speed						
RGT-1 Volume	160	160	160	160	120	160
R1 Stirrer Speed	Medium	Medium	NA	Medium	Medium	Medium
RGT-2 Volume	40	40	40	40	30	40
R2 Stirrer Speed	High	High	High	High	High	High
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
Reference Ranges						
Test	ALP	ALP	ALP	ALP	ALP	ALP
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
Category Male						
Normal-Lower Limit	53	53	53	53	53	53
Normal-Upper Limit	128	128	128	128	128	128
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Category Female						
Normal-Lower Limit	42	42	42	42	42	42
Normal-Upper Limit	98	98	98	98	98	98
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Revision Number						
Revision	<A-100-ALP-2 01.03.2016>	<A-200-ALP-2 01.03.2016>	<A-300/600-ALP-2 01.03.2016>	<A-640-ALP-2 01.03.2016>	<A-1000-ALP-2 01.03.2016>	<A-180-ALP-2 01.03.2016>

ASSAY PARAMETERS (SI units)

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details						
Test	ALP	ALP	ALP	ALP	ALP	ALP
Test Code	4	4	4	4	4	4
Report Name	Alkaline Phosphatase	Alkaline Phosphatase	Alkaline Phosphatase	Alkaline Phosphatase	Alkaline Phosphatase	Alkaline Phosphatase
Unit	µkat/l	µkat/l	µkat/l	µkat/l	µkat/l	µkat/l
Decimal Places	2	2	2	2	2	2
Wavelength-Primary	405	405	415	415	405	405
Wavelength-Secondary	700	700	700	700	700	700
Assay type	Rate A	Rate A	Rate A	Rate A	Rate A	Rate A
Curve type	Linear	Linear	Linear	Linear	Linear	Linear
M1 Start	0	0	0	0	0	0
M1 End	0	0	0	0	0	0
M2 Start	21	21	21	32	16	21
M2 End	31	31	39	52	26	25
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	2.4	2.4	2.4	2.4	2.4	2.4
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	0.075	0.075	0.075	0.075	0.075	0.075
Technical Maximum	22	22	22	22	22	22
Y=aX+b						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	0	0	0	0	0	0
Reagent Abs Max	1.4	1.4	1.4	1.4	1.4	1.4
Auto Rerun	No	No	No	No	No	No
Total Reagents	2	2	2	2	2	2
Reagent R1	ALP R1	ALP R1	ALP R1	ALP R1	ALP R1	ALP R1
Reagent R2	ALP R2	ALP R2	ALP R2	ALP R2	ALP R2	ALP R2
Reagent R3	NA	NA	NA	NA	NA	NA
Test Volumes						
Test	ALP	ALP	ALP	ALP	ALP	ALP
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Sample Volumes						
Normal	4	4	4	4	3	4
Dilution Ratio	1	1	1	1	1	1
Increase	8	8	8	8	6	8
Dilution Ratio	1	1	1	1	1	1
Decrease	2	2	2	2	3	2
Dilution Ratio	5	5	5	5	5	5
Standard volume	4	4	4	4	3	4
Reagent Volumes and Stirrer speed						
RGT-1 Volume	160	160	160	160	120	160
R1 Stirrer Speed	Medium	Medium	NA	Medium	Medium	Medium
RGT-2 Volume	40	40	40	40	30	40
R2 Stirrer Speed	High	High	High	High	High	High
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
Reference Ranges						
Test	ALP	ALP	ALP	ALP	ALP	ALP
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
Category Male						
Normal-Lower Limit	0.90	0.90	0.90	0.90	0.90	0.90
Normal-Upper Limit	2.18	2.18	2.18	2.18	2.18	2.18
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Category Female						
Normal-Lower Limit	0.71	0.71	0.71	0.71	0.71	0.71
Normal-Upper Limit	1.67	1.67	1.67	1.67	1.67	1.67
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
Revision Number						
Revision	<ASI-100-ALP-2 01.03.2016>	<ASI-200-ALP-2 01.03.2016>	<ASI-300/600-ALP-2 01.03.2016>	<ASI-640-ALP-2 01.03.2016>	<ASI-1000-ALP-2 01.03.2016>	<ASI-180-ALP-2 01.03.2016>

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1. Zilva JF, Pannall PR, "Plasma Enzymes in Diagnosis" in Clinical Chemistry in Diagnosis and Treatment. Lloyd London 1979: Chapter 15 : 343.
2. IFCC method for the measurement of ALP J. Clin. Chem. Clin. Biochem. 1983: 21: 731-48.
3. Young DS. Effects of Drugs on Clinical Laboratory Tests. Third Edition 1990 : 3 : 19-25.
4. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.
5. Kaplan and Pesce (Eds.) Clinical Chemistry, Theory analysis and correlation. Second Edition. CV Mosby Co. 1989.

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