AST/GOT

| Cat. No. | Pack Name | Packaging (Content) |
|----------|---------------------|------------------------------|
| XSYS0016 | AST/GOT 330 | R1: 6 × 44 ml, R2: 6 × 11 ml |
| XSYS0073 | AST/GOT 564 XL-1000 | R1: 6 × 72 ml, R2: 6 × 22 ml |

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of AST/GOT (Aspartate Aminotransferase) in human serum and plasma.

CLINICAL SIGNIFICANCE

AST/GOT is widely distributed with high concentrations in the heart, liver, skeletal muscle, kidney and erythrocytes. Damage or disease to any of these tissues such as myocardial infarction, viral hepatitis, liver necrosis, cirrhosis and muscular dystrophy may result in raised levels of AST/GOT.

PRINCIPLE

This reagent is based on IFCC recommendations, without pyridoxal phosphate. The series of reactions involved in the assay system is as follows:

1. AST/GOT present in the sample catalyses the transfer of the amino group from L-aspartate to 2-oxoglutarate forming oxaloacetate and L-glutamate.
2. Oxaloacetate in the presence of NADH and Malate dehydrogenase (MDH) is reduced to L-malate. In this reaction NADH is oxidized to NAD⁺. The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm due to the oxidation of NADH to NAD⁺.
3. Addition of Lactate dehydrogenase (LDH) to the reagent is necessary to achieve rapid and complete reduction of endogenous pyruvate so that it does not interfere with the assay.

REAGENT COMPOSITION

R1

| | |
|----------------------|------------|
| Tris buffer (pH 7.8) | 110 mmol/l |
| L-aspartic acid | 340 mmol/l |
| LDH | ≥ 4000 U/l |
| MDH | ≥ 750 U/l |

R2

| | |
|----------------|-------------|
| CAPSO | 20 mmol/l |
| 2-oxoglutarate | 85 mmol/l |
| NADH | 1.05 mmol/l |

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum or plasma (heparin, EDTA).

It is recommended to follow NCCLS procedures (or similar standardized conditions).

Loss activity:

| | | |
|---------------|-------------|--------|
| within 3 days | at 2–8 °C | < 8 % |
| within 3 days | at 15–25 °C | < 10 % |

Stability

at least 3 months at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL 4x3 (Cat.No. XSYS0034) or XL MULTICAL 10x3 (Cat.No. XSYS0122) is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to IFCC formulation (without PDP).

QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5, Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

CALCULATION

The XL Results are calculated automatically by the instrument.

UNIT CONVERSION

U/l x 0.017 = µkat/l

EXPECTED VALUES ⁴

At 37 °C

Men up to 35 U/l

Women up to 31 U/l

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

| | |
|---------------------------------|--------------|
| Limit of quantification: | 3.84 U/l |
| Linearity: | 390 U/l |
| Measuring range: | 3.84–390 U/l |

| Intra-assay precision Within run (n=20) | Mean (U/l) | SD (U/l) | CV (%) |
|--|---------------|-------------|-----------|
| Sample 1 | 103.2 | 0.60 | 0.54 |
| Sample 2 | 313.2 | 1.68 | 0.54 |

| Inter-assay precision Run to run (n=20) | Mean (U/l) | SD (U/l) | CV (%) |
|--|---------------|-------------|-----------|
| Sample 1 | 43.8 | 0.60 | 1.37 |
| Sample 2 | 115.2 | 1.08 | 0.92 |

COMPARISON

A comparison between XL-Systems AST/GOT (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.941 x - 3.96 U/l

r = 0.996

INTERFERENCES

Following substances do not interfere:

bilirubin up to 30 mg/dl, triglycerides up to 2000 mg/dl, haemolysis interferes due to AST activity from erythrocytes.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent 1 of the kit contains Aspartic acid, Tris(hydroxymethyl) aminomethane, Sodium hydroxide.

UFI: YMDF-7EN4-1H7J-XF4E

Warning

Hazard statement:

H315 Causes skin irritation.

H319 Causes serious eye irritation.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P302 + P352 IF ON SKIN: Wash with plenty of water and soap.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337 + P313 If eye irritation persists: Get medical advice/attention.

Reagent 2 is not classified as hazardous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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REFERENCES / ЛИТЕРАТУРА / ЛІТЕРАТУРА / LITERATURA / LITERATÚRA / REFERENCIAS

1. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 55-65.
2. Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.
3. Schumann G, Bonora R, Ceriotti F, Férard G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002;40:725-33.
4. Tietz Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.

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