# UREA

Cat. No.	Pack Name	Packaging (Content)
XSYS0020	UREA 275	R1: 5 x 44 ml, R2: 5 x 11 ml
XSYS0075	UREA 564 XL-1000	R1: 6 x 72 ml, R2: 6 x 22 ml

## INTENDED USE


Diagnostic reagent for quantitative *in vitro* determination of Urea in human serum, plasma and urine.

## CLINICAL SIGNIFICANCE

Urea is the major end product of protein nitrogen metabolism in humans. It constitutes the largest fraction of the nonprotein nitrogen component of the blood. Urea is produced in the liver and excreted through the kidneys in the urine. Consequently, the circulating levels of urea depend upon protein intake, protein catabolism and kidney function. Elevated urea levels can occur with dietary changes, diseases which impair kidney function, liver diseases, congestive heart failure, diabetes and infections.

## PRINCIPLE

The enzyme methodology employed in this reagent is based on the reaction first described by Talke and Schubert. To shorten and simplify the assay, the calculations are based on the discovery of Tiffany et al. that urea concentration is proportional to absorbance change over a fixed time interval.

Urea is hydrolysed in the presence of water and Urease to produce ammonia and carbon dioxide.

In the presence of Glutamate Dehydrogenase (GLDH) and reduced Nicotinamide Adenine Dinucleotide (NADH), ammonia combines with  $\alpha$ -ketoglutarate ( $\alpha$ -KG) to produce L-Glutamate.

The reaction is monitored by measuring the rate of decrease in absorbance at 340 nm as NADH is converted to  $\text{NAD}^+$ . The rate of oxidation of NADH is proportional to the urea concentration in the sample.

## REAGENT COMPOSITION

### R1

Tris Buffer	100 mmol/l
$\alpha$ -Ketoglutarate	5.49 mmol/l
Urease (Jack Bean)	$\geq 10$ KU/l
GLDH (Microorganism)	$\geq 3.8$ KU/l

### R2

NADH	1.66 mmol/l
------	-------------

Also contains non-reactive fillers and stabilisers.

## REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

## STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

## SPECIMEN COLLECTION AND HANDLING

Use unheamolytic serum, plasma or urine.

**Do not** use samples preserved with fluoride and heparin ammonium salt as anti-coagulant.

It is recommended to follow NCCLS procedures (or similar standardized conditions). **Urine:** Dilute urine 1 + 99 with distilled water and multiply results by 100.

## Stability

in serum/plasma:	7 days	at 20–25 °C
	7 days	at 4–8 °C
in urine:	1 year	at -20 °C
	2 days	at 20–25 °C
	2 days	at 4–8 °C
	1 month	at -20 °C

Discard contaminated specimens.

## CALIBRATION

Calibration with calibrator XL MULTICAL 4x3 (Cat. No. XSYS0034) or XL MULTICAL 10x3 (Cat. No. XSYS0122) is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

## Traceability:

This calibrator has been standardized using the SRM 909c reference material.

## QUALITY CONTROL

For quality control ERBA NORM 4x5 (Cat. No. BLT00080) or ERBA NORM 10x5 (Cat. No. XSYS0123) and ERBA PATH 4x5 (Cat. No. BLT00081) or ERBA PATH 10x5 (Cat. No. XSYS0124) are recommended.

## CALCULATION

Results are calculated automatically by the instrument.

## UNIT CONVERSION

mg/dl x 0.1665 = mmol/l  
Urea (mg/dl) x 0.467 = BUN (mg/dl)  
BUN (mg/dl) x 2.14 = Urea (mg/dl)

## EXPECTED VALUES <sup>1,2</sup>

### In Serum / Plasma <sup>1</sup>

Adults	(mg/dl)	(mmol/l)
Global	17–43	2.8–7.2
Women < 50 years	15–40	2.6–6.7
Women > 50 years	21–43	3.5–7.2
Men < 50 years	19–44	3.2–7.3
Men > 50 years	18–55	3.0–9.2

### Children

1–3 years	11–36	1.8–6.0
4–13 years	15–36	2.5–6.0
14–19 years	18–45	2.9–7.5

### Urea / Creatinine ratio <sup>1</sup>

20–35 [(mg/dl)/(mg/dl)]

### Urea in Urine <sup>2</sup>

26–43 g/24 h (0.43–0.72 mol/24 h)

**It is recommended that each laboratory verify this range or derives reference interval for the population it serves.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

**Limit of quantification:** 1.5 mg/dl  
**Linearity:** 300 mg/dl of Urea or 140 mg/dl of BUN

**Measuring range:** 1.5–300 mg/dl

Intra-assay precision Within run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	28.08	0.287	1.02
Sample 2	27.49	0.240	0.94

Inter-assay precision Run to run (n=20)	Mean (mg/dl)	SD (mg/dl)	CV (%)
Sample 1	45.09	0.719	1.61
Sample 2	203.59	2.395	1.58

## COMPARISON

A comparison between XL-Systems Urea (y) and a commercially available test (x) using 40 samples gave following results:

y = 1.034 x - 0.295 mg/dl  
r = 0.994

## INTERFERENCES

Following substances do not interfere:

haemoglobin up to 7.5 g/l, bilirubin up to 30 mg/dl, triglycerides up to 2000 mg/dl.

## WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents of the kit are not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

Erba Lachema s.r.o., Karásek 2219/1d, 621 00 Brno, CZ  
e-mail: diagnostics@erba.com, www.erbamannheim.com

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- Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 374-7.
- Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 1838.
- Talka, H. Schubert, G.E. Klin. Wochschr. 19; 43: 174.
- Tiffany, T.O. Jansen, J.M.Burtis CA, Overton JB, Scott CD. Clin Chem 1972; 18: 829-40.
- Kaplan, LA. in "Clinical Theory, Analysis and Correlation." Kaplan LA, Pesce AJ. (Ed) C V Mosby Company St Louis 1984; 1257-61.
- Shephard, MD, Mezzachi, RD. Clin. Biochem. Revs. 1983; 4: 61-7.
- Young, D.S. Effects of Drugs on Clinical Laboratory Tests. Third Edition. 1990; 21: 5.
- Wachtel, M. et al, Creation and Verification of Reference Intervals. Laboratory Medicine 1995; 26: 593-7.
- National Committee for Clinical Laboratory Standards. User evaluation of Precision Performance of Clinical Chemistry Devices. NCCLS 1984, NCCLS, Publication EP5-T.
- Tietz Textbook of Clinical Chemistry. Burtis CA and Ashwood ER, Fifth Edition, 2012.

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