LDL DIRECT

| Cat. No. | Pack Name | Packaging (Content) |
|----------|-----------|------------------------------|
| XSYS0044 | LDL C 80 | R1: 2 x 30 ml, R2: 2 x 10 ml |

INTENDED USE


Diagnostic reagent for quantitative *in vitro* determination of LDL Cholesterol in human serum and plasma.

CLINICAL SIGNIFICANCE

Low Density Lipoproteins (LDL) are synthesized in the liver by the action of various lipolytic enzymes on triglyceride-rich Very Low Density Lipoproteins (VLDLs). Specific LDL receptors exist to facilitate the elimination of LDL from plasma by liver parenchymal cells. It has been shown that most of the cholesterol stored in atherosclerotic plaques originates from LDL. For this reason the LDL Cholesterol concentration is considered to be the most important clinical predictor, of all single parameters, with respect to coronary atherosclerosis.²⁻⁸ Accurate measurement of LDL Cholesterol is of vital importance in therapies which focus on lipid reduction to prevent atherosclerosis or reduce its progress and to avoid plaque rupture. Can be applied on automated analyzers.

PRINCIPLE

The assay is based on a modified polyvinyl sulfonic acid (PVS) and polyethylene-glycol methyl ether (PEGME) coupled classic precipitation method with the improvements in using optimized quantities of PVS/PEGME and selected detergents. LDL, VLDL, and chylomicron (CM) react with PVS and PEGME and the reaction results in inaccessibility of LDL, VLDL and CM by cholesterol oxidase (CHOD) and cholesterol esterase (CHER), whereas HDL reacts with the enzymes. Addition of R2 containing a specific detergent releases LDL from the PVS/PEGME complex. The released LDL reacts with the enzymes to produce H₂O₂ which is quantified by the Trinder reaction.

REAGENT COMPOSITION

| | | |
|--|--|-----------|
| R1 | | |
| MES buffer (pH 6.5) | | 50 mmol/l |
| Polyvinylsulfonic acid | | 50 mg/l |
| Polyethyleneglycolmethylether | | 30 ml/l |
| 4-aminoantipyrine | | 0.9 g/l |
| Cholesterol esterase | | 5 kU/l |
| Cholesterol oxidase | | 20 kU/l |
| Peroxidase | | 5 kU/l |
| R2 | | |
| MES buffer (pH 6.5) | | 50 mmol/l |
| Detergent | | |
| TODB N,N-Bis(4-sulfobutyl)-3-methylaniline | | 3 mmol/l |

REAGENT PREPARATION

Reagents R1 and R2 are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C. On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated. Once opened both reagents R1 & R2 are stable for 60 days at 2–8 °C, when protected from contamination. Reagents are light-sensitive. Do not let bottles remain open. Keep containers tightly closed.

SPECIMEN COLLECTION AND HANDLING

Use serum or heparin plasma. It is recommended to follow NCCLS procedures (or similar standardized conditions).

Stability in serum/plasma: 12 hours at 20–25 °C
10 days at 4–8 °C
12 weeks at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with HDL/LDL CAL, Cat. No. XSYS0061 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized using the NIST SRM 1951b reference material.

QUALITY CONTROL

For quality control ERBA NORM, Cat. No. BLT00080 and ERBA PATH, Cat. No. BLT00081 are recommended.

CALCULATION

Results are calculated automatically by the instrument.

UNIT CONVERSION

mg/dl x 0.026 = mmol/l

EXPECTED VALUES ¹¹

| | |
|---------------------|----------------------|
| Less than 100 mg/dl | – optimal |
| 100–129 mg/dl | – near/above optimal |
| 130–159 mg/dl | – borderline high |
| 160–189 mg/dl | – high |
| ≥ 190 mg/dl | – very high |

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 2.60 mg/dl
Linearity: 263 mg/dl
Measuring range: 2.60–263 mg/dl

PRECISION

| Intra-assay precision Within run (n=20) | Mean (mg/dl) | SD (mg/dl) | CV (%) |
|--|-----------------|---------------|-----------|
| Sample 1 | 50.346 | 0.885 | 1.72 |
| Sample 2 | 82.308 | 1.808 | 2.21 |

| Inter-assay precision Run to run (n=20) | Mean (mg/dl) | SD (mg/dl) | CV (%) |
|--|-----------------|---------------|-----------|
| Sample 1 | 47.00 | 0.885 | 1.91 |
| Sample 2 | 92.69 | 1.500 | 1.61 |

COMPARISON

A comparison between XL-Systems LDL DIRECT (y) and a commercially available test (x) using 40 samples gave following results:

y = 0.964x -1.615 mg/dl
r = 0.995

INTERFERENCES

Following substances do not interfere: haemoglobin up to 10 g/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl. Interference by N-acetylcysteine (NAC), acetaminophen and metamizole causes falsely low results. To carry out the test, blood withdrawal should be performed prior to administration of drugs.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person. Reagents of the kit are not classified as dangerous.

WASTE MANAGEMENT

Please refer to local legal requirements.

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ASSAY PARAMETERS (conventional units)

| Instrument | XL-100 EM-100 | XL-200 EM-200 | XL-300/600 EM-360 | XL-640 | XL-1000 | XL-180 |
|----------------------|------------------|------------------|----------------------|------------|------------|------------|
| Test Details | | | | | | |
| Test | LDL | LDL | LDL | LDL | LDL | LDL |
| Test Code | 50 | 50 | 50 | 50 | 50 | 50 |
| Report Name | LDL Direct | LDL Direct | LDL Direct | LDL Direct | LDL Direct | LDL Direct |
| Unit | mg/dl | mg/dl | mg/dl | mg/dl | mg/dl | mg/dl |
| Decimal Places | 1 | 1 | 1 | 1 | 1 | 1 |
| Wavelength-Primary | 600 | 600 | 600 | 600 | 600 | 600 |
| Wavelength-Secondary | 700 | 700 | 700 | 700 | 700 | 700 |
| Assay type | 2-Point | 2-Point | 2-Point | 2-Point | 2-Point | 2-Point |
| Curve type | Linear | Linear | Linear | Linear | Linear | Linear |
| M1 Start | 16 | 16 | 12 | 24 | 14 | 16 |
| M1 End | 16 | 16 | 12 | 24 | 14 | 16 |
| M2 Start | 34 | 34 | 48 | 61 | 29 | 34 |
| M2 End | 34 | 36 | 51 | 63 | 31 | 34 |
| Sample replicates | 1 | 1 | 1 | 1 | 1 | 1 |
| Standard replicates | 3 | 3 | 3 | 3 | 3 | 3 |
| Control replicates | 1 | 1 | 1 | 1 | 1 | 1 |
| Control interval | 0 | 0 | 0 | 0 | 0 | 0 |
| Reaction Direction | Increasing | Increasing | Increasing | Increasing | Increasing | Increasing |
| React. Abs. Limit | 0.8 | 0.8 | 0.9 | 0.9 | 0.9 | 0.9 |
| Prozone Limit % | 0 | 0 | 0 | 0 | 0 | 0 |
| Prozone Check | Lower | Lower | Lower | Lower | Lower | Lower |
| Linearity Limit % | 0 | 0 | 0 | 0 | 0 | 0 |
| Delta Abs/Min | 0 | 0 | 0 | 0 | 0 | 0 |
| Technical Minimum | 2.60 | 2.60 | 2.60 | 2.60 | 2.60 | 2.60 |
| Technical Maximum | 263 | 263 | 263 | 263 | 263 | 263 |
| Y=aX+b | | | | | | |
| a= | 1 | 1 | 1 | 1 | 1 | 1 |
| b= | 0 | 0 | 0 | 0 | 0 | 0 |
| Reagent Abs Min | 0 | 0 | 0 | 0 | 0 | 0 |
| Reagent Abs Max | 0.3 | 0.3 | 0.3 | 0.3 | 0.3 | 0.3 |
| Auto Rerun | No | No | No | No | No | No |
| Total Reagents | 2 | 2 | 2 | 2 | 2 | 2 |
| Reagent R1 | LDL R1 | LDL R1 | LDL R1 | LDL R1 | LDL R1 | LDL R1 |
| Reagent R2 | LDL R2 | LDL R2 | LDL R2 | LDL R2 | LDL R2 | LDL R2 |
| Reagent R3 | NA | NA | NA | NA | NA | NA |

Test Volumes

| | | | | | | |
|-----------------------|-------|-------|-------|-------|-------|-------|
| Test | LDL | LDL | LDL | LDL | LDL | LDL |
| Sample Type | SERUM | SERUM | SERUM | SERUM | SERUM | SERUM |
| Sample Volumes | | | | | | |
| Normal | 2 | 2 | 3 | 2 | 2 | 2 |
| Dilution Ratio | 1 | 1 | 1 | 1 | 1 | 1 |
| Increase | 4 | 4 | 4 | 4 | 4 | 4 |
| Dilution Ratio | 1 | 1 | 1 | 1 | 1 | 1 |
| Decrease | 2 | 2 | 2 | 2 | 2 | 2 |
| Dilution Ratio | 5 | 5 | 1 | 5 | 5 | 5 |
| Standard volume | 2 | 2 | 3 | 2 | 2 | 2 |

Reagent Volumes and Stirrer speed

| | | | | | | |
|------------------|------|------|-----|------|------|------|
| RGT-1 Volume | 180 | 180 | 210 | 180 | 180 | 180 |
| R1 Stirrer Speed | High | High | NA | High | High | High |
| RGT-2 Volume | 60 | 60 | 70 | 60 | 0 | 60 |
| R2 Stirrer Speed | High | High | NA | High | NA | High |
| RGT-3 Volume | 0 | 0 | 0 | 0 | 60 | 0 |
| R3 Stirrer Speed | NA | NA | NA | NA | High | NA |

Reference Ranges

| | | | | | | |
|------------------------|---------|---------|---------|---------|---------|---------|
| Test | LDL | LDL | LDL | LDL | LDL | LDL |
| Sample Type | SERUM | SERUM | SERUM | SERUM | SERUM | SERUM |
| Reference Range | Default | Default | Default | Default | Default | Default |
| Category Male | | | | | | |
| Normal-Lower Limit | 0 | 0 | 0 | 0 | 0 | 0 |
| Normal-Upper Limit | 130 | 130 | 130 | 130 | 130 | 130 |
| Panic-Lower Limit | NA | NA | NA | NA | NA | NA |
| Panic-Upper Limit | NA | NA | NA | NA | NA | NA |
| Category Female | | | | | | |
| Normal-Lower Limit | 0 | 0 | 0 | 0 | 0 | 0 |
| Normal-Upper Limit | 130 | 130 | 130 | 130 | 130 | 130 |
| Panic-Lower Limit | NA | NA | NA | NA | NA | NA |
| Panic-Upper Limit | NA | NA | NA | NA | NA | NA |

Revision Number

| | | | | | | |
|----------|---------------------------------|---------------------------------|-------------------------------------|---------------------------------|----------------------------------|---------------------------------|
| Revision | <A-100- LDL-2 26.09.2013> | <A-200- LDL-2 26.09.2013> | <A-300/600- LDL-2 26.09.2013> | <A-640- LDL-2 26.09.2013> | <A-1000- LDL-2 26.09.2013> | <A-180- LDL-1 12.12.2013> |
|----------|---------------------------------|---------------------------------|-------------------------------------|---------------------------------|----------------------------------|---------------------------------|

ASSAY PARAMETERS (SI units)

| Instrument | XL-100 EM-100 | XL-200 EM-200 | XL-300/600 EM-360 | XL-640 | XL-1000 | XL-180 |
|----------------------|------------------|------------------|----------------------|------------|------------|------------|
| Test Details | | | | | | |
| Test | LDL | LDL | LDL | LDL | LDL | LDL |
| Test Code | 50 | 50 | 50 | 50 | 50 | 50 |
| Report Name | LDL Direct | LDL Direct | LDL Direct | LDL Direct | LDL Direct | LDL Direct |
| Unit | mmol/l | mmol/l | mmol/l | mmol/l | mmol/l | mmol/l |
| Decimal Places | 2 | 2 | 2 | 2 | 2 | 2 |
| Wavelength-Primary | 600 | 600 | 600 | 600 | 600 | 600 |
| Wavelength-Secondary | 700 | 700 | 700 | 700 | 700 | 700 |
| Assay type | 2-Point | 2-Point | 2-Point | 2-Point | 2-Point | 2-Point |
| Curve type | Linear | Linear | Linear | Linear | Linear | Linear |
| M1 Start | 16 | 16 | 12 | 24 | 14 | 16 |
| M1 End | 16 | 16 | 12 | 24 | 14 | 16 |
| M2 Start | 34 | 34 | 48 | 61 | 29 | 34 |
| M2 End | 34 | 36 | 51 | 63 | 31 | 34 |
| Sample replicates | 1 | 1 | 1 | 1 | 1 | 1 |
| Standard replicates | 3 | 3 | 3 | 3 | 3 | 3 |
| Control replicates | 1 | 1 | 1 | 1 | 1 | 1 |
| Control interval | 0 | 0 | 0 | 0 | 0 | 0 |
| Reaction Direction | Increasing | Increasing | Increasing | Increasing | Increasing | Increasing |
| React. Abs. Limit | 0.8 | 0.8 | 0.9 | 0.9 | 0.9 | 0.9 |
| Prozone Limit % | 0 | 0 | 0 | 0 | 0 | 0 |
| Prozone Check | Lower | Lower | Lower | Lower | Lower | Lower |
| Linearity Limit % | 0 | 0 | 0 | 0 | 0 | 0 |
| Delta Abs/Min | 0 | 0 | 0 | 0 | 0 | 0 |
| Technical Minimum | 0.07 | 0.07 | 0.07 | 0.07 | 0.07 | 0.07 |
| Technical Maximum | 6.84 | 6.84 | 6.84 | 6.84 | 6.84 | 6.84 |
| Y=aX+b | | | | | | |
| a= | 1 | 1 | 1 | 1 | 1 | 1 |
| b= | 0 | 0 | 0 | 0 | 0 | 0 |
| Reagent Abs Min | 0 | 0 | 0 | 0 | 0 | 0 |
| Reagent Abs Max | 0.3 | 0.3 | 0.3 | 0.3 | 0.3 | 0.3 |
| Auto Rerun | No | No | No | No | No | No |
| Total Reagents | 2 | 2 | 2 | 2 | 2 | 2 |
| Reagent R1 | LDL R1 | LDL R1 | LDL R1 | LDL R1 | LDL R1 | LDL R1 |
| Reagent R2 | LDL R2 | LDL R2 | LDL R2 | LDL R2 | LDL R2 | LDL R2 |
| Reagent R3 | NA | NA | NA | NA | NA | NA |

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1. "Executive Summary of the Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (Adult treatment Panel III)", JAMA, 285:2486 (2001).
2. Crouse JR et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26; 566 (1985).
3. Castelli WP et al., Cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
4. Barr DP, Russ EM, Eder HA, Protein-lipid relationships in human plasma, Am. J. Med., 11;480 (1951).
5. Badimon JJ, Badimon L., Fuester V., Regression of Atherosclerotic Lesions by High-Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit, Journal of Clinical Investigation, 85:1234-41(1990).
6. Gordon T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am.J. Med., 62;707 (1977).
6. Kannel WB, Castelli WP, Gordon T., Cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Am. Intern. Med., 90; 85 (1979).
7. William P., Robinson D., Baily A., High density lipoprotein and coronary risk factor, Lancet, 1;72 (1979).
8. Castelli, W. P., et al, Cholesterol and other lipids in coronary heart disease.
9. Third Report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of the High Blood Cholesterol in Adults (Adult Treatment Panel III).
10. Pisani T, Gebiski CP, Leary Et, et al. Accurate Direct Determination of Low-Density Lipoprotein Cholesterol Assay. Arch Pathol Lab Med 1995; 119:1127)
11. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Burtis, C.A., Ashwood, E.R., Bruns, D.E.; 5th edition, WB Saunders Comp., 2012.

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