# RHEUMATOID FACTOR

Cat. No.	Pack Name	Packaging (Content)
XSYS0048	RF	R1: 2 x 40 ml (Buffer), R2: 2 x 8 ml (RF reagent)

EN

## INTENDED USE

Quantitative determination of Rheumatoid Factor (RF) in human serum by turbidimetric immunoassay (Aggregated human IgG method).

## CLINICAL SIGNIFICANCE

The diagnosis of rheumatoid arthritis (RA) is based largely on clinical examination, but laboratory tests (e.g. RF Test) are useful to support the clinical diagnosis and to evaluate the severity and course of the disease in the individual patient. RF is a term used to describe a variety of antibodies (in most cases of the IgM type) that will react with modified human IgG (e.g. IgG in circulating immune complexes, IgG adsorbed to latex, etc.) and IgG of animal origin. RF is highly associated with rheumatoid arthritis, as high as 90 % of patients with RA have RF titers of more than 20 IU/ml.

## METHODOLOGY

Measurement of antigen-antibody reaction by the end-point method.

## REAGENT COMPOSITION

### R1 (Buffer)

Good's buffer (pH 7.4) 50 mmol/l

Sodium azide (< 0.1 %)

### R2 (RF reagent)

Heat-aggregated human IgG < 0.5 mg/ml

Sodium azide (< 0.1 %)

## REAGENT PREPARATION

Liquid reagents are ready to use.

## STABILITY AND STORAGE

The reagents are stable until expiry date when kept at 2–8 °C. Stability in the instrument is at least 4 weeks if contamination is avoided. Do not freeze.

## SAMPLE COLLECTION

Use fresh serum. If the test can not be carried out on the same day, the serum may be stored at 2–8 °C for 48 hours. If stored for a longer period, the sample should be frozen.

## MATERIALS REQUIRED BUT NOT PROVIDED

- Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance accurately at 340 nm may be used.
- Analyser specific consumables such as sample cups.
- Controls.
- Saline (9 g/l NaCl)

## ASSAY PROCEDURE

Refer to the assay parameters for details.

## CALIBRATION

Blank: Saline

Cat. No.:	Product name	Pack name	Content
XSYS0052	RF CALIBRATOR	RF CAL SH	1 x 1 ml

Calibration curve: generate a 6 point calibration curve by diluting the calibrator 1:32, 1:16, 1:8, 1:4, 1:2 and undiluted in saline.

## Calibration frequency

Calibration verification: Not necessary.

## Quality control

For quality control use

Cat. No.:	Product name	Pack name	Content
BLT20039	RF CONTROL	RF CON	1 x 1 ml
BLT20034	MULTICONTROL LEVEL 1	MULTICON L1	1 x 1 ml
BLT20035	MULTICONTROL LEVEL 2	MULTICON L2	1 x 1 ml

## CALCULATION

Results are calculated automatically by the instrument.

## EXPECTED VALUES

### Reference Values

Normal values: 0–20 IU/ml (WHO)

**This range is given for orientation only. Each Laboratory should establish its own reference values.**

## PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

### Accuracy

A comparison between XL-Systems RF (y) and a commercially available test (x) using 50 samples gave following results:

$y = 1.189x - 3.690$  IU/ml

$r = 0.996$

Measuring Range: 0–500 IU/ml

Detection Limit: 1.61 IU/ml

Hook Effect: No risk

## PRECISION

Intra-assay precision Within run	Mean (IU/ml)	SD (IU/ml)	CV (%)
Sample 1	54.5	0.85	1.56
Sample 2	146.5	3.99	2.72

Inter-assay precision Run to run	Mean (IU/ml)	SD (IU/ml)	CV (%)
Sample 1	24.95	1.57	6.56
Sample 2	50.08	2.36	4.71

**Specificity:** Monospecific

**Interferences:** No interference for: Hemoglobin (250 mg/dl), Bilirubin (40 mg/dl), Triacylglycerols (1 g/dl).

**Limitations:** None

**Stability at 4 °C:** At least 3 years after production

## WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Components of the kit are not classified as dangerous but contain less than 0.1 % sodium azide - classified as very toxic and dangerous substance for the environment.

2. Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion. Flush drains with water thoroughly after disposing of fluids containing sodium azide.

3. Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.

## Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagents R1 and R2 are not classified as dangerous.

## WASTE MANAGEMENT

Please refer to local legal requirements.

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# ASSAY PARAMETERS

Instrument	XL-100 EM-100	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
<b>Test Details</b>						
Test	RF	RF	RF	RF	RF	RF
Test Code	38	38	38	38	38	38
Report Name	Rheumatoid Factor	Rheumatoid Factor	Rheumatoid Factor	Rheumatoid Factor	Rheumatoid Factor	Rheumatoid Factor
Unit	IU/ml	IU/ml	IU/ml	IU/ml	IU/ml	IU/ml
Decimal Places	1	1	1	1	1	1
Wavelength-Primary	340	340	340	340	340	340
Wavelength-Secondary	0	0	0	0	0	0
Assay type	2-Point	2-Point	2-Point	2-Point	2-Point	2-Point
Curve type	Cubic Spline	Cubic Spline	Cubic Spline	Cubic Spline	Cubic Spline	Cubic Spline
M1 Start	16	16	12	24	10	16
M1 End	16	16	12	24	10	16
M2 Start	34	36	50	62	31	34
M2 End	34	36	50	62	31	34
Sample replicates	1	1	1	1	1	1
Standard replicates	3	3	3	3	3	3
Control replicates	1	1	1	1	1	1
Control interval	0	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	NA	NA	NA	NA	NA	NA
Prozone Limit %	0	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0	0
Technical Minimum	NA	NA	NA	NA	NA	NA
Technical Maximum	NA	NA	NA	NA	NA	NA
<b>Y=aX+b</b>						
a=	1	1	1	1	1	1
b=	0	0	0	0	0	0
Reagent Abs Min	NA	NA	NA	NA	NA	NA
Reagent Abs Max	0	0	0	0	0	0
Auto Rerun	No	No	No	No	No	No
Total Reagents	2	2	2	2	2	2
Reagent R1	RF R1	RF R1	RF R1	RF R1	RF R1	RF R1
Reagent R2	RF R2	RF R2	RF R2	RF R2	RF R2	RF R2
Reagent R3	NA	NA	NA	NA	NA	NA
<b>Test Volumes</b>						
Test	RF	RF	RF	RF	RF	RF
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
<b>Sample Volumes</b>						
Normal	10	10	12	10	10	10
Dilution Ratio	1	1	1	1	1	1
Increase	20	20	24	20	20	20
Dilution Ratio	1	1	1	1	1	1
Decrease	5	5	6	5	5	5
Dilution Ratio	1	1	1	1	1	1
Standard volume	10	10	12	10	10	10
<b>Reagent Volumes and Stirrer speed</b>						
RGT-1 Volume	180	180	200	180	180	180
R1 Stirrer Speed	High	High	NA	High	High	High
RGT-2 Volume	36	36	40	36	36	36
R2 Stirrer Speed	High	High	NA	High	High	High
RGT-3 Volume	0	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA	NA
<b>Reference Ranges</b>						
Test	RF	RF	RF	RF	RF	RF
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default	Default
<b>Category Male</b>						
Normal-Lower Limit	0	0	0	0	0	0
Normal-Upper Limit	20	20	20	20	20	20
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
<b>Category Female</b>						
Normal-Lower Limit	0	0	0	0	0	0
Normal-Upper Limit	20	20	20	20	20	20
Panic-Lower Limit	NA	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA	NA
<b>Revision Number</b>						
Revision	<A-100- RF-1 20.08.2013>	<A-200- RF-1 20.08.2013>	<A-300/600- RF-1 20.08.2013>	<A-640- RF-1 20.08.2013>	<A-1000- RF-1 20.08.2013>	<A-180- RF-1 12.12.2013>

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