HbA1c

| Cat. No. | Pack Name | Packaging (Content) |
|----------|-----------|--|
| XSYS0054 | HBA1C | R1: 1 x 24 ml, R2a: 1 x 8 ml, R2b: 1 x 4 ml, R3: 2 x 50 ml |

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of hemoglobin A1c in whole blood on photometric systems.

CLINICAL SIGNIFICANCE

Hemoglobin A1c (HbA1c) is a glycosylated hemoglobin which is formed by the non-enzymatic reaction of glucose with native hemoglobin. This process runs continuously throughout the circulatory life of the red cell (average life time 100–120 days). The rate of glycation is directly proportional to the concentration of glucose in the blood. The blood level of HbA1c represents the average blood glucose level over the preceding 6 to 8 weeks (due to the kinetics of erythrocyte turnover this period is more affected by the blood glucose level than the preceding weeks). Therefore, HbA1c is suitable for retrospective long-term monitoring of blood glucose concentration in individuals with diabetes mellitus. Clinical studies have shown that lowering of HbA1c level can help to prevent or delay the incidence of late diabetic complications.

As the amount of HbA1c also depends on the total quantity of hemoglobin the reported HbA1c value is indicated as a percentage of the total hemoglobin concentration.

Falsely low values (low HbA1c despite high blood glucose) may occur in people with conditions with shortened red blood cell survival (hemolytic diseases) or significant recent blood loss (higher fraction of young erythrocytes). Falsely high values (high HbA1c despite normal blood glucose) have been reported in iron deficiency anemia (high proportion of old erythrocytes). These circumstances have to be considered in clinical interpretation of HbA1c values.

METHODOLOGY

Particle enhanced immunoturbidimetric test.

HbA1c is determined directly without measurement of total hemoglobin.

PRINCIPLE

Total Hb and HbA1c in hemolyzed blood bind with the same affinity to particles in R1. The amount of binding is proportional to the relative concentration of both substances in the blood. Mouse anti-human HbA1c monoclonal antibody (R2a) binds to particle bound HbA1c. Goat anti-mouse IgG polyclonal antibody (R2b) interacts with the monoclonal mouse anti-human HbA1c antibody and agglutination takes place. The measured absorbance is proportional to the HbA1c bound to particles, which in turn is proportional to the percentage of HbA1c in the sample.

STANDARDIZATION

The assay is standardized according to the approved IFCC reference method [3]. Calibration according to DCCT/NGSP is also possible. Corresponding calibrator values are listed in the package insert of the calibrator set HbA1c liquid.

DCCT/NGSP and IFCC values show a linear relationship and can therefore be calculated from each other. Also new IFCC recommended units mmol/mol (mmol HbA1c / mol Hb) can be easily calculated from IFCC values. Recalculation equations are following:

IFCC = (NGSP – 2.15) / 0.915

NGSP = 0.915 x IFCC + 2.15

mmol/mol = 10 x IFCC

IFCC: International Federation of Clinical Chemistry [3,4]

DCCT: Diabetes Control and Complications Trial [5]

NGSP: National Glycohemoglobin Standardization Program [6]

REAGENT COMPOSITION

| | | |
|-------------|--|-----------|
| R1: | Buffer | 20 mmol/l |
| | Latex | 1.5 % |
| R2a: | Buffer | 10 mmol/l |
| | Mouse anti-human HbA1c monoclonal antibody | 5.5 mg/dl |
| R2b: | Buffer | 1 mmol/l |
| | Goat anti-mouse IgG polyclonal antibody | 67 mg/dl |
| | Stabilizers | |
| R3: | Hemolyzing solution | |

REAGENT PREPARATION

Transfer 4 ml of R2b into bottle R2a and mix well immediately.

Ratio between R2a and R2b must be 2/1. Stability of premixed R2a/R2b: One month stored at 2–8 °C.

STABILITY AND STORAGE

The reagents are stable until expiry date when kept at 2–8 °C. Stability in the instrument is at least 4 weeks if contamination is avoided. Do not freeze.

SPECIMEN COLLECTION

Whole blood collected with EDTA.

Sample preparation:

Hemolyzing Solution (R3) 500 µl

Sample/calibrator/control 10 µl

Mix and allow to stand for 5 minutes or until complete lysis is apparent.

In the case of instruments which can process 3 reagents, sample preparation can be performed on board.

Specimen stability:

Whole blood 1 week at 2–8 °C

Hemolysate 10 hours at 15–25 °C

Hemolysate 10 days at 2–8 °C

ASSAY PROCEDURE

Application sheets for automated systems are available on request.

Wavelength 660 nm

Optical path 1 cm

Temperature 37 °C

Measurement Against air

3-component system – ready-to-use

| | |
|---|--------|
| Sample or calibrator | 20 µl |
| Reagent 1 | 750 µl |
| Mix, incubate for 2 min., then add | |
| Reagent 2a | 250 µl |
| Mix, incubate for 3 min., then add: | |
| Reagent 2b | 125 µl |
| Mix, read absorbance after exactly 2 min. | |

2-component system - premixed R2a/R2b

| | |
|---|---------|
| Sample or calibrator | 30 µl |
| Reagent 1 | 1000 µl |
| Mix, incubate for 5 min., then add: | |
| Reagent 2a/2b | 500 µl |
| Mix, read absorbance after exactly 5 min. | |

MATERIALS REQUIRED BUT NOT PROVIDED

General laboratory equipment.

CALIBRATION

The concentration of HbA1c in unknown samples is derived from a calibration curve using an appropriate mathematical model such as spline. The calibration curve is obtained with 4 calibrators at different levels and NaCl solution (9 g/l) for determination of the zero value.

Stability of calibration: 3-component system 8 weeks
2-component system 2 weeks

For calibration use the HBA1C CAL SET with 4 different concentrations (XSYS0057).

CONTROLS

For internal quality control use the HBA1C CON L (XSYS0055), HBA1C CON H (XSYS0056).

CALCULATION

Results are calculated automatically by the instrument.

EXPECTED VALUES ⁷

Reference intervals should be established or verified by the laboratory based on an appropriate non-diabetic patient population.

| | % NGSP | % IFCC | mmol/mol |
|-------------------|--------|--------|----------|
| Non-diabetics | 4–6 | 3–4 | 30–40 |
| Target of therapy | < 7 | < 5 | < 50 |
| Change of therapy | > 8 | > 6 | > 60 |

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Accuracy

A comparison between XL-Systems HbA1c (y) and commercially available test (x) using 70 samples were carried. Serum samples were assayed in parallel and the results compared by linear regression analysis according to Passing-Bablok. The following statistics were obtained:
r = 0.973
y = 0.935 x + 0.507 % NGSP

Measuring Range

The test has been developed to determine concentrations of HbA1c within a measuring range from 4.90–15.87 % DCCT/NGSP, 3–15 % IFCC, 30–150 mmol/mol.

The assay is applicable for hemoglobin concentrations in blood from 6.6 to 26 g/dl.

Specificity / Interferences

Due to its antibodies, HbA1c is a specific immunoassay for human HbA1c. No interference was observed by ascorbic acid up to 60 mg/dl, conjugated and unconjugated bilirubin up to 60 mg/dl, lipemia up to 2000 mg/dl triglycerides, RF up to 250 IU/ml, carbamylated Hb up to 10 mmol/l, acetylated Hb up to 10 mmol/l and urea up to 300 mg/dl.

No interference is observed by Hemoglobin variants HbS and elevated HbA2. Elevated levels of HbF may lead to false HbA1c values. Alcoholism and ingestion of large doses of aspirin may lead to inconsistent results [1].

Sensitivity /Limit of Detection

The limit of detection is 1.84 % NGSP HbA1c.

PRECISION

Values according to NGSP

| Within-run precision n = 20 | Mean (% NGSP) | SD (% NGSP) | CV (%) |
|--------------------------------|------------------|----------------|-----------|
| Sample 1 | 6.98 | 0.04 | 0.60 |
| Sample 2 | 10.31 | 0.06 | 0.58 |

| Between day precision n = 20 | Mean (% NGSP) | SD (% NGSP) | CV (%) |
|---------------------------------|------------------|----------------|-----------|
| Sample 1 | 6.74 | 0.05 | 0.80 |
| Sample 2 | 9.98 | 0.19 | 1.87 |

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Components of the kit are not classified as dangerous but reagents R1, R2a and R2b contain less than 0.1 % sodium azide - classified as very toxic and dangerous substance for the environment.

In very rare cases, samples of patients with gammopathy might give falsified results.

Immediately after HbA1c measurement cleaning of cuvettes is necessary. Use the alkaline cuvette washing solution which is recommended by the analyzer manufacturer.

Take necessary precautions for use of laboratory reagent.

WASTE MANAGEMENT

Please refer to local legal requirements.

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ASSAY PARAMETERS (conventional units)

| Instrument | XL-100 EM-100 | XL-200 EM-200 | XL-300/600 EM-360 | XL-640 | XL-1000 | XL-180 |
|----------------------|------------------|------------------|----------------------|---------------|-------------|---------------|
| Test Details | | | | | | |
| Test | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C |
| Test Code | 5 | 5 | 5 | 5 | 5 | 5 |
| Report Name | HbA1c | HbA1c | HbA1c | HbA1c | HbA1c | HbA1c |
| Unit | % | % | % | % | % | % |
| Decimal Places | 1 | 1 | 1 | 1 | 1 | 1 |
| Wavelength-Primary | 660 | 660 | 660 | 660 | 660 | 660 |
| Wavelength-Secondary | 0 | 0 | 0 | 0 | 0 | 0 |
| Assay type | 2-Point | 2-Point | 2-Point | 2-Point | 2-Point | 2-Point |
| Curve type | Exponential | Exponential | Exponential | Exponential | Exponential | Exponential |
| M1 Start | 20 | 18 | 17 | 27 | 16 | 18 |
| M1 End | 20 | 18 | 17 | 27 | 16 | 18 |
| M2 Start | 34 | 36 | 50 | 62 | 31 | 34 |
| M2 End | 34 | 36 | 50 | 62 | 31 | 34 |
| Sample replicates | 1 | 1 | 1 | 1 | 1 | 1 |
| Standard replicates | 3 | 3 | 3 | 3 | 3 | 3 |
| Control replicates | 1 | 1 | 1 | 1 | 1 | 1 |
| Control interval | 0 | 0 | 0 | 0 | 0 | 0 |
| Reaction Direction | Increasing | Increasing | Increasing | Increasing | Increasing | Increasing |
| React. Abs. Limit | 0 | 0 | 0 | 0 | 0 | 0 |
| Prozone Limit % | 0 | 0 | 0 | 0 | 0 | 0 |
| Prozone Check | Lower | Lower | Lower | Lower | Lower | Lower |
| Linearity Limit % | 0 | 0 | 0 | 0 | 0 | 0 |
| Delta Abs/Min | 0 | 0 | 0 | 0 | 0 | 0 |
| Technical Minimum | NA | NA | NA | NA | NA | NA |
| Technical Maximum | NA | NA | NA | NA | NA | NA |
| Y=aX+b | | | | | | |
| a= | 1 | 1 | 1 | 1 | 1 | 1 |
| b= | 0 | 0 | 0 | 0 | 0 | 0 |
| Reagent Abs Min | NA | NA | NA | NA | NA | NA |
| Reagent Abs Max | 0 | 0 | 0 | 0 | 0 | 0 |
| Auto Rerun | No | No | No | No | No | No |
| Total Reagents | 2 | 2 | 2 | 2 | 3 | 2 |
| Reagent R1 | HBA1C R1 | HBA1C R1 | HBA1C R1 | HBA1C R1 | HBA1C R1 | HBA1C R1 |
| Reagent R2 | HBA1C R2a+R2b | HBA1C R2a+R2b | HBA1C R2a+R2b | HBA1C R2a+R2b | HBA1C R2a | HBA1C R2a+R2b |
| Reagent R3 | NA | NA | NA | NA | HBA1C R2b | NA |

Test Volumes

| | | | | | | |
|-----------------------|-------|-------|-------|-------|-------|-------|
| Test | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C |
| Sample Type | BLOOD | BLOOD | BLOOD | BLOOD | BLOOD | BLOOD |
| Sample Volumes | | | | | | |
| Normal | 3.5 | 3.5 | 3.5 | 3 | 3 | 3.5 |
| Dilution Ratio | 1 | 1 | 1 | 1 | 1 | 1 |
| Increase | 7 | 7 | 7 | 6 | 6 | 7 |
| Dilution Ratio | 1 | 1 | 1 | 1 | 1 | 1 |
| Decrease | 2 | 2 | 2 | 2 | 2 | 2 |
| Dilution Ratio | 1 | 1 | 1 | 1 | 1 | 1 |
| Standard volume | 3.5 | 3.5 | 3.5 | 3 | 3 | 3.5 |

Reagent Volumes and Stirrer speed

| | | | | | | |
|------------------|------|------|------|------|------|------|
| RGT-1 Volume | 150 | 140 | 140 | 120 | 120 | 150 |
| R1 Stirrer Speed | High | High | High | High | High | High |
| RGT-2 Volume | 70 | 70 | 70 | 60 | 40 | 70 |
| R2 Stirrer Speed | High | High | High | High | High | High |
| RGT-3 Volume | 0 | 0 | 0 | 0 | 20 | 0 |
| R3 Stirrer Speed | NA | NA | NA | NA | High | NA |

Reference Ranges

| | | | | | | |
|------------------------|---------|---------|---------|---------|---------|---------|
| Test | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C |
| Sample Type | BLOOD | BLOOD | BLOOD | BLOOD | BLOOD | BLOOD |
| Reference Range | Default | Default | Default | Default | Default | Default |
| Category Male | | | | | | |
| Normal-Lower Limit | 4 | 4 | 4 | 4 | 4 | 4 |
| Normal-Upper Limit | 6 | 6 | 6 | 6 | 6 | 6 |
| Panic-Lower Limit | NA | NA | NA | NA | NA | NA |
| Panic-Upper Limit | NA | NA | NA | NA | NA | NA |
| Category Female | | | | | | |
| Normal-Lower Limit | 4 | 4 | 4 | 4 | 4 | 4 |
| Normal-Upper Limit | 6 | 6 | 6 | 6 | 6 | 6 |
| Panic-Lower Limit | NA | NA | NA | NA | NA | NA |
| Panic-Upper Limit | NA | NA | NA | NA | NA | NA |

Revision Number

| | | | | | | |
|----------|------------------------------------|------------------------------------|--|------------------------------------|-------------------------------------|------------------------------------|
| Revision | <A-100- HBA1C-4 29. 1. 2019> | <A-200- HBA1C-3 29. 1. 2019> | <A-300/600- HBA1C-3 29. 1. 2019> | <A-640- HBA1C-3 29. 1. 2019> | <A-1000- HBA1C-3 29. 1. 2019> | <A-180- HBA1C-2 29. 1. 2019> |
|----------|------------------------------------|------------------------------------|--|------------------------------------|-------------------------------------|------------------------------------|

ASSAY PARAMETERS (SI units)

| Instrument | XL-100 EM-100 | XL-200 EM-200 | XL-300/600 EM-360 | XL-640 | XL-1000 | XL-180 |
|----------------------|------------------|------------------|----------------------|---------------|-------------|---------------|
| Test Details | | | | | | |
| Test | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C |
| Test Code | 5 | 5 | 5 | 5 | 5 | 5 |
| Report Name | HbA1c | HbA1c | HbA1c | HbA1c | HbA1c | HbA1c |
| Unit | mmol/mol | mmol/mol | mmol/mol | mmol/mol | mmol/mol | mmol/mol |
| Decimal Places | 1 | 1 | 1 | 1 | 1 | 1 |
| Wavelength-Primary | 660 | 660 | 660 | 660 | 660 | 660 |
| Wavelength-Secondary | 0 | 0 | 0 | 0 | 0 | 0 |
| Assay type | 2-Point | 2-Point | 2-Point | 2-Point | 2-Point | 2-Point |
| Curve type | Exponential | Exponential | Exponential | Exponential | Exponential | Exponential |
| M1 Start | 20 | 18 | 17 | 27 | 16 | 18 |
| M1 End | 20 | 18 | 17 | 27 | 16 | 18 |
| M2 Start | 34 | 36 | 50 | 62 | 31 | 34 |
| M2 End | 34 | 36 | 50 | 62 | 31 | 34 |
| Sample replicates | 1 | 1 | 1 | 1 | 1 | 1 |
| Standard replicates | 3 | 3 | 3 | 3 | 3 | 3 |
| Control replicates | 1 | 1 | 1 | 1 | 1 | 1 |
| Control interval | 0 | 0 | 0 | 0 | 0 | 0 |
| Reaction Direction | Increasing | Increasing | Increasing | Increasing | Increasing | Increasing |
| React. Abs. Limit | 0 | 0 | 0 | 0 | 0 | 0 |
| Prozone Limit % | 0 | 0 | 0 | 0 | 0 | 0 |
| Prozone Check | Lower | Lower | Lower | Lower | Lower | Lower |
| Linearity Limit % | 0 | 0 | 0 | 0 | 0 | 0 |
| Delta Abs/Min | 0 | 0 | 0 | 0 | 0 | 0 |
| Technical Minimum | NA | NA | NA | NA | NA | NA |
| Technical Maximum | NA | NA | NA | NA | NA | NA |
| Y=aX+b | | | | | | |
| a= | 1 | 1 | 1 | 1 | 1 | 1 |
| b= | 0 | 0 | 0 | 0 | 0 | 0 |
| Reagent Abs Min | NA | NA | NA | NA | NA | NA |
| Reagent Abs Max | 0 | 0 | 0 | 0 | 0 | 0 |
| Auto Rerun | No | No | No | No | No | No |
| Total Reagents | 2 | 2 | 2 | 2 | 3 | 2 |
| Reagent R1 | HBA1C R1 | HBA1C R1 | HBA1C R1 | HBA1C R1 | HBA1C R1 | HBA1C R1 |
| Reagent R2 | HBA1C R2a+R2b | HBA1C R2a+R2b | HBA1C R2a+R2b | HBA1C R2a+R2b | HBA1C R2a | HBA1C R2a+R2b |
| Reagent R3 | NA | NA | NA | NA | HBA1C R2b | NA |

Test Volumes

| | | | | | | |
|-----------------------|-------|-------|-------|-------|-------|-------|
| Test | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C | HBA1C |
| Sample Type | BLOOD | BLOOD | BLOOD | BLOOD | BLOOD | BLOOD |
| Sample Volumes | | | | | | |
| Normal | 3.5 | 3.5 | 3.5 | 3 | 3 | 3.5 |
| Dilution Ratio | 1 | 1 | 1 | 1 | 1 | 1 |
| Increase | 7 | 7 | 7 | 6 | 6 | 7 |
| Dilution Ratio | 1 | 1 | 1 | 1 | 1 | 1 |
| Decrease | 2 | 2 | 2 | 2 | 2 | 2 |
| Dilution Ratio | 1 | 1 | 1 | 1 | 1 | 1 |
| Standard volume | 3.5 | 3.5 | 3.5 | 3 | 3 | 3.5 |

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1. Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 142-48.
2. Sacks DB. Carbohydrates. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B. Saunders Company; 1999.p. 790-6.
3. Jeppsson JO, Kobold U, Barr J, Finke A et al. Approved IFCC reference method for the measurement of HbA1c in human blood. Clin Chem Lab Med 2002;40:78-89.
4. Hoelzel W, Weykamp C et al. IFCC Reference System for Measurement of Hemoglobin A1c in Human Blood and the National Standardization Schemes in the United States, Japan, and Sweden: A Method-Comparison Study. Clin Chem 2004; 50:1:166-74.
5. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes in the development and progression of long-term complications in insulin- -dependent diabetes mellitus. N Engl J Med.1993;329:977-86.
6. Little RR, Rohlfing CL, Wiedmeyer HM, Myers GL et al. The National Glycohemoglobin Standardization Program: A Five-Years Progress Report. Clin Chem 2001;47:1985-92.
7. Miedema K. Standardization of HbA1c and Optimal Range of Monitoring. Scand J Clin Lab Invest 2005;65 (Suppl 240):61-72.

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