HDL/LDL CAL

For *in vitro* diagnostic use only.

Cat. No.:	Pack Name:	Packaging (Content):
XSYS0061	HDL/LDL CAL	2 × 1 ml

INTENDED USE

HDL/LDL Calibrator is intended for quantitative determination of HDL and LDL cholesterol using a HDL and LDL Direct reagent.

CERTIFIED VALUES

Calibrator LOT	Expiry	Parameter	Concentration	
xxxxxxxxx	rrrr-mm	HDL	xx.x mg/dl	x.xx mmol/l
		LDL	xxx mg/dl	x.xx mmol/l

UNIT CONVERSION

mg/dl × 0.026 = mmol/l

REAGENTS SUPPLIED

HDL/LDL CAL for 2 × 1 ml

REAGENT STORAGE AND STABILITY

Unopened reagents are stable till the expiry date mentioned on the label when stored at 2–8 °C. The reconstituted calibrator is stable for 6 days at 2–8 °C.

CALIBRATOR RECONSTITUTION

Reconstitute with 1ml of deionised water at 20–25 °C and mix gently (avoid foaming).

Allow to stand for at least 30 minutes until complete reconstitution before use.

Store reconstituted calibrator at 2–8 °C.

WARNING:

Serum used in the manufacture of this product has been tested and found non reactive for hepatitis B surface antigen (HbsAg), antibody to Hepatitis C (HCV), HIV-1 p24 antigen, and antibody to HIV-1/2. The test procedures do not guarantee that all infectious agents will be detected. Because no test method can offer complete assurance that Hepatitis B virus Hepatitis C virus and HIV-1/2 or other infectious agents are absent, the material should be handled as potentially infectious.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent is not classified as dangerous.

USED SYMBOLS

Catalogue number

Manufacturer

See instructions for use

Lot number


In vitro diagnostic medical device

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