LIPASE

Cat. No.	Pack Name	Packaging (Content)
XSYS0081	LIP 110	R1: 2 x 44 ml, R2: 2 x 11 ml

INTENDED USE

Diagnostic reagent for quantitative *in vitro* determination of Lipase in human serum and plasma.

CLINICAL SIGNIFICANCE

Lipases are enzymes which hydrolyze glycerol ester of long fatty acids. The enzyme and its cofactor colipase is produced in the pancreas, lipase being also secreted in small amounts by the salivary glands as well as by gastric, pulmonary and intestinal mucosa. Bile acids and colipase form micellar complexes with the lipids and bind lipase on the substrate / water interface. Determination of lipase is used for investigation of pancreatic disorders. In acute pancreatitis the lipase concentrations rise to 2–50 fold to upper reference limit within 4–8 hours after begin of abdominal pain peaking at 24 hours and decreasing within 8 to 14 days. Elevated lipase values can also be observed in chronic pancreatitis and obstruction of the pancreatic duct.

PRINCIPLE

Enzymatic color test.

The colorimetric substrate 1,2-o-dilauryl-rac-glycero-3-glutaric acid-(6-methylresorufin)-ester is cleaved by pancreatic lipase and the resulting dicarboxylic acid ester is hydrolysed under the alkaline test conditions to yield the chromophore methylresorufin. The kinetic of colour formation at 580 nm is monitored and it is proportional to lipase activity in sample.

REAGENT COMPOSITION

R1

Good's Buffer pH 8.0	
Taurodesoxycholate	≥ 1 mmol/l
Desoxycholate	≥ 1 mmol/l
Calcium ions	≥ 1 mmol/l
Colipase	≥ 1 mg/l

R2

Tartrate Buffer pH 4.0	
Lipase Color Substrate	≥ 0.1 mmol/l

REAGENT PREPARATION

Reagents are liquid, ready to use.

STABILITY AND STORAGE

The unopened reagents are stable till the expiry date stated on the bottle and kit label when stored at 2–8 °C.

On board stability: min. 30 days if refrigerated (2–10 °C) and not contaminated.

Caution: reagent R2 is a microemulsion. Therefore, a slight apparent precipitation could occur, showing a light red deposit on the bottom of vial. It is a normal behaviour and it is recommended to resuspend solution before analysis with a mild shaking.

SPECIMEN COLLECTION & HANDLING

It is recommended to follow NCCLS procedures (or similar standardized conditions) Use serum, plasma (heparin).

Stability

in serum / plasma: 7 days at 20–25 °C
7 days at 4–8 °C
1 year at -20 °C

Discard contaminated specimens.

CALIBRATION

Calibration with calibrator XL MULTICAL 4x3, Cat. No. XSYS0034 or XL MULTICAL 10x3, Cat. No. XSYS0122 is recommended.

Calibration frequency: it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

Traceability:

This calibrator has been standardized to IFCC formulation.

QUALITY CONTROL

For quality control ERBA NORM 4x5, Cat. No. BLT00080 or ERBA NORM 10x5, Cat. No. XSYS0123 and ERBA PATH 4x5, Cat. No. BLT00081 or ERBA PATH 10x5, Cat. No. XSYS0124 are recommended.

CALCULATION

The XL Results are calculated automatically by the instrument

UNIT CONVERSION

U/l x 0.017 = µkat/l

EXPECTED VALUES ³

Serum

≤ 60 U/l (≤ 1.0 µkat/l)

It is recommended that each laboratory verify this range or derives reference interval for the population it serves.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Limit of quantification: 8.2 U/l

Linearity: 300 U/l

Measuring range: 8.2–300 U/l

PRECISION

Intra-assay precision Within run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	209.4	4.38	2.09
Sample 2	85.2	2.82	3.29

Inter-assay precision Run to run (n=20)	Mean (U/l)	SD (U/l)	CV (%)
Sample 1	52.8	1.92	3.64
Sample 2	171	2.76	2.5

COMPARISON

A comparison between XL-Systems Lipase (y) and a commercially available test (x) using 40 samples gave following results:

r = 1.000

y = 1.065 x + 1.98 U/l

INTERFERENCES

Following substances do not interfere:

haemoglobin up to 4.5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 1000 mg/dl. N-acetyl-p-benzoquinone imine (metabolite of paracetamol) could generate erroneously low results in samples for patients that have taken toxic doses of paracetamol.

WARNING AND PRECAUTIONS

For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.

Hazards identification in accordance with Regulation (EC) No 1272/2008

Reagent 1 is not classified as dangerous.

Reagent 2 of the kit contains propan-1-ol.

Danger

Hazard statement:

H318 Causes serious eye damage.

Precautionary statement:

P280 Wear protective gloves/protective clothing/eye protection.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

WASTE MANAGEMENT

Please refer to local legal requirements.

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1. Lorentz K Lipase. In: Thomas L, editor. Clinical laboratory diagnostics. 1st ed. Frankfurt: TH- Books Verlagsgesellschaft; 1998. p. 95-7.
2. Moss DW, Henderson AR. Digestive enzymes of pancreatic origin. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 689-708.
3. Tietz N, Shuey DF. Lipase in serum – the elusive enzyme: an overview. Clin Chem 1993; 39:746- 56.
4. Lott J, Patel ST, Sawhney AK, Kazmierczak SC, LoveJE. Assays of serum lipase: analytical and clinical considerations. Clin Chem 1986;32:1290-1302.
5. Leybold A, Junge W. Importance of colipase for the measurement of serum lipase activity. Adv Clin Enzymol 1986;4:60-7.
6. Borgström B. The action of bile salts and other detergents on pancreatic lipase and the interaction with colipase. Biochimica et Biophysica Acta 1977;488:381-91.
7. Garguori Y, Julien R, Bois A, Verger R, Sarda L. Studies on the detergent inhibition of pancreatic lipase activity. J of Lipid Research 1983;24:1336-42.
8. Guder WG, Zafta B et al. The quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag;2001 p.36-7.
9. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2 Washington, DC: The American Association for Clinical Chemistry Press 2000.
10. Junge W, Abicht K, Goldman J. Evaluation of the colorimetric liquid assay for pancreatic lipase on Hitachi analysers in 7 clinical centres in Europe. Clin Chem Lab Med 1999; 37,Special suppl:469.

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