FERRITIN

Cat. No.	Pack Name	Packaging (Content)
XSYS0101	FRTN	R1: 2 × 14.5 ml (Buffer), R2: 2 × 7.7 ml (Latex)

EN

INTENDED USE

Quantitative determination of Ferritin in human serum and plasma by immunoturbidimetry.

CLINICAL SIGNIFICANCE

The plasma Ferritin concentration declines very early in the development of iron deficiency. On the other hand, a large number of chronic diseases result in increased serum Ferritin concentrations. These diseases include chronic infections, chronic inflammatory disorders such as rheumatoid arthritis or renal disease, Gaucher's disease, and numerous types of malignancies, especially lymphomas, leukaemias, breast cancer and neuroblastoma. Increase in plasma Ferritin concentration also occurs in viral hepatitis or following toxic liver injury because of release of Ferritin from damaged liver cells. Plasma Ferritin concentration is also increased with increases of iron stores, as seen in patients with hemosiderosis or hemochromatosis. Besides the use of Ferritin as an iron metabolism parameter, Ferritin as also gained importance as a tumour marker for therapeutic drug monitoring and follow-up.

PRINCIPLE

The assay is based on a latex agglutination reaction. The measured difference of absorbances at 600 nm by the 2-point method is proportional to the concentration of ferritin in the sample.

REAGENT COMPOSITION

R1 (Buffer)
TRIS Buffer (pH7.2) 120mM
Non reactant components and preservatives

R2 (Latex)
Latex particles coated with rabbit antibodies against human Ferritin.
Non reactant components and preservatives.

REAGENT PREPARATION

Liquid reagents, ready to use. R2 needs to be mixed via inversion 5–10 times before being placed on the instrument and the stirring must be repeated on a weekly basis.

STABILITY AND STORAGE

The reagents are stable until expiry date when kept at 2–8 °C. Stability in the instrument is at least 8 weeks if contamination is avoided. Do not freeze.

SAMPLE COLLECTION

Serum, heparinised or EDTA-collected plasma. Ferritin is stable in serum/plasma up to 7 days at 2–8 °C and up to 6 months at -20 °C.

MATERIALS REQUIRED BUT NOT PROVIDED

- Any instrument with temperature control of 37 ± 0.5 °C that is capable of reading absorbance accurately at 600 nm may be used
- Analyser specific consumables such as sample cups
- Calibrator
- Controls
- Saline (9 g/l NaCl)

ASSAY PROCEDURE

Refer to the assay parameters for details.

CALIBRATION

Blank: saline

Cat. No.	Product name	Pack name	Content
XSYS0104	FERRITIN CALIBRATOR	FRTN CAL	1 × 1 ml

Calibration curve: generate a 5 point calibration curve by diluting the calibrator 1:16, 1:8, 1:4, 1:2 and undiluted in saline.

Calibration frequency it is recommended to do a calibration

- after reagent lot change
- as required by internal quality control procedures

QUALITY CONTROL

For quality control use

Cat. No.	Product name	Pack name	Content
BLT20042	FERRITIN CONTROL LOW	FRTN CON L	1 × 1 ml
BLT20041	FERRITIN CONTROL HIGH	FRTN CON H	1 × 1 ml

CALCULATION

Results are calculated automatically by the instrument.

EXPECTED VALUES

Reference Values

Newborn	25–200 ng/ml (µg/l)
1 month	200–600 ng/ml (µg/l)
2–5 months	50–200 ng/ml (µg/l)
6 months – 15 years	7–140 ng/ml (µg/l)
Adult	
Male	20–250 ng/ml (µg/l)
Female	10–120 ng/ml (µg/l)

Reference values are given for orientation only. Each laboratory should establish its own reference values.

PERFORMANCE DATA

Data contained within this section is representative of performance on ERBA XL systems. Data obtained in your laboratory may differ from these values.

Measuring Range:	5.0–420 ng/ml (µg/l)
Limit of quantification:	5.0 ng/ml (µg/l)
Hook Effect:	No risk

PRECISION

Intra-assay precision Within run (n = 20)	Mean (ng/ml)	SD (ng/ml)	CV (%)
Sample 1	62.1	1.9	3.11
Sample 2	118.0	4.3	3.66

Inter-assay precision Run to run (n = 20)	Mean (ng/ml)	SD (ng/ml)	CV (%)
Sample 1	47.1	1.6	3.51
Sample 2	143.1	6.7	4.71

COMPARISON

A comparison between XL-Systems FERRITIN (y) and a commercially available test (x) using 53 samples gave following results:
 $y = 0.923 x - 3.140$ ng/ml
 $r = 0.996$

INTERFERENCES

Following substances do not interfere:
haemoglobin up to 5 g/l, bilirubin up to 40 mg/dl, triglycerides up to 3000 mg/dl

WARNING AND PRECAUTIONS

1. For *in vitro* diagnostic use. To be handled by entitled and professionally educated person.
2. Product (R2) contains sodium azide < 0.1 %.
3. Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV1 and HIV2 antibodies, as well as for the hepatitis B surface antigen and anti-hepatitis C antibodies, using a method approved by the FDA.

R1

Warning

Contains:
reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one (3:1)

Hazard statement:

H317 May cause an allergic skin reaction.

Precautionary statement:

P280 Wear protective gloves.
P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

R2

Supplemental information:

EUH 208 Contains 2-chloroacetamide. May produce an allergic reaction

WASTE MANAGEMENT

Please refer to local legal requirements.

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ASSAY PARAMETERS (conventional units)

Instrument	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details					
Test	FERR	FERR	FERR	FERR	FERR
Test Code	3	3	3	3	3
Report Name	Ferritin	Ferritin	Ferritin	Ferritin	Ferritin
Unit	ng/ml	ng/ml	ng/ml	ng/ml	ng/ml
Decimal Places	1	1	1	1	1
Wavelength-Primary	600	600	600	600	600
Wavelength-Secondary	0	0	0	0	0
Assay type	2-Point	2-Point	2-Point	2-Point	2-Point
Curve type	Cubic Spline	Cubic Spline	Cubic Spline	Cubic Spline	Cubic Spline
M1 Start	19	17	29	14	19
M1 End	19	17	29	14	19
M2 Start	36	48	63	31	34
M2 End	36	48	63	31	34
Sample replicates	1	1	1	1	1
Standard replicates	3	3	3	3	3
Control replicates	1	1	1	1	1
Control interval	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	NA	NA	NA	NA	NA
Prozone Limit %	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0
Technical Minimum	NA	NA	NA	NA	NA
Technical Maximum	NA	NA	NA	NA	NA
Y=aX+b					
a=	1	1	1	1	1
b=	0	0	0	0	0
Reagent Abs Min	NA	NA	NA	NA	NA
Reagent Abs Max	0	0	0	0	0
Auto Rerun	No	No	No	No	No
Total Reagents	2	2	2	2	2
Reagent R1	FERR R1	FERR R1	FERR R1	FERR R1	FERR R1
Reagent R2	FERR R2	FERR R2	FERR R2	FERR R2	FERR R2
Reagent R3	NA	NA	NA	NA	NA

Test Volumes					
Test	FERR	FERR	FERR	FERR	FERR
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM
Sample Volumes					
Normal	20	20	20	16.7	25
Dilution Ratio	1	1	1	1	1
Increase	30	30	30	25	37
Dilution Ratio	1	1	1	1	1
Decrease	2	2	2	2	2
Dilution Ratio	1	1	1	1	1
Standard volume	20	20	20	16.7	25
Reagent Volumes and Stirrer speed					
RGT-1 Volume	120	120	120	100	150
R1 Stirrer Speed	High	NA	High	High	High
RGT-2 Volume	60	60	60	50	75
R2 Stirrer Speed	High	NA	High	High	High
RGT-3 Volume	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA

Reference Ranges					
Test	FERR	FERR	FERR	FERR	FERR
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default
Category Male					
Normal-Lower Limit	20	20	20	20	20
Normal-Upper Limit	250	250	250	250	250
Panic-Lower Limit	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA
Category Female					
Normal-Lower Limit	10	20	10	10	10
Normal-Upper Limit	120	250	120	120	120
Panic-Lower Limit	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA

Revision Number					
Revision	<A-200-FERR-1 18.02.2020>	<A-300/600-FERR-1 25.01.2021>	<A-640-FERR-1 18.02.2020>	<A-1000-FERR-1 18.02.2020>	<A-180-FERR-2 21.05.2020>

ASSAY PARAMETERS (SI units)

Instrument	XL-200 EM-200	XL-300/600 EM-360	XL-640	XL-1000	XL-180
Test Details					
Test	FERR	FERR	FERR	FERR	FERR
Test Code	3	3	3	3	3
Report Name	Ferritin	Ferritin	Ferritin	Ferritin	Ferritin
Unit	µg/l	µg/l	µg/l	µg/l	µg/l
Decimal Places	1	1	1	1	1
Wavelength-Primary	600	600	600	600	600
Wavelength-Secondary	0	0	0	0	0
Assay type	2-Point	2-Point	2-Point	2-Point	2-Point
Curve type	Cubic Spline	Cubic Spline	Cubic Spline	Cubic Spline	Cubic Spline
M1 Start	19	17	29	14	19
M1 End	19	17	29	14	19
M2 Start	36	48	63	31	34
M2 End	36	48	63	31	34
Sample replicates	1	1	1	1	1
Standard replicates	3	3	3	3	3
Control replicates	1	1	1	1	1
Control interval	0	0	0	0	0
Reaction Direction	Increasing	Increasing	Increasing	Increasing	Increasing
React. Abs. Limit	NA	NA	NA	NA	NA
Prozone Limit %	0	0	0	0	0
Prozone Check	Lower	Lower	Lower	Lower	Lower
Linearity Limit %	0	0	0	0	0
Delta Abs/Min	0	0	0	0	0
Technical Minimum	NA	NA	NA	NA	NA
Technical Maximum	NA	NA	NA	NA	NA
Y=aX+b					
a=	1	1	1	1	1
b=	0	0	0	0	0
Reagent Abs Min	NA	NA	NA	NA	NA
Reagent Abs Max	0	0	0	0	0
Auto Rerun	No	No	No	No	No
Total Reagents	2	2	2	2	2
Reagent R1	FERR R1	FERR R1	FERR R1	FERR R1	FERR R1
Reagent R2	FERR R2	FERR R2	FERR R2	FERR R2	FERR R2
Reagent R3	NA	NA	NA	NA	NA

Test Volumes					
Test	FERR	FERR	FERR	FERR	FERR
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM
Sample Volumes					
Normal	20	20	20	16.7	25
Dilution Ratio	1	1	1	1	1
Increase	30	30	30	25	37
Dilution Ratio	1	1	1	1	1
Decrease	2	2	2	2	2
Dilution Ratio	1	1	1	1	1
Standard volume	20	20	20	16.7	25
Reagent Volumes and Stirrer speed					
RGT-1 Volume	120	120	120	100	150
R1 Stirrer Speed	High	NA	High	High	High
RGT-2 Volume	60	60	60	50	75
R2 Stirrer Speed	High	NA	High	High	High
RGT-3 Volume	0	0	0	0	0
R3 Stirrer Speed	NA	NA	NA	NA	NA

Reference Ranges					
Test	FERR	FERR	FERR	FERR	FERR
Sample Type	SERUM	SERUM	SERUM	SERUM	SERUM
Reference Range	Default	Default	Default	Default	Default
Category Male					
Normal-Lower Limit	20	20	20	20	20
Normal-Upper Limit	250	250	250	250	250
Panic-Lower Limit	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA
Category Female					
Normal-Lower Limit	10	10	10	10	10
Normal-Upper Limit	120	120	120	120	120
Panic-Lower Limit	NA	NA	NA	NA	NA
Panic-Upper Limit	NA	NA	NA	NA	NA

Revision Number					
Revision	<ASI-200-FERR-1 18.02.2020>	<ASI-300/600-FERR-1 25.01.2021>	<ASI-640-FERR-1 18.02.2020>	<ASI-1000-FERR-1 18.02.2020>	<ASI-180-FERR-2 21.05.2020>

REFERENCES / ЛИТЕРАТУРА / LITERATURA

1. Tietz, NW, ed. Clinical Guide to Laboratory Tests. 3rd. ed. Philadelphia: W.B. Saunders Company Ltd., 1995.
2. Burtis, C.A., Ashwood, E.R., Bruns, D.E. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 5th edition: W.B. Saunders Company Ltd., 2012.
3. Jacobs, DJ, Demott, WR, Grady, HJ, Horvat, RJ, Huestis, DW and Kasten, BL, JR, eds. Laboratory Test Handbook. 4th. ed. Ohio, Hudson: Lexi-Comp Inc., 1996.
4. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests. 2nd. ed. Washington, DC: The American Association for Clinical Chemistry Press, 1997.
5. Bernard A., Lauwerys R. Turbidimetric latex immunoassay for serum ferritin. Journal of Immunological methods 1984; 71:141-147.

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